Soluciones Dentales

3M™ Filtek™ P60

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7 Agosto 2025

Ficha Técnica

Descripción:

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- Prepare Tooth and Isolate.**
Apply liner/base if desired.
- Mix a level scoop of Vitrebond powder and a drop of Vitrebond liquid on a mixing pad.
 - Apply a thin layer of the liner/base to dentin surfaces using a ball applicator.
 - Light cure for 30 seconds.
- Etch:**
- Apply 3M™ Scotchbond™ etchant to enamel and dentin. Wait 15 seconds. Etchant on Vitrebond base is not deleterious.
 - Rinse.
 - Blot excess water, leaving tooth moist.

- Bond:**
- Using a fully saturated brush tip for each coat, apply 2 consecutive coats of 3M Single Bond adhesive to enamel and dentin.
 - Dry gently for 2-5 seconds.
 - Light cure for 10 seconds.

- Place Restorative:**
- Place 3M Filtek P60 restorative in increments less than 2.5mm.
 - Light cure each increment for 20 seconds.

- Finish and Polish:**
- Finish occlusal surface using an appropriate finishing instrument.
 - Finish interproximal surfaces with 3M™ Sof-Lex™ Pop-on™ (extra-thin discs and Sof-Lex strips).

- Check Occlusion:**
- Check lateral and centric occlusion.
 - Adjust if necessary.

Please refer to instructions for more detailed information as well as precautionary and safety information.
3M Customer Review 1-800-454-2299 © 1999 3M

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Cod Munnich	Producto	Presentación
527816	FILTEK P60 A3	JERINGA DE 4 GR
527817	FILTEK P60 B2	JERINGA DE 4 GR
527818	FILTEK P60 C2	JERINGA DE 4 GR

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