Soluciones Dentales

3M™ Filtek™ P60

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7 Agosto 2025

Ficha Técnica

Descripción:

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3M™ Direct Posterior Restorations


3M™ Vitrebond™ Light Cure Glass Ionomer Liner/Base

3M™ Single Bond Dental Adhesive System


3M™ Filtek™ P60 Posterior Restorative

Prepare Tooth and Isolate.
Apply liner/base if desired.
• Mix a level scoop of Vitrebond powder and a drop of Vitrebond liquid on a mixing pad.
• Apply a thin layer of the liner/base to dentin surfaces using a ball applicator.
• Light cure for 30 seconds.

Etch:
• Apply 3M™ Scotchbond™ etchant to enamel and dentin. Wait 15 seconds. Etchant on Vitrebond base is not deleterious.
• Rinse.
• Blot excess water, leaving tooth moist.

Bond:
• Using a fully saturated brush tip for each coat, apply 2 consecutive coats of 3M Single Bond adhesive to enamel and dentin.
• Dry gently for 2-5 seconds.
• Light cure for 10 seconds.

Place Restorative:
• Place 3M Filtek P60 restorative in increments less than 2.5mm.
• Light cure each increment for 20 seconds.

Finish and Polish:
• Finish occlusal surface using an appropriate finishing instrument.
• Finish interproximal surfaces with 3M™ Sof-Lex™ Pop-on™ (extra-thin discs and Sof-Lex strips).

Check Occlusion:
• Check lateral and centric occlusion.
• Adjust if necessary.

Please refer to instructions for more detailed information as well as precautionary and warranty information.
3M Customer Review 1-800-454-2299 © 1999 3M

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Cod Munnich	Producto	Presentación
527816	FILTEK P60 A3	JERINGA DE 4 GR
527817	FILTEK P60 B2	JERINGA DE 4 GR
527818	FILTEK P60 C2	JERINGA DE 4 GR

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