**Boston
Scientific**

Hurricane[™] RX

BILIARY

Balloon Dilatation Catheter

Directions for Use	2
Instrucciones de uso	7
Mode d'emploi	12
Gebrauchsanweisung	17
Istruzioni per l'uso	23
Gebruiksaanwijzing	28
Instruções de Utilização	33

91055637-01

2015-06

TABLE OF CONTENTS

WARNING 3

DEVICE DESCRIPTION..... 3

 Balloon Construction 3

 Catheter Construction 3

INTENDED USE/INDICATIONS FOR USE 4

CONTRAINDICATIONS..... 4

PRECAUTIONS..... 4

ADVERSE EVENTS 4

HOW SUPPLIED..... 4

HANDLING AND STORAGE 4

OPERATIONAL INSTRUCTIONS..... 4

 Inspection and Preparation..... 4

 Catheter Insertion 5

 Balloon Inflation 5

WARNINGS..... 5

 Catheter Withdrawal 6

WARRANTY 6

Hurricane™ RX

BILIARY

Balloon Dilatation Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Hurricane RX Biliary Balloon Dilatation Catheter, is a double lumen catheter with a balloon mounted at the distal tip. Dilatation balloon catheters are used to exert radial force to dilate narrow duct segments, as well as the Sphincter of Oddi.

Balloon Construction

Balloons are constructed with specially treated polymers and processes, which provide maximum strength. Each balloon inflates to its stated diameter and length at a specific pressure. Maximum pressure is dependent on the balloon size and is identified on the label.

CAUTION: DO NOT EXCEED THE MAXIMUM PRESSURE PRINTED ON THE LABEL.

Catheter Construction

The catheter is multi-lumen. The guidewire lumen begins at the distal tip of the catheter, travels under the black sleeve where it exits at the proximal end of the black sleeve, which is 25 cm from the catheter tip. This lumen is used to pass the catheter over a 0.035 in (0.89 mm) guidewire to provide for the rapid exchange feature. The lumen marked “**Balloon**” is the balloon inflation lumen. The catheter shaft tapers beneath the balloon segment to achieve the lowest possible deflated profile. The lumen marked “**Injection**”, at the proximal end of the catheter, can be used for contrast medium injection for fluoroscopic visualization of the bile duct. There is a FloSwitch™ Adapter connected to the hub of the “Injection” lumen. The FloSwitch Adapter is to remain closed until distal dye contrast injections are required.

Two radiopaque markers are placed under the balloon segment of the catheter to provide visual reference points for balloon positioning within the stricture or duct. A stainless steel stiffening stylet is inserted in the proximal portion of the catheter shaft to augment catheter pushability through the endoscope and should be left in place. This stylet ends under the black tunnel.

INTENDED USE/INDICATIONS FOR USE

The Hurricane™ RX Biliary Balloon Dilatation Catheter is recommended for endoscopic dilatation of strictures of the biliary tree and the Sphincter of Oddi.

The Hurricane RX Biliary Balloon Dilatation Catheter may be used for injection of contrast medium for fluoroscopic visualization of the bile ducts.

CONTRAINDICATIONS

None known.

PRECAUTIONS

A thorough understanding of the technical principles, clinical applications and risks associated with balloon dilatation of the biliary tract is necessary before using this product.

Any use of this product for procedures other than those indicated in these instructions is not recommended.

ADVERSE EVENTS

The adverse events that may result from a balloon dilatation procedure include but may not be limited to:

- perforation
- hemorrhage
- pain and tenderness
- sepsis/infection
- allergic reaction to contrast medium
- inflammation

HOW SUPPLIED

Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible. Do not reuse.

HANDLING AND STORAGE

Store in a cool, dry, dark place. Rotate inventory so that catheters are used prior to the expiration date on the package label.

OPERATIONAL INSTRUCTIONS

Inspection and Preparation

1. Prior to use, carefully examine the unit to verify that the catheter or sterile package has not been damaged in shipment.
2. Remove plastic mandrel from distal end of catheter.

Note: Do not remove stylet from proximal injection hub.

3. Close FloSwitch™ Adapter attached to injection hub.
4. **Prepare balloon inflation syringe with contrast medium*.**
5. Attach balloon inflation hub to inflation syringe with a gauge to monitor pressure.
6. To maximize endoscope passage, apply a vacuum to the catheter before removing protective sleeve from the balloon section of the catheter. Close stopcock, disconnect syringe from stopcock, remove air from syringe and reconnect to balloon hub.
7. Remove protective sleeve from balloon prior to insertion.
8. **For Contrast Injection:** Prime second syringe (10 ml (cc) recommended) with contrast medium* and attach to FloSwitch Adapter, open FloSwitch Adapter and flush if needed.

* **Note:** Commercially available contrast fluids have different viscosities that may affect balloon inflation and deflation performance as well as distal contrast injection performance. The higher the viscosity (thicker) the slower the fluid will move through the catheter. Dilute the contrast fluid with saline as required to improve inflate-deflate time as well as distal contrast injection performance.

PRECAUTION: Do not pre-inflate, pretest balloon, or attempt to re-fold balloon into protective sleeve.

Catheter Insertion

The Hurricane™ RX Biliary Balloon Dilatation Catheter is designed to pass through a duodenoscope with a minimum 3.2 mm working channel. Each catheter's working length is 180 cm. The balloon should be advanced over the guidewire and through the scope channel using short, deliberate 2-3 cm movements to prevent inadvertent kinking of the catheter shaft. Fluoroscopy can be used to confirm proper placement.

1. The Hurricane RX Biliary Balloon Dilatation Catheters are to be used with a Boston Scientific RX locking device.
2. Ensure the RX locking device is firmly seated onto the endoscope against the side of the working channel port.
3. With the 0.035 in (0.89 mm) guidewire positioned and locked into place, backload the Hurricane RX Biliary Balloon Dilatation Catheter over the guidewire. Slide the Hurricane RX Biliary Balloon Dilatation Catheter onto the guidewire and through the first 25 cm until the guidewire exits the proximal end of the black sleeve. Grasp the guidewire and slide the catheter over the guidewire until it reaches the RX locking device.
4. Unlock the guidewire from guidewire locking device. Pass the distal 25 cm section of the catheter past the guidewire locking device while maintaining tension on the guidewire. Lock the guidewire into place and continue advancing the Hurricane RX Biliary Balloon Dilatation Catheter through the scope until cannulation is achieved.
5. Two radiopaque markers located under the dilatation section of the balloon will aid in positioning the catheter while using fluoroscopy.
6. Contrast medium can be injected through the FloSwitch™ Injection Port to visualize the bile duct if needed. The FloSwitch Adapter must be opened prior to injection. The guidewire tip may be pulled back into the catheter to help facilitate contrast injection.

Note: Use caution if retracting wire for aid in contrast injection. Retracting the guidewire tip more than 20 cm into the catheter tip may facilitate loss of guidewire cannulation through the Hurricane RX Biliary Balloon Dilatation Catheter.

Balloon Inflation

Under fluoroscopic guidance, the balloon should be advanced and positioned across the stricture to be dilated. The balloon can be inflated by means of a hand-held syringe. It is recommended that a 10 ml (cc) or larger syringe with a gauge to monitor balloon pressure be used when inflating a Hurricane RX Biliary Balloon Dilatation Catheter. The typical inflation medium is a 1:1 mixture of contrast medium and normal saline.

WARNINGS

- Do not use air or any gaseous substances as a balloon inflation medium.
 - To prevent balloon burst, do not exceed the Maximum Pressure printed on the label.
 - Monitor pressure using an inflation device and pressure gauge. As dilatation takes place, the pressure reading may fluctuate. Adjust the balloon inflation pressure, as necessary to maintain desired pressure.
 - If the balloon does rupture or a significant loss of pressure within the balloon occurs, deflate the balloon and carefully remove the balloon and endoscope together as a unit. Do not attempt to withdraw a ruptured balloon through the endoscope.
-

Catheter Withdrawal

1. When dilatation has been completed, deflate the balloon by applying suction to the balloon lumen.
2. Do not pull back on the catheter until the balloon is completely deflated.
3. Slowly remove the catheter from the endoscope.
4. The catheter should be disposed of in accordance with accepted hospital guidelines.

Note: The larger the syringe diameter, the greater the suction that is applied. For maximum deflation, a 50 ml (cc) syringe is recommended.

PRECAUTION: If resistance is felt when removing a catheter from the scope, stop and remove the endoscope and balloon as one unit.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC’s control directly affect the instrument and the results obtained from its use. BSC’s obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

FloSwitch is a trademark of Argon Medical Devices, Inc. Corporation.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2015 Boston Scientific Corporation or its affiliates.