LightTrail™ Reusable

Laser Fiber

Directions for Use	2
Instrucciones de uso	8
Mode d'emploi	15
Gebrauchsanweisung	22
Istruzioni per l'uso	29
Gebruiksaanwijzing	36
Instruções de Utilização	43

TABLE OF CONTENTS

WARNING3

DEVICE DESCRIPTION3

INTENDED USE/INDICATIONS FOR USE3

CONTRAINDICATIONS.....3

WARNINGS.....3

PRECAUTIONS.....3

ADVERSE EVENTS4

PRODUCT -SPECIFIC SIZE INFORMATION4

HOW SUPPLIED.....4

DIRECTIONS FOR USE.....4

INTRA-OPERATIVE INSTRUCTIONS4

POST-OPERATIVE INSTRUCTIONS5

REPROCESSING (CLEANING, DISINFECTION AND STERILIZATION).....5

WARRANTY7

LightTrail™ Reusable

Laser Fiber

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

After each use, clean, sterilize and reuse LightTrail Reusable Laser Fibers in accordance with hospital, administrative and/or local government policy.

The utilization of cleaning and sterilization practices that do not meet globally accepted reprocessing and sterilization industry standards may create a risk of contamination of the LightTrail Reusable Laser Fibers, potentially resulting in patient infection or cross-infection.

For any equipment used in LightTrail Reusable Laser Fibers device cleaning and sterilization, use only the materials and parameters specified per the equipment manufacturer's manual.

Using a contaminated LightTrail Reusable Laser Fiber may lead to injury, illness, infection and/or sepsis that may result in death of the patient.

DEVICE DESCRIPTION

The LightTrail Reusable laser fibers are supplied EO sterilized. They consist of an SMA-905 connector and a silica core fiber with a polymer jacket. The laser fiber transmits laser energy from the laser console to the treatment site through the fiber tip. An automatic identification feature enables the laser system to recognize the connected laser fiber.

For use with Boston Scientific (StarMedTec) laser systems. Please refer to the laser system User Manual for complete information regarding applications, contraindications, potential adverse events, precautions and warnings.

User Information

Before using the LightTrail Laser Fibers, the physician should fully understand the use of the Holmium (Ho:YAG) and/or Thulium (Tm:YAG) laser safety considerations, tissue interaction and proper technique specific to the treatment for which the physician intends to use this product.

Contents

Reusable laser fiber with SM-905 connector.

INTENDED USE/INDICATIONS FOR USE

The LightTrail Reusable laser fibers are compatible with Boston Scientific (StarMedTec) lasers for indications that are cleared for these laser systems.

This may include but is not limited to surgical procedures involving endoscopic, laparoscopic and open surgical procedures involving vaporization, ablation, fragmentation of calculi, coagulation, hemostasis, excision, resection and incision of soft tissues.

CONTRAINDICATIONS

Refer to the laser User Manual for contraindications related that may be specific to each surgical specialty.

WARNINGS

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these and other instructions relevant to the procedure. Failure to do so may result in complications.

- Improper use of the device or use of a damaged device may result in severe eye or tissue damage, fire in the treatment room and accidental laser exposure to the treatment room personnel or patient.
- Baskets, guidewires and other Ureteroscope accessories may be damaged by direct contact with the laser treatment beam.
- Do not bend fiber at sharp angles. If visible light (aiming beam) can be seen leaking from the fiber, fiber failure may result when therapeutic energy is applied as the fiber is deflected beyond the optical limits of total internal reflection.
- Fiber should not be clamped with forceps or other securing instruments as it could result in fiber damage or breakage.
- Ensure that all procedure room personnel wear appropriate protective eyewear during the delivery of laser energy. Failure to do so may result in injury.
- Residues of blood and tissue and/or damage can lead to overheating of the LightTrail Reusable Laser Fiber tip. This bears a significant risk for the patient and the attending physician. In this case the LightTrail Reusable Laser Fiber tip needs to be cleaned while the laser is OFF; LightTrail Reusable Laser Fiber tip can also be shortened with an adequate tool.

Perforation may occur while executing laser treatment near sensitive areas (blood vessels, intestine). The laser radiation should only be applied for the time to achieve the desired effect of ablation, coagulation, lithotripsy or vaporization.

Gas irrigation especially with endoscopic or interstitial procedures may cause gas embolism.

DO NOT APPLY TO THE CENTRAL NERVOUS SYSTEM AND THE CIRCULATORY SYSTEM

PRECAUTIONS

- All personnel within the treatment room shall wear protective laser eyewear appropriate for Ho:YAG and Tm:YAG lasers. See laser User Manual for requirements concerning protective eyewear.
- The use of irrigation is recommended throughout the procedure to absorb any heat produced.
- Care should be exercised with the glass tip to avoid severe impact or side stresses that may fracture the tip.
- Observe all safety information given in the user manual.
- The general regulations and information for the safe handling of laser radiation must be applied. Apply LightTrail Reusable Laser Fiber only under operating room conditions.
- When using the LightTrail Reusable Laser Fiber in contact surgery treatments avoid forceful movements on tissue surface and lateral pressure to the LightTrail Reusable Laser Fiber tip.
- In the event of no output energy fiber connector may be hot to touch.

ADVERSE EVENTS

The potential adverse effects associated with Ho:YAG and Tm:YAG laser and fibers may include but are not limited to:

- Adhesions
- Discomfort
- Edema
- Embolism
- Extravasation
- Fever
- Fistula
- Healing, Impaired
- Hematoma
- Hematuria
- Hemorrhage
- Hydronephrosis
- Hypertension
- Incontinence
- Infection
- Inflammation
- Irritative voiding symptoms (e.g. dysuria, urgency)
- Laceration
- Nerve injury
- Obstruction
- Pain
- Perforation
- Reflux
- Sexual Dysfunction
- Stenosis
- Sterility/infertility
- Thermal damage
- Vasovagal response

PRODUCT -SPECIFIC SIZE INFORMATION

PRODUCT DESCRIPTION	REFERENCE NUMBER	CORE SIZE	OUTSIDE DIAMETER	WORKING LENGTH	MAXIMUM POWER
LIGHTTRAIL™ REUSABLE 230	M0068F64510	230µm	450µm	3.0m	HOLMIUM: 10W THULIUM: 30W
LIGHTTRAIL REUSABLE 365	M0068F64530	365µm	600µm	3.0m	HOLMIUM: 30W THULIUM: 80W
LIGHTTRAIL REUSABLE 600	M0068F64550	600µm	860µm	3.0m	HOLMIUM: 50W THULIUM: 120W
LIGHTTRAIL REUSABLE 800	M0068F64570	800µm	1180µm	3.0m	HOLMIUM: 50W THULIUM: 120W
LIGHTTRAIL REUSABLE 270	M0068F64630	272µm	420µm	3.0m	HOLMIUM: 22.5W THULIUM: N/A

HOW SUPPLIED

- The LightTrail Reusable Laser Fiber is provided sterilized and must be cleaned, disinfected, and resterilized after each use.
- Do not use if the package is opened or damaged.
- Do not use if labeling is incomplete or illegible.
- Do not use past the expiration date displayed on product label.
- Store in a cool, dry, dark place.

DIRECTIONS FOR USE

Safety information is shown on the labeling of the laser devices and their operating manuals.

Prior to activation of the laser, review the laser User Manual for instructions for proper set-up and operation of the laser system.

Device Preparation

1. Ensure that the sterile packaging is not open, torn or punctured and the product is not damaged. If necessary, return the device to Boston Scientific for replacement.
2. Open the pouch. Slide fiber on shipping card out of pouch into sterile field.
3. Carefully remove the fiber from the shipping card by removing the CONNECTOR end first. Uncoil the remainder of the fiber from the shipping card. Handle the fiber with care as damage may occur if struck or bent sharply. Inspect and treat the output tip with particular care as it represents the most delicate, easily damaged portion of assembly.
4. Ensure that the distal end is intact and that the SMA connector is clean.
CAUTION: Do not use any LightTrail Reusable Laser Fiber with damaged distal end or damaged SMA connector. Avoid touching the exposed connector end.
5. Make sure the laser is in “OFF” mode before inserting the fiber into the SMA port.
6. Attach the connector to the laser and finger-tighten until a secure connection is achieved.
7. Ensure the laser aiming beam is on.
8. Carefully inspect the fiber for kinks, punctures, fractures or other damage. If the fiber appears damaged, do not use the device. Return it to Boston Scientific for replacement.
9. Hold the fiber tip to a non-reflective surface, and ensure that a circular green spot appears. If the spot is weak or not visible, do not use and return the device to Boston Scientific for replacement.
10. Set the laser treatment parameters as instructed in the laser User Manual and within the appropriate power level for the fiber.

INTRA-OPERATIVE INSTRUCTIONS

1. Insert the fiber through the endoscopic delivery system until the fiber tip is in the field of view.
2. Ensure that the fiber tip and aiming beam are in clear view and at a safe distance from the tip of the endoscopic delivery system when laser energy is being applied. Do not use the laser before precise examination and targeting of the tissue to be treated. The LightTrail Reusable Laser Fiber tip must be in the correct position.
3. Position the aiming beam on the target site.
4. Place the laser in “READY” mode.
5. Depress the footswitch to deliver the treatment beam.

POST-OPERATIVE INSTRUCTIONS

The LightTrail™ Laser Fibers are reusable devices. After use, dispose of the device properly or follow the reprocessing steps below.

REPROCESSING (CLEANING, DISINFECTION AND STERILIZATION)

General information

The LightTrail Reusable Laser Fiber is provided sterile.

After application, the LightTrail Reusable Laser Fiber must be disposed or cleaned, disinfected and sterilized according to this procedure.

A reprocessing includes:

- pre-cleaning and preparation
- cleaning and disinfection

as well as

- packaging and sterilization

Note that reprocessing has to be performed in correct order according to this instruction.

Only after preparation and reprocessing - within the limits of maximum 10 application cycles - the LightTrail Reusable Laser Fiber may be used again.

Using the LightTrail Reusable Laser Fiber beyond the allowed 10 application cycles is not permitted and can lead to unforeseeable risks for patient, operator and laser device.

Please note, that devices for reprocessing (washer/disinfector, sterilizer) must be maintained and checked periodically. Please check the particular service manual.

Furthermore any valid legal provisions as well as valid hygienic provisions of the operator (e.g. clinic or hospital) must be kept.

Ensure all personal protective equipment as required by hospital policy is used.

CAUTION: During the reprocessing process do not bend the LightTrail Reusable Laser Fiber below 15 cm. A smaller bending diameter can lead to damage on LightTrail Reusable Laser Fiber and put patients and operators at risk for the next application.

CAUTION: Clean, prepare and check the distal end of LightTrail Reusable Laser Fiber immediately after application.

CAUTION: Put on the plug protection cap before cleaning and disinfection of LightTrail Reusable Laser Fiber. Remove the plug protection cap of LightTrail Reusable Laser Fiber before packaging and sterilization.

CAUTION: Any pathogenic organism polluted medical device can be a source of infections. The use of such medical devices assumes appropriate cleaning, disinfection and sterilization.

CAUTION: To avoid the spreading of pathogenic organisms, the connection plug of LightTrail Reusable Laser Fiber as well as the laser device should not be touched with contaminated gloves or instruments.

CAUTION: Machine thermal disinfection processes are preferred to chemical or chemical-thermal processes.

CAUTION: Chemical disinfection involves the danger of disinfectant residues in openings and on instrument surfaces.

CAUTION: Manual disinfection procedures have lower effects and reproducibility and are only acceptable in case of unavailability of a mechanical procedure.

CAUTION: Avoid extreme temperature differences and rapid temperature changes during cleaning and rinsing cycles.

CAUTION: Do not reuse the LightTrail Reusable Laser Fiber after use with patients with supposed or diagnosed prion disease. Dispose of it in this case immediately after use.

Please note: Defective products must be reprocessed prior to reshipment for complaint check.

Please note: In order to avoid TASS (Toxic Anterior Segment Syndrome) assure that any used alkaline or enzymatic purifier or disinfectant is rinsed and neutralized completely.

Please note: Ensure that during pre-cleaning, cleaning and disinfection there are no affixed residues, as these affect the cleaning, disinfection and sterilization performance.

Please note: Due to the intended application and the type of construction for the LightTrail Reusable Laser Fiber, the classification „critical A“ (=“without special requirements concerning the conditioning”) is advised.

Please note: Boston Scientific is not responsible for deviations from this instruction.

Pre-cleaning:

1. Clean the application area of the fiber several times with a clean, lint-free damp cloth with water or pH neutral enzymatic cleaning detergent to remove visible blood and tissue residue

Note: To avoid the fixation of protein do not use a cleaning agent based on aldehyde, alcohol or peracetic acid.

Note: Water or Cleaning solution should not exceed 43 °C (109 °F)

Preparation:

1. After pre-cleaning, strip fiber to remove the jacket of the LightTrail Reusable Laser Fiber for a length of 40 mm.
2. Carefully score the optical fiber about 5 mm below the end of coating (scratch tool e.g. ceramic knife or diamond scribe) and break the LightTrail Reusable Laser Fiber by pulling in longitudinal direction of the fiber.
3. To check the new tip, point the fiber end at a plain white surface. The spot of the pilot beam should be round with no stray edges.

Rinsing:

1. After fiber preparation step above, remove fiber from laser device and place protective cap on SMA connector (plug).
Note: Protective cap must remain on connector for cleaning and disinfection.
2. Coil the fiber with a diameter of at least 15 cm (6 inches) and secure fiber in the position by looping the distal end through the circle 2-3 times so fiber does not uncoil.
3. Rinse the LightTrail Reusable Laser Fiber with water or pH neutral enzymatic detergent cleaning agent. If necessary, remove any remaining contamination with a lint free and damp cloth. No visible contamination should remain on the LightTrail Reusable Laser Fiber.

Machine cleaning and disinfection (Disinfector, RDG)

General information on the application of cleaning / disinfection devices

Use of disinfection methods other than the described procedures of disinfection of the LightTrail Reusable Laser Fiber are not permitted.

The device has to comply with the requirements of ISO 15883 / AAMI ST15883 and must be validated to these requirements.

The device should have an approved program for thermal disinfection in conformity with EN 15883 / AAMI ST15883 (A0 value > 3000).

The program used must be adequate for the instruments and has to have a sufficient number of rinsing cycles. For rinsing, only use water, that is either sterile, distilled, deionized or filtered by reverse osmosis.

Note: Use filtered air for drying.

General Information on the application of cleaning agent or disinfectant

The cleaning agent must be suitable for the cleaning of instruments made of metal and plastics. Unless thermal disinfection is used, make sure that adequate disinfectant with approved efficiency (e.g. DGHM / VAH or FDA certified or respectively CE approved) is applied and that the disinfectant is compatible with the used cleaning agent.

If an alkaline cleaning agent is used (pH <12.2), an adequate neutralization has to be performed in accordance to the cleaning agent manufacturers' instructions.

Make sure that the cleaning agents and chemicals used are compatible with the instruments.

Comply with the manufacturers' specified concentrations for the cleaning agent and disinfectant solutions used.

Thermal reconditioning

The Vario-TD procedure is validated for the thermal reconditioning of LightTrail™ Reusable Laser Fiber. The process must not be modified:

1. Place a maximum of three LightTrail Reusable Laser Fibers one above the other in a separate instrument basket (grid box) and secure them to minimize movements of the fibers caused by of the water pressure during cleaning and possible damage of the fibers.
2. Position the instrument basket in the cleaning device at an appropriate location.
3. Start the program.

Procedure (parameter):

	Operation	Temperature (°C)	Length (Minimum in minutes)	use of cleaning agent
1	pre-rinsing	5 - 26	5	no
2	cleaning	50 - 55	5	Alkaline detergent (i.e. Neodisher Mediclean Forte)
3	intermediate rinsing	25 - 30	5	no
4	disinfection	90 - 136	5	no
5	drying	15 - 135	5	no

4. When the program is complete, remove the LightTrail Reusable Laser Fiber.
5. Inspect the LightTrail Reusable Laser Fiber in accordance with the Inspection section of this instruction. Repeat the thermal reconditioning process if impurities are still visible on the product.

Please note: If necessary re-dry with a lint-free, sterile cloth. Grooves and cavities in the plug area have to be re-dried with sterile filtered compressed air (< 4 bar).

Chemical reconditioning

A chemical-thermal reconditioning process for the LightTrail Reusable Laser Fiber may be used as an alternative to thermal reconditioning. The procedure must not be modified:

1. Place a maximum of three LightTrail Reusable Laser Fibers one about the other in a separate instrument basket (grid box) and secure them to minimize

movements of the fibers caused by the water pressure during cleaning and possible damage of the fibers.

2. Position the instrument basket in the cleaning device at an appropriate location.
3. Start the program.

Procedure (parameter):

	Operation	Temperature (°C)	Length (Minimum in minutes)	use of cleaning agent
1	pre-rinsing	8 - 12	1	no
2	cleaning	45 - 55	5	pH neutral, enzymatic detergent (i.e. Neodisher Medizym)
3	intermediate rinsing	25 - 30	5	no
4	disinfection	53 - 55	5	Glutaral <10%, Glyoxal <10% (i.e. Neodisher Septo DN)
5	rinsing	8 - 12	2	no
6	drying	15 - 135	5	no

4. Remove the LightTrail Reusable Laser Fiber at the end of the program off the cleaning device.
5. Inspect the LightTrail Reusable Laser Fiber in accordance with the Inspection section of this instruction. Repeat the chemical reconditioning process if impurities are still visible on the product.

Please note: If necessary re-dry by means of a lint-free, sterile cloth. Grooves and cavities in the plug area have to be re-dried with sterile filtered compressed air (< 4 bar).

Manual cleaning and disinfection

Please clean LightTrail Reusable Laser Fiber according to following described procedure.

Cleaning process

1. Immerse the coiled LightTrail Reusable Laser Fiber for at least 5 minutes in room temperature, fresh distilled or deionized water.
2. Remove all blood or cellular residue with a soft brush or a sterile, lint-free cloth, paying special attention to the cavities, threads and holes in the connector (plug) area.
3. Remove the LightTrail Reusable Laser Fiber from the bath.
4. Flush the connector (plug) a minimum of 5 times with a least 10 mL pH neutral, enzymatic solution (i.e. Neodisher Medizym) using a water pistol, syringe or equivalent.
5. Immerse the entire LightTrail Reusable Laser Fiber in a bath of pH neutral, enzymatic solution for at least 5 minutes.
6. Flush the connector (plug) again a minimum of 5 times with a least 10mL pH neutral, enzymatic solution (i.e. Neodisher Medizym).
7. Remove the LightTrail Reusable Laser Fiber from the bath.
8. Rinse the fiber thoroughly at least 2 times, at least 15 seconds each, with sterile distilled or deionized water.

- Rinse the connector at least 1 more time, at least 15 seconds, with sterile distilled or deionized water with the connector (plug) downward so water can run off.
- Check visually if the cleaning of the LightTrail™ Reusable Fiber was successful. Repeat cleaning process steps 1-9 again if necessary.

Disinfection process

- Flush the connector (plug) a minimum of 5 times with a least 10 mL sterile distilled or deionized water.
- Immerse the LightTrail Reusable Laser Fiber in solution of <1% of Ortho-Phthalaldehyde (i.e. Cidex OPA) disinfectant for at least 12 minutes.
- Flush connector (plug) again 5 times with at least 10 mL solution of <1% of Ortho-Phthalaldehyde (i.e. Cidex OPA) disinfectant.
- Remove LightTrail Reusable Laser Fiber from disinfection solution.
- Rinse the fiber thoroughly at least 2 times, 15 seconds each, with sterile distilled or deionized water.
- Rinse the connector at least 1 more time, at least 15 seconds, with sterile distilled or deionized water with the connector (plug) downward so water can run off.
- Check visually if the LightTrail Reusable Laser Fiber is free of soil and contamination. Repeat Cleaning and Disinfection procedures if impurities are still present.
- Dry the LightTrail Reusable Laser Fiber with sterile, filtered compressed air with pressure <4 bar (58 psi).
- Inspect the LightTrail Reusable Laser Fiber per Inspection section of this instruction.

Inspection

Macroscopic visual inspection


- Check the LightTrail Reusable Laser Fiber after cleaning and disinfection for corrosion, damaged surfaces, splintered materials, and contamination.
- Repeat cleaning and disinfection if viewable contamination is present.
- Dispose the LightTrail Reusable Laser Fiber if damage is present.

Microscopic inspection

- Remove the plug protective cap.
- Check the proximal fiber end of the plug with a microscope (fiberscope) with sufficient resolution.
- The proximal fiber end may show excessive scratches, residues, or respectively discolorations due to cleaning agents or disinfectants.
- Dispose the LightTrail Reusable Laser Fiber in case of damage.

Packaging

- If a scale like the plastic scale shown in this section is available, mark the performed reconditioning cycle by cutting off the next cycle number. For that purpose, use a sterile scissor.
- Any use of LightTrail Reusable Laser Fiber is counted electronically and shown on the laser device prepared for this; see User Manual of particular laser device.
- Ensure that the plug protective cap is off.
- Put the coiled LightTrail Reusable Laser Fiber into a suitable disposable sterilization package (single- or double wrapping) that conforms with the following requirements:

- DIN EN ISO / ANSI AAMI ISO 11607-1/2
- Approved for steam sterilization (temperature resist up to at least 135 °C (275 °F) sufficient steam transmission).
- Pay attention to sufficient protection from mechanical damages. Double packaging is preferred to single packaging.

Sterilization

Ensure the protective cap is removed prior to sterilization. Sterilize the LightTrail Reusable Laser Fiber using the following procedures:

Steam sterilization

Steam sterilization is subject to following criteria:

- Steam sterilizer according to DIN EN 13060-2004, DIN EN 285, AAMI / ANSI ST55
- Validated according to ISO / ANSI AAMI ISO 17665 (valid commissioning and product-specific performance evaluation)
- Maximal sterilization temperature 135 °C (275 °F); (plus tolerance according to EN ISO / ANSI AAMI ISO 17665)
- Sterilization time (exposure time at sterilization temperature)
 - Gravity Cycle: 30 minutes at 121 °C (250 °F) with 15 minute dry time
 - Pre-vacuum Cycle: 4 minutes at 132 °C (270 °F) with 20 minute dry time
 - Pre-vacuum Cycle*: 3 minutes at 134 °C (273 °F) with 30 minute dry time
 - Pre-vacuum Cycle: 3 minutes at 135 °C (275 °F) with 16 minute dry time
 - Pre-vacuum Cycle*: 18 minutes at 134 °C (273 °F) with 30 minute dry time

*This cycle is not recognized by the FDA for use in the United States.

Storage

After sterilization, place the LightTrail Reusable Laser Fiber in the sterilization package and store in a cool, dry, dark place.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2019 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd