EndoVive™ One Step Button™

INITIAL PLACEMENT PUSH AND PULL

Low Profile Percutaneous Endoscopic
Gastrostomy Kit

50652590-01

2020-04

Directions for Use	2
Instrucciones de uso	7
Mode d'emploi	12
Gebrauchsanweisung	17
Istruzioni per l'uso	22
Gebruiksaanwijzing	27
Instruções de Utilização	32

EndoVive™ One Step Button™

INITIAL PLACEMENT PUSH AND PULL

Low Profile Percutaneous Endoscopic Gastrostomy Kit

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The EndoVive One Step Button Device is a single-use, sterile, disposable device used during initial placement for direct feeding. The Kit contains the components needed to perform the procedure. The One Step Button Device is to be placed with a guidewire.

The catheter is Magnetic Resonance (MR) Safe.

Contents

- (1) Guidewire or Insertion Wire
- (1) 5 g Lubricating Jelly (if applicable)
- (1) One Step Button Device
- (1) Packet Antiseptic Ointment (if applicable)
- (4) Drainage Sponges and/or Gauze Pads
- (1) 1 mm Spacer Disc
- (1) 2 mm Spacer Disc

- (1) 3 mm Spacer Disc
- (1) 4 mm Spacer Disc
- (1) 5 mm Spacer Disc
- (1) Right Angle Feeding Set
- (1) Bolus Feeding set
- (1) 60 cc Syringe
- (1) Button Decompression Tube
- (1) Three Pack Iodine Swabs (if applicable)
- (1) Fenestrated Drape
- (1) 5.0 ml Syringe
- (1) 19 ga (1.1 mm) x 1.5 in (38.1 mm) Filter Needle (if applicable)
- (1) 25 ga (0.5 mm) x 1 in (25.4 mm) injection Needle (if applicable)
- (1) Surgical Scalpel
- (1) Percutaneous Stoma Measuring Device
- (1) Bumper
- (4) Gauze Pads
- (1) 5 g Lubricating Jelly (if applicable)
- (1) Retrieval Snare

INTENDED USE/INDICATIONS FOR USE

The One Step Button Device is intended to provide nutrition to a patient directly into the stomach through a stoma. It is indicated for use on patients who are unable to consume nutrition by conventional means.

CONTRAINDICATIONS

The One Step Button Device is not indicated for the following:

- Children under 10 Kilos/22 pounds.
- The 24F is not indicated in children under 15 kilos/33 pounds.
- Not recommended in individuals with abnormally high intra-abdominal pressures.
- Obstruction of the esophagus which may prevent the introduction or removal of the feeding tube.
- Inability to identify transabdominal illumination or needle placement.
- Multiple surgical procedures near the gastrostomy site.
- Patients on anticoagulant drugs.
- High medical risk patients.

PRECAUTIONS

Do not use this product without reading and understanding the complete instructions enclosed herein.

ADVERSE EVENTS

The following complications may occur with the use of direct feeding devices: fever, gastric distention, infection, blockage/occlusion, tissue necrosis, migration, peritonitis, sepsis, erosion/embedding into the gastric wall ("Buried Bumper Syndrome"), aspiration, bleeding, fistula, gastroparesis, GE reflux, pain, perforation,

ulceration, tube clogging, malposition, leakage, kinking, inadvertent removal, small bowel obstruction, granulation tissue, and pneumoperitoneum.

OPERATIONAL INSTRUCTIONS

Prior to Use

Caution: Carefully examine the unit to verify that the contents have not been damaged in shipment. DO NOT USE if damaged. Immediately return damaged product to Boston Scientific Corporation.

Note: In order to select a One Step Button™ Device of the appropriate size, use the Percutaneous Stoma Measuring Device (P.S.M.D.) included with this kit. See P.S.M.D. Instructions For Use.

Caution: This device is intended to connect to enteral feeding devices. Care should be taken to avoid the potential harm of misconnection with non-enteral feeding connectors.

Directions

Push Method

1. Placement of Guidewire

- Pass the guidewire down the already placed P.S.M.D. cannula. Upon visualization of the guidewire, grasp the wire with the snare.
- Withdraw the endoscope, snare, and the attached snared wire.
- Unlock the P.S.M.D. and completely compress the Thumb Ring until it lies flat against the top of the spool's barrel.

Warning: DO NOT REMOVE MEASURING DEVICE IN THE 'BLOOM' POSITION. Failure to compress the P.S.M.D. to eliminate the extended "BLOOM" could cause damage to the gastric and abdominal walls when removed.

2. One Step Button Catheter Preparation

- Remove One Step Button Catheter from the package and apply water soluble lubricant to both ends of the catheter.
- Keep the P.S.M.D. cannula in the abdomen, perpendicular to the abdominal wall.

3. Advancement of the One Step Button Catheter over the Guidewire

- Slide the One Step Button Catheter over the guidewire exiting the patient's mouth (Figure 1).


Figure 1

- FIRM TENSION IS REQUIRED AT BOTH ENDS OF THE GUIDEWIRE to allow for passage of the One Step Button Catheter into the alimentary tract and into the stomach.
- Slowly advance the One Step Button Device over the guidewire until the P.S.M.D. cannula is displaced by the catheter.
- Remove and discard the P.S.M.D.
- Gently pull the catheter through the abdominal wall. Once the second dilating portion of the catheter begins to exit the abdominal wall, caution must be taken to advance the remaining portion of the catheter slowly.
- When the red strip is visible, advance the catheter very slowly. With your other hand, provide countertraction on the abdomen around the stoma site. Advance the catheter in small increments until the folded encased wings of the One Step Button Device have exited the abdominal wall.
- After the One Step Button Catheter is pulled through the abdominal wall, the endoscopist should re-scope the patient to view the internal dome placement.
- Extending out of the stomach is the dilating tube tip and encased Button flange.
- Keep the catheter perpendicular to the stomach with care not to allow the encased Button flange to retract back into the abdominal wall.

4. Removal of Sheath

- Grasp the loose end of the tear strip (black suture material located at the distal end of the red strip), and pull up along red strip (away from patient) (Figure 2).
- Peel away sheath and release external Button flange.
- Discard sheath, tubing and dilating tip.
- Remove the guidewire either through the mouth or through the Button catheter.
- Discard guidewire.


Figure 2

Pull Method

1. Placement of Insertion Wire

- Pass the insertion wire down the already placed P.S.M.D. cannula. Upon visualization of the insertion wire, grasp the wire with the snare.
- Withdraw the endoscope, snare, and the attached snared wire.
- Unlock the P.S.M.D. and completely compress the Thumb Ring until it lies flat against the top of the spoon's barrel.

Warning: DO NOT REMOVE MEASURING DEVICE IN THE 'BLOOM' POSITION. Failure to compress the P.S.M.D. to eliminate the extended 'BLOOM' could cause damage to the gastric and abdominal walls when removed.

2. Attachment of the One Step Button™ Catheter to the Insertion Wire

- Pass the insertion wire through the One Step Button Device wire loop (Figure 3).
- Pass the Button end of the One Step catheter through the insertion wire loop and pull the entire catheter through (Figure 4).
- Grasp the dilating tip between the thumb and forefinger, tighten the loops together to form attachment (Figure 5).
- Apply water soluble lubricant to both ends of the catheter.
- Keep the locked P.S.M.D. in the abdomen, perpendicular to the abdominal wall.


Figure 3


Figure 4


Figure 5

3. One Step Button Device Catheter Placement


Figure 6

- Slowly pull the catheter through the alimentary tract by applying tension on the insertion wire at the stoma site. Continue to pull the catheter into the stomach. The P.S.M.D. will be displaced by the catheter and should be removed.
- Gently pull the catheter through the abdominal wall. Once the second dilating portion (proximal end) of the catheter begins to exit the abdominal wall, caution must be taken to advance the remaining portion of the catheter slowly.
- When the red strip is visible, advance the catheter very slowly. With your hand, provide countertraction on the abdomen around the stoma site. Advance the catheter in small increments until the folded encased wings of the One Step Button Device have exited the abdominal wall.
- After the One Step Button Device Catheter is pulled through the abdominal wall, the endoscopist should re-scope the patient to view the internal dome placement.
- Extending out of the stomach is the dilating tube tip and encased Button flange.
- Keep the catheter perpendicular to the stomach with care not to allow the encased Button flange to retract back into the abdominal wall.

4. Removal of Sheath

- Grasp the loose end of the tear strip (black suture material located at the distal end of the red strip), and pull up along red strip (away from patient) (Figure 7).
- Peel away sheath and release external Button flange.
- Discard sheath, tubing and dilating tip.


Figure 7

5. One Step Button™ Device Placement Confirmation

- With the One Step Button Device in place:
 - Irrigate with 5-10 mL of sterile saline.
 - ROTATE the Button to ensure comfortable, stable placement. The Button should rotate freely.

WARNING: DO NOT INJECT AIR INTO THE BUTTON TO DETERMINE DEVICE LOCATION, THIS COULD LEAD TO GASTRIC DISTENSION.

WARNING: IT IS RECOMMENDED THAT 24 HOURS ELAPSE BEFORE THE FIRST FEEDING, AS GASTRIC DISTENSION MAY OCCUR.

6. Use of Button Spacers (Optional)

Indications: If the Button shaft is too long, spacers can be used between the stoma and external flange of the Button. The purpose of the spacers is to provide the optimal Button length and fit. Spacers are available in 1 mm, 2 mm, 3 mm, 4 mm, and 5 mm discs. If the stoma tract is between sizes, always choose the larger size Button and utilize the spacers.

- Spacers can stay in place until a stoma tract has been developed.
- Selection of spacers is based on desired tension required to compress fascia layers during the formation of the stoma.

Placement of Spacers: The Button must be in place prior to placement of the spacer. Carefully stretch spacer disc over the longer side of the Button's external flange, then through the shorter side so that it is seated in between the Button's external flange and the skin's surface. The spacer disc should fit loosely and should not prevent the Button from being easily rotated.

WARNING: SPACERS SHOULD BE USED ONLY AFTER THE STOMACH HAS BEEN DECOMPRESSED AND THE BUTTON CHECKED FOR EASE OF ROTATION. TOO TIGHT A FIT CAN CAUSE PRESSURE NECROSIS.

- To confirm proper fit, the Button should rotate with ease.

Removal of Spacers: Carefully cut the appropriate spacer disc, being careful NOT TO CUT the Button. Gently pull away from the button and discard.

7. Tube Care

- Once the feeding has been completed, flush the device with 5-10 mL (3-5 mL for children) of water to rinse any remaining particulate.
- Cleanse stoma site and external portion of gastrostomy tube with soap and water. The bolster should be rotated daily for site hygiene. The stoma site should be clean and dry at all times.
- Monitor the tube for possible migration. Tube migration could result in the following: inability to feed, peritonitis, infection and associated sequelae.
- Prior to discharge, instruct the patient on device cleaning, monitoring, and medication administration (when applicable). Ensure that the patient is provided with a Patient Care guide for additional device monitoring instructions.

8. Removal of Button

- Grasp the tabs and, using a constant traction, carefully pull the button out of the stoma. NOTE: If desired, a viscous topical anesthetic gel may be applied to the skin surrounding the stoma and into the stoma itself. Wipe the skin dry if gastric content of enteral feeding product spills.
- Removal can also be done by placing an obturator rod through the catheter to distend/elongate the Button dome for ease of removal. The Button can also be removed endoscopically if desired. Insert the endoscope, cut the shaft and retrieve the dome.

Warning: DO NOT let the dome pass through the gastrointestinal tract. Obstruction can result in associated sequelae.

- Following removal and before spontaneous closure of the stoma, insert a new Button into the stoma if enteral feeding by this route is still desired.
- If spontaneous closure of stoma is desired, apply a dressing over the stoma site.

Warning: Removal should not occur until stoma is fully established which may take two months or more, otherwise loss of patency of stoma tract could occur. Removal is recommended to be performed by a qualified clinician.

Note: The silicone feeding tubes should remain in place until the stoma is mature. The healthcare provider should decide when a replacement device is necessary, as precise device life cannot be predicted due to several factors. These factors may include medications, nutrition, gastric pH and tube care. Clinical literature suggests that silicone PEG implant replacement be evaluated approximately one year after implantation.

HOW SUPPLIED

The EndoVive™ One Step Button Device, Initial Placement Gastrostomy Kit is separately packaged as one (1) unit, sterilized by ethylene oxide gas. The unit is included as a component in the EndoVive One Step Button Device Tray. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the sterilization expiration date on package label.

Kits that contain antibiotics or antiseptics must be stored between 15-30°C

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
© 2020 Boston Scientific Corporation or its affiliates.