Extractor™ Pro DL

Retrieval Balloon Catheter

Directions for Use	2
Instrucciones de uso	6
Mode d'emploi	10
Gebrauchsanweisung	14
Istruzioni per l'uso	19
Gebruiksaanwijzing	24
Instruções de Utilização	29

90965436-01

2014-11

TABLE OF CONTENTS

WARNING 3

DEVICE DESCRIPTION..... 3

 Catheter Shaft 3

 Balloon/Radiopaque Marker 3

 Contrast Injection..... 3

 Contents..... 3

INTENDED USE/INDICATIONS FOR USE..... 3

CONTRAINDICATIONS..... 3

WARNINGS..... 4

PRECAUTIONS..... 4

ADVERSE EVENTS 4

HOW SUPPLIED..... 4

 Handling and Storage..... 4

DEVICE PREPARATION..... 4

DIRECTIONS FOR USE..... 5

 Stone Retrieval 5

 Duct Occlusion 5

DEVICE REMOVAL..... 5

WARRANTY 5

Extractor™ Pro DL

Retrieval Balloon Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Extractor Pro DL Retrieval Balloon Catheter is a double lumen retrieval balloon catheter. The catheter may be placed with or without the aid of a guidewire. The Extractor Pro DL Retrieval Balloon Catheter is capable of accepting a 0.035 in (0.89 mm) guidewire in the guidewire lumen. The balloon lumen is used to inflate the balloon using the syringe packaged with the device. The different balloon sizes are marked on the syringe for your convenience.

Catheter Shaft

The Extractor Pro DL Retrieval Balloon Catheter shaft tapers from 7F (2.3 mm) to 6F (2.0 mm) at the distal 7.5 cm portion. The balloon lumen of the catheter allows for inflation of the balloon only. The guidewire channel of the catheter accommodates a 0.035 in (0.89 mm) guidewire through the length of the Extractor Pro DL Retrieval Balloon Catheter.

Balloon/Radiopaque Marker

A latex balloon is mounted at the distal tip of the catheter with a radiopaque marker located on the proximal end of the balloon. A pre-staked syringe is enclosed in the package. The syringe should be used for inflating the Extractor Pro DL Retrieval Balloon Catheter. The syringe is marked with 2 numeric markers, which indicate balloon diameters.

Contrast Injection

Contrast media injection may be achieved through the guidewire lumen by removing the guidewire and attaching a contrast filled syringe to the proximal luer.

Warning: DO NOT EXCEED THE LARGER O.D (OUTSIDE DIAMETER) OF THE BALLOON INDICATED ON THE ENCLOSED PRE-STAKED SYRINGE WHEN INFLATED.

Contents

- Balloon Catheter
- Stopcock
- Syringe

INTENDED USE/INDICATIONS FOR USE

The Extractor Pro DL Retrieval Balloon Catheter is used endoscopically to 1) remove stones from the biliary system; or 2) to facilitate injection of contrast medium while occluding the duct with the balloon.

CONTRAINDICATIONS

Contraindications for this device are those specific to endoscopic retrograde cholangiopancreatography (ERCP) and endoscopic sphincterotomy (ES).

WARNINGS

- Check for proper position of the balloon catheter using contrast medium injection and fluoroscopy. Balloon inflation in an improper location may lead to patient injury.
 - If resistance is met during procedure, do not advance the balloon without first determining the cause or resistance and taking remedial action.
-

PRECAUTIONS

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

- The Extractor™ Pro DL Retrieval Balloon Catheter family of catheters should only be used by or under the supervision of physicians thoroughly trained in endoscopic retrograde cholangiopancreatography.
- Monitor balloon inflation under fluoroscopy.
- The Extractor Pro DL Retrieval Balloon Catheter family of catheters are designed for use with endoscopes which have a working channel of at least 3.2 mm (or larger).
- Due to various guidewire tolerances among manufacturers, it is recommended the Extractor Pro DL Retrieval Balloon Catheter family of catheters be used with a Boston Scientific 0.035 in (0.89 mm) guidewire.
- The Extractor Pro DL Retrieval Balloon Catheter family of catheters are supplied sterile. If the catheter package is opened or damaged, do not use catheter.

ADVERSE EVENTS

Adverse events include, but may not be limited to:

- pancreatitis
- perforation
- hemorrhage
- hematoma
- cholangitis
- stone impaction
- septicemia/infection
- allergic reaction to contrast medium

HOW SUPPLIED

The Extractor Pro DL Retrieval Balloon Catheter is supplied sterile and intended for single use only. Sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that catheters are used prior to the sterilization expiration date on package label. Store in a cool, dry, dark place.

DEVICE PREPARATION

Precaution: A thorough understanding of the technical principles, clinical applications, and risks associated with endoscopic sphincterotomy (ES) and endoscopic retrograde cholangiopancreatography (ERCP) is necessary before using this product.

Precaution: The Extractor Pro DL Retrieval Balloon Catheter is supplied in a sealed sterile pouch. Prior to use, inspect the pouch for any breach to ensure a sterile device.

1. Prior to use, carefully examine the unit to verify that the sterile package or product has not been damaged in shipment.
 2. Remove the Extractor Pro DL Retrieval Balloon Catheter from its package and inspect for any damage.
 3. Prior to clinical use, test the catheter by flushing it with sterile water or saline.
-

Caution: Kinks in the catheter will hinder injection capability. Do not use the Extractor Pro DL Retrieval Balloon Catheter if any defect is found during inspection.

4. Use the pre-staked syringe enclosed in the package to inflate the balloon. The syringe is marked with 2 numeric markers, which indicate balloon diameters. Prior to connecting the syringe to the catheter, pull

the plunger back until it stops. Prior to use, inflate the balloon to the maximum balloon diameter to ensure there are no leaks. Do not use the catheter if leaks are detected. To completely deflate the balloon, pull back the plunger until it stops. The catheter is now ready for use.

5. If contrast medium will be injected through the guidewire lumen, use a 30 ml (30 cc) syringe to inject the contrast or saline through the lumen to expel any air bubbles before inserting the catheter through the endoscope.

Precaution: The Boston Scientific RX biopsy cap is not recommended to be used with the Extractor™ Pro DL 15-18 mm Retrieval Balloon Catheter.

DIRECTIONS FOR USE

With the endoscope in the proper position, insert the Extractor Pro DL Retrieval Balloon Catheter independently or over a guidewire and through the instrument channel to the appropriate location. A radiopaque marker located just proximal to the balloon will aid in positioning the catheter while using fluoroscopy.

Caution: The Extractor Pro DL Retrieval Balloon Catheter should be advanced through the endoscope using short, deliberate 2-3 cm movements to prevent inadvertent damage to the device.

Stone Retrieval

1. Position the catheter past the stone(s) being retrieved from the appropriate duct.
2. Check for proper position of the balloon catheter using contrast medium injection and fluoroscopy. When the balloon is in the desired position, inflate the balloon with the enclosed syringe.
3. The plunger should line up with the corresponding marker on the syringe in order to inflate to the desired balloon diameter.

Note: For the 9 mm size of the 9-12 mm balloon, inflate to the 12 mm marker first, then to the 9 mm marker on the syringe.

4. With the balloon inflated, withdraw the catheter from the duct to retrieve the stone(s).

Duct Occlusion

1. To facilitate the injection of contrast medium, position the catheter in the appropriate location to occlude the duct.
2. Check for proper position of the balloon catheter using contrast medium injection and fluoroscopy. Balloon inflation in an improper location may lead to patient injury.
3. When the catheter is in the desired position, inflate the balloon with the syringe. The syringe is marked with 2 numeric markers, which indicate balloon diameters.
4. Inject the contrast medium using a 30 ml (30 cc) syringe.

DEVICE REMOVAL

Completely deflate the Extractor Pro DL Retrieval Balloon Catheter prior to withdrawing the device into the scope.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2014 Boston Scientific Corporation or its affiliates.