Obtryx™ II System

CURVED

Transobturator Sling with
PrecisionBlue™ Design

Directions for Use	2
Instrucciones de uso	9
Mode d'emploi	15
Gebrauchsanweisung	21
Istruzioni per l'uso	27
Gebruiksaanwijzing	33
Instruções de Utilização	39

50839563-01

2019-07

TABLE OF CONTENTS

WARNING 3

DEVICE DESCRIPTION..... 3

INDICATIONS FOR USE..... 3

CONTRAINDICATIONS 4

HOW SUPPLIED..... 4

DIRECTIONS FOR USE..... 4

TENSION MESH/SLEEVE REMOVAL..... 5

GENERAL WARNING..... 6

POST PROCEDURAL WARNING 6

ADVERSE EVENTS 6

PRECAUTIONS..... 7

WARRANTY 8

Obtryx™ II System

CURVED

Transobturator Sling with PrecisionBlue™ Design

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician trained in use of surgical mesh for repair of stress urinary incontinence.

WARNING

Contents supplied STERILE using ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

WARNING

This product is intended for use only by clinicians with adequate training and experience in the surgical treatment of stress urinary incontinence (SUI). The physician is advised to consult the medical literature regarding techniques, complications and hazards associated with the intended procedures.

DEVICE DESCRIPTION

The Obtryx II System is a sterile, single use system consisting of two (2) delivery devices and one (1) mesh assembly. The mesh assembly is comprised of a polypropylene knitted mesh with dilator legs and a center tab. At the distal ends of the dilator legs there are association loops designed to be placed in the needle slot of the distal end of the delivery device. The disposable delivery device consists of a handle with a stainless steel needle. The needle is designed to facilitate the passage of the mesh assembly through bodily tissues for placement through the obturator foramen.

The mid-urethral slings are composed of materials that are electrically nonconductive, nonmetallic, and nonmagnetic, and pose no known hazards resulting from exposure to any MR (Magnetic Resonance) environment. Therefore, the suburethral slings are considered MR Safe.

INDICATIONS FOR USE

The mesh implant is intended for use as a suburethral sling for the treatment of stress urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

The Obtryx II (Curved) Delivery Device is intended for use as an aid in insertion, placement, fixation, and anchoring of Obtryx II surgical mesh during urogynaecological procedures.

CONTRAINDICATIONS

The mesh suburethral sling implant is contraindicated in the following patients:

- Pregnant patients, patients with potential for future growth or patients that are considering future pregnancies.
- Any patients with soft tissue pathology into which the implant is to be placed.
- Patients with any pathology which would compromise implant placement.
- Patients with any pathology, such as blood supply limitations or infections that would compromise healing.

HOW SUPPLIED

The device is supplied sterile. Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store at a controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the expiration date on package label.

DIRECTIONS FOR USE

Prior to Use

Carefully examine the system to verify that neither the contents nor the sterilized package has been damaged in shipment. **DO NOT USE** if sterile barrier on product is damaged. Immediately return damaged product to Boston Scientific.

The design of the Obtryx™ II System allows the operator a percutaneous approach utilizing a transobturator technique. See Figure 1 for parts description.


Figure 1: Parts Description

Prepare and drape the patient using standard surgical practice.

WARNING

Assure that the bladder is empty prior to initiating the use of this product. Ensure that the bladder, urethra and other important landmarks are properly identified.

Steps to Use

1. Prepare the skin lateral to the inferior pubic ramus and vaginal operative sites.
2. Make a 1.0 cm to 1.5 cm vertical midline incision on the anterior vaginal wall at the level of the mid-urethra. Dissect bilaterally to the interior portion of the inferior pubic ramus at the 45° angle off the midline creating a pathway for delivery device placement.
3. Create a vertical skin incision large enough to insert tip of needle just lateral to the edge of the inferior pubic ramus at the junction where the inferior pubic ramus and the adductor longus muscle meet. Repeat on the contralateral side.

WARNING

If excessive force is encountered during advancement/withdrawal, stop and determine remedial action prior to proceeding.

4. Grasp the device handle and insert one (1) needle through one (1) skin incision, piercing through the obturator muscle and obturator membrane. Turn the handle at the 45° angle medial towards the midline. Place the opposite hand's forefinger into the lateral dissection of the vaginal incision, placing the fingertip on the distal end of the needle. Guide the distal end of the needle around the inferior pubic ramus through the vaginal incision, maintaining contact with the finger.

WARNING

Pay careful attention to avoid the adductor longus tendon with the delivery device.

WARNING

Make sure the delivery device and mesh assembly pass sufficiently lateral to the urethra in order to avoid urethral injury.

5. Engage one (1) association loop to the distal end of the needle (see Figure 2) protruding through the vagina.


Figure 2: Association Loop Engagement

6. Pull the needle out through the skin incision. Be sure that the mesh assembly is not twisted and lies flat under the urethra, with the blue center tab positioned suburethrally, facing outward.
7. Remove the association loop from the needle (see Figure 3).


Figure 3: Association Loop Removal

8. Repeat Steps 4-7 on the contralateral side with the second needle.
9. Cystoscopy may be performed at this time, to be determined at the physician's discretion.
10. Next see section "Tension Mesh/Sleeve Removal."

TENSION MESH/SLEEVE REMOVAL

1. Adjust the mesh/sleeve by pulling outwards on the dilators so that the blue center tab is centered below the urethra.
2. Appropriately tension the mesh/sleeve according to physician preference.
3. Once proper tension is achieved, cut the leader loop that is on the outside of the sleeve that is connecting the dilator leg and sleeve to the mesh. Pull outward on the dilator to remove the sleeve leaving the mesh in place. Repeat on the other side. (See Figure 4).


Figure 4: Tension Mesh/Sleeve Removal

4. Grasp the blue center tab and cut the center tab lead located on the side of the center tab to release the tab from the mesh. Remove the center tab and center tab lead from the vaginal canal.
5. Gently pushing downward on the skin incisions, cut the distal ends of the mesh and confirm that those ends retract into the skin incisions.
6. Close all incisions per standard practice.

GENERAL WARNING

The risks and benefits of performing a suburethral sling procedure in the following should be carefully considered:

- Careful consideration should be given to performing this procedure for patients with untreated coagulopathies or who are being treated with either anticoagulants or antiplatelet agents.
- Patients with hypertonic bladders or vesico ureteral reflux.
- Take special care in cases of bladder prolapse because of anatomical distortion. If the patient requires a cystocele repair, it should be done prior to the suburethral sling procedure.
- Vaginal and urinary tract infection should be treated prior to a suburethral sling implantation procedure.
- User should be familiar with surgical procedures and techniques involving nonabsorbable meshes.
- Good surgical practices should be followed for management of contamination or infected wounds.
- Mesh is considered a permanent implant. Removal of mesh or correction of mesh related complications may involve multiple surgeries.
- Complete removal of mesh may not be possible and additional surgeries may not always fully correct the complications.

POST PROCEDURAL WARNING

- If subsequent infection occurs, follow appropriate medical intervention practices.
- The patient should be advised that future pregnancies may negate the effects of this procedure and the patients may again become incontinent.

ADVERSE EVENTS

The following adverse events have been reported due to suburethral sling placement, but are not limited to:

- As with all implants, local irritation at the wound site and/or a foreign body response may occur.
- Tissue responses to the mesh implant could include:
 - erosion/exposure/extrusion of the mesh through the vaginal or urethral mucosa, bladder wall, or other surrounding tissue, which may result in stone formation
 - scarring/scar contracture
 - device migration
 - fistula formation and inflammation

The occurrences of these events may require surgical intervention and possible removal of the entire mesh.

- Like all foreign bodies, the mesh may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.
- Allergic reaction has been reported.
- Known risks of surgical procedures for the treatment of incontinence include:
 - pain, ongoing pain (pelvic, vaginal, groin/thigh, dyspareunia)
 - infection
 - detrusor instability
 - complete failure of the procedure
 - voiding dysfunction (incontinence, mild to moderate incontinence due to incomplete urethral support or due to overactive bladder)
 - bruising, bleeding (vaginal, hematoma formation)
 - abscess
 - vaginal discharge
 - dehiscence of vaginal incision
 - edema and erythema at the wound site
 - perforation or laceration of vessels, nerves, bladder, urethra, bowel or other tissues may occur during placement.

The occurrence of these events may require surgical intervention. In some instances the response to these events may persist as a permanent condition after the intervention.

PRECAUTIONS

- The use of polypropylene mesh in urogynecologic procedures such as the treatment of stress urinary incontinence regardless of the route of delivery (transvaginal, suprapubic or transobturator), has been associated with cases of erosion. Erosion has been reported in bladder, vagina, urethra, ureter, and bowel. Treatment of the erosion may require surgical removal.
- As with all surgical procedures, certain risk factors are known to impact patient outcomes in the pelvic floor which include, but are not limited to, impaired vascularity (e.g. diabetes, smoking status, estrogen status, pelvic floor radiation exposure, etc.), age, pelvic floor myalgia, impaired wound healing (e.g. diabetes, steroid usage, etc.), or active infection in or near the surgical site. The above pathophysiologic conditions must be considered when determining whether the patient is an appropriate candidate for mesh implantation, either by transvaginal, suprapubic or transobturator route.
- Standard surgical practices should be followed for the suburethral sling procedure as well as for the management of contaminated or infected wounds.
- Bleeding can occur. Check carefully before releasing patient from the hospital.
- Do not remove the protective plastic sleeve covering mesh implant until proper position has been confirmed.
- Ensure the mesh is placed without tension under the mid-urethra.
- Use of this device should be done with the understanding that subsequent infection may require removal of the mesh.
- Physician should determine when it is suitable for each patient to return to normal activities.
- Patients should be counseled when to resume vigorous activities (heavy lifting, exercise), and intercourse after the procedure.

- Should dysuria, bleeding or other problems occur, the patient should be instructed to immediately contact the physician.
- Do not use any mechanical means of contact with the mesh (such as clips, staples etc.) within the urethral support region of the mesh as mechanical damage to the mesh may occur.
- Avoid excessive tension on the mesh during handling.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.**
**Boston Scientific Limited
**Boston Scientific (Australia) Pty Ltd
**© 2019 Boston Scientific Corporation or its affiliates.