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zh-TW

This document is valid only on the date printed. If unsure of the print date, please re-print to ensure use of the latest revision of the IFU (available at www.e-ifu.com). The onus resides with the user to ensure that the most up-to-date IFU is used.

en

ENVOY™ Distal Access Guiding Catheter

IMPORTANT INFORMATION

Please read before use

Some products may not be licensed in all jurisdictions.

STERILE EO

 Rx Only

MD

Description

The ENVOY™ Distal Access (DA) Guiding Catheter is a variable stiffness, braided catheter with a large non-tapered lumen that facilitates the intravascular passage of interventional devices. The distal segment is flexible for navigation into distal vasculature. The catheter has an outer hydrophilic coating that reduces friction during manipulation in the vessel. The lubricious PTFE lined inner lumen is designed to facilitate movement of the guide wires and other devices. A luer fitting located on the end of the catheter hub can be used to attach accessories. The distal section of the catheter is radiopaque to aid visualization under fluoroscopy, and the distal tip is clearly distinguished by a radiopaque marker band which is approximately 3 mm from the distal tip. The catheter is available with preshaped tips to facilitate positioning. A peel away introducer is included to facilitate insertion into the sheath.

Intended Use

The Envoy Distal Access Guiding Catheter is intended as a conduit device for the intravascular introduction of interventional/diagnostic devices to the neuro vasculature.

Indications for Use

The ENVOY Distal Access Guiding Catheter is indicated for use in the neuro vasculature for the intravascular introduction of interventional/diagnostic devices to facilitate neurovascular procedures.

Patient Population

The ENVOY Distal Access Guiding Catheter is intended for use in patients undergoing diagnostic or interventional procedures that require neuro intravascular access for delivery of diagnostic or interventional devices.

Intended User

The ENVOY Distal Access Guiding Catheter should be used only by physicians trained in interventional percutaneous procedures at medical facilities with the appropriate imaging equipment.

Clinical Benefits (Applicable only in markets subject to EU Regulation 2017/745)

As an accessory device comprising part of the endovascular platform for neurovascular index procedures, the clinical benefit of the ENVOY Distal Access Guiding Catheter relates to the provision of a minimally invasive diagnostic or treatment option to patients through facilitating the delivery of diagnostic or interventional devices endovascularly as an alternative to open surgery.

A Summary of Safety and Clinical Performance (SSCP) can be found at the following link:
<https://ec.europa.eu/tools/eudamed>.

This is the SSCP location after the launch of EUDAMED.

Note: The Summary of Safety and Clinical Performance and the EUDAMED website are not reviewed by USFDA.

Contraindications

None known.

Warnings

Discard the guiding catheter after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Do not use with Ethiodol™ or Lipiodol™ contrast media, or other such contrast media which incorporates the components of these agents.

Precautions

Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged,
- contents appear damaged, or
- the expiry date has passed.
- Store in a dark, dry place between 15° and 30°C (59° and 86°F).
- Exposure to temperatures above 54°C (130°F) may damage the catheter.
- Do not expose to organic solvents.
- Inspect the guiding catheter before use to verify that its size, shape, and condition are suitable for the specific procedure.
- If strong resistance is met during manipulation, discontinue the procedure and determine the cause for the resistance before proceeding. If the cause of the resistance cannot be determined, withdraw the catheter.
- Limit the exposure to X-ray radiation doses to patients and physicians by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors when possible.
- Extreme care must be taken to avoid damage to the vasculature through which the catheter passes. The guiding catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.

- Torquing the catheter excessively while kinked may cause damage which could result in possible separation along the catheter shaft. Should the guiding catheter shaft become severely kinked, withdraw the entire system (guiding catheter, guidewire, and catheter sheath introducer).
- Advancement, manipulation, and withdrawal of the guiding catheter should always be performed under fluoroscopic guidance.
- Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

Potential Adverse Events and Reporting

Possible complications include, but are not limited to the following:

- acute occlusion
- air embolism
- death
- distal embolization
- emboli
- false aneurysm formation
- hematoma or hemorrhage at puncture site
- infection
- intracranial hemorrhage
- ischemia
- neurological deficits including stroke
- vessel spasm, thrombosis, dissection, or perforation

Use of this device requires fluoroscopy which presents potential risks to physicians and patients associated with x-ray exposure. Possible risks include, but are not limited, to the following:

- alopecia
- burns ranging in severity from skin reddening to ulcers
- cataracts
- delayed neoplasia

Note: If a serious incident occurs during the use of this device or as a result of its use, please report it to the manufacturer and/or its authorized representative and to your national authority.

Recommended Procedure

1. Gently remove the ENVOY DA Guiding Catheter from its packaging.
2. Inspect the catheter upon removal from packaging to verify that it is undamaged.
WARNING: Do not use a catheter that has been damaged in any way. If damage is detected, replace with another ENVOY DA Guiding Catheter that is not damaged.
3. Prior to use, flush the lumen of the Guiding Catheter with a heparinized saline solution. Immerse 10 cm of the distal end of the catheter in the heparinized saline to wet the hydrophilic coating.
4. Use the supplied peel away introducer to aid in the insertion of the catheter through the catheter sheath introducer.
5. Introduce the guiding catheter into the vasculature through the catheter sheath introducer, and over a .038" (0.97 mm) or .035" (0.89 mm) guidewire using a percutaneous entry technique of choice.
6. Set up a continuous heparinized saline flush through the sidearm of a rotating hemostatic valve attached to the Luer lock fitting of the guiding catheter hub.
Note: It is recommended that continuous heparinized saline flush be maintained between the ENVOY DA Guiding Catheter and any intraluminal device passed coaxially through it.
7. Under fluoroscopic guidance, advance the guiding catheter over the guidewire until the desired position is attained.
8. Remove the guidewire prior to the introduction of other intravascular devices or infusion of contrast agents.

How Supplied

②

The ENVOY Distal Access Guiding Catheter is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Cerenovus will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner unit is not opened or damaged, the product is sterile and nonpyrogenic.

Disposal

Dispose of all used devices in accordance with hospital policy. Do not discard this product as unsorted municipal waste. Contact your local authorities for disposal instructions or return it to Cerenovus for disposal.

Warranty

Cerenovus warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

TM ENVOY is a trademark of Cerenovus or affiliated companies.

TM Ethiodol and Lipiodol are trademarks of Guerbet S.A.


Visit www.e-ifu.com for the latest version of this document. For technical assistance, contact your local sales representative.

Interpretation of symbols

	Radiopaque
	Do not Autoclave
	Do not resterilize
	Do not reuse
	Do not use if package is damaged
	Manufacturer
	Nonpyrogenic
	Consult instructions for use
	Caution
	Quantity

Medical Device

Temperature limits 15-30°C
(59-86°F)

Packaging Unit

Authorized European
Representative

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© Cerenovus 2021. All rights reserved.