Cerenovus Microcatheters

[Cerenovus Microcatheters](#)

en

[Cerenovus mikrokatetri](#)

lv

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[„Cerenovus“ mikrokateteriai](#)

lt

This document is valid only on the date printed. If unsure of the print date, please re-print to ensure use of the latest revision of the IFU (available at www.e-ifu.com). The onus resides with the user to ensure that the most up-to-date IFU is used.

en

Cerenovus Microcatheters

IMPORTANT INFORMATION

Please read before use

Some products may not be licensed in all jurisdictions.

STERILE	EO
---------	----

Rx Only

MD

Description

The Cerenovus Microcatheters are variable stiffness, single lumen catheters designed to access small, tortuous vasculature. They are available in a variety of outer and inner diameters. Each configuration has a hydrophilic coating to provide lubricity for navigation of vessels. The inner lumen is lined with lubricious PTFE to facilitate movement of guidewires and other devices. The distal sections of the catheter bodies are radiopaque to aid visualization under fluoroscopy, and the distal tips are clearly distinguished by a radiopaque marker.

Product Features

RAPIDTRANSIT™ Microcatheter

PROWLER™ PLUS Microcatheter

PROWLER™ SELECT™ PLUS Microcatheter

	French	Inch	mm
Proximal O.D.	2.8	--	.95
Distal O.D.	2.3	--	.75
Proximal I.D.		.021	.5
Distal I.D.		.021	.5
Min. I.D. Guiding Catheter		.042	1.1
Max. O.D. Guidewire		.018	.46
Max. PVA Compatibility			1000 µm

PROWLER™ 14 Microcatheter

PROWLER™ SELECT™ LP-ES Microcatheter

	French	Inch	mm
Proximal O.D.	2.3	--	.75
Distal O.D.	1.9	--	.65
Proximal I.D.		.0165	.4
Distal I.D.		.0165	.4
Min. I.D. Guiding Catheter		.035	.9
Max. O.D. Guidewire		.014	.36
Max. PVA Compatibility			500 µm

PROWLER™ 10 Microcatheter

	French	Inch	mm
Proximal O.D.	2.3	--	.75
Distal O.D.	1.7	--	.55
Proximal I.D.		.015	.4
Distal I.D.		.015	.4
Min. I.D. Guiding Catheter		.035	.9
Max. O.D. Guidewire		.012	.30
Max. PVA Compatibility			500 µm

Intended Use

The Cerenovus Microcatheter is intended as a conduit for endovascular use during diagnostic and/or interventional procedures in the peripheral and neuro vasculature.

Indications

The Cerenovus Microcatheter is indicated for use in the infusion/introduction of interventional agents/devices in the peripheral and neuro vasculature.

Patient Population

The Cerenovus Microcatheter is intended for use in patients undergoing diagnostic and/or interventional procedures.

Intended User

The device should be used only by physicians trained in endovascular procedures at medical facilities with the appropriate imaging equipment.

Clinical Benefits (Applicable only in markets subject to EU Regulation 2017/745)

The Cerenovus Microcatheter provides conduit for endovascular catheter-based interventions for less invasive alternatives to open surgery in patient diagnostic and treatment procedures.

A Summary of Safety and Clinical Performance (SSCP) can be found at the following link:

<https://ec.europa.eu/tools/eudamed>.

This is the SSCP location after the launch of EUDAMED.

Note: The Summary of Safety and Clinical Performance and the EUDAMED website are not reviewed by USFDA.

Contraindications

The Cerenovus Microcatheters are contraindicated for delivery of infusates that contain dimethyl sulfoxide (DMSO).

Warnings

THIS DEVICE IS INTENDED FOR ONE USE ONLY. Discard the microcatheter after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Never advance or withdraw an intraluminal device against resistance until the cause of resistance is determined by fluoroscopy. If the cause cannot be determined, withdraw the catheter. Movement against resistance may result in damage to the vessel.

The infusion pressure should not exceed the maximum listed pressure for each catheter, as indicated in the flowrate charts. Pressure in excess of the recommended range may result in catheter rupture or tip severance.

Precautions

- Store in a dark, dry place between 15° and 30°C (59° and 86°F).
- Do not use if package is open or damaged.
- Use prior to the “Use By” date.
- Read and follow the Instructions for Use of all agents or contrast media used with the microcatheters.

Potential Adverse Events and Reporting

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications, which may occur during or after the procedure.

Possible complications include, but are not limited to the following:

- Access site complications
- Allergic reaction
- Embolism
- Hemorrhage
- Infection
- Neurological deficit including stroke and death
- Vascular injury including but not limited to pseudoaneurysm, vasospasm, vessel dissection, and vessel perforation

Use of this device requires fluoroscopy which presents potential risks to physicians and patients associated with x-ray exposure. Possible risks include, but are not limited, to the following:

- Alopecia
- Burns ranging in severity from skin reddening to ulcers
- Cataracts
- Delayed neoplasia

Note: If a serious incident occurs during the use of this device or as a result of its use, please report it to the manufacturer and/or its authorized representative and to your national authority.

RAPIDTRANSIT™ Microcatheter

FLOWRATE	100%	50/50	100%
CHART	Contrast-76%	Contrast-76%/	Saline
.021" I.D.	(cc/sec)	Saline (cc/sec)	(cc/sec)

Approximate Values

Total Catheter Length (cm)	Dead Space (cc)	@ 100 psi (689 kPa)	@ 200 psi (1379 kPa)	@ 100 psi	@ 200 psi	@ 100 psi	@ 200 psi
175	0.55	0.1	0.1	0.3	0.5	0.7	1.1
155	0.50	0.1	0.2	0.4	0.8	0.8	1.3
140	0.46	0.1	0.2	0.5	0.8	0.9	1.4
105	0.39	0.1	0.2	0.6	1.1	1.1	1.5
75	0.30	0.2	0.3	0.8	1.3	1.3	1.8

PROWLER™ PLUS Microcatheter
PROWLER™ SELECT™ PLUS Microcatheter

FLOWRATE	100%	50/50	100%
CHART	Contrast-76%	Contrast-76%/	Saline
.021" I.D.	(cc/sec)	Saline (cc/sec)	(cc/sec)

Approximate Values

Total Catheter Length (cm)	Dead Space (cc)	@ 100 psi (689 kPa)	@ 300 psi (2068 kPa)	@ 100 psi	@ 300 psi	@ 100 psi	@ 300 psi
175	0.55	0.1	0.1	0.3	0.8	0.7	1.5
165	0.53	0.1	0.1	0.4	1.0	0.7	1.6
155	0.50	0.1	0.2	0.4	1.1	0.8	1.7
140	0.46	0.1	0.2	0.5	1.2	0.9	1.9
115	0.40	0.1	0.2	0.6	1.5	0.9	1.9
105	0.39	0.1	0.3	0.6	1.5	1.1	1.9
90	0.34	0.1	0.4	0.7	1.7	1.2	2.2
75	0.30	0.2	0.5	0.8	1.8	1.3	2.4

PROWLER™ 14 Microcatheter
PROWLER™ SELECT™ LP-ES Microcatheter

FLOWRATE	100%	50/50	100%
CHART	Contrast-76%	Contrast-76%/	Saline
.0165" I.D.	(cc/sec)	Saline (cc/sec)	(cc/sec)

Approximate Values

Total Catheter Length (cm)	Dead Space (cc)	@ 100 psi (689 kPa)	@ 300 psi (2068 kPa)	@ 100 psi	@ 300 psi	@ 100 psi	@ 300 psi
175	0.38	0.02	0.1	0.2	0.4	0.3	0.7
155	0.35	0.04	0.1	0.2	0.5	0.4	0.7
140	0.35	0.05	0.2	0.2	0.6	0.4	0.9
75	0.25	0.1	0.2	0.3	0.7	0.5	1.0

PROWLER™ 10 Microcatheter

FLOWRATE	100%	50/50	100%
CHART	Contrast-76%	Contrast-76%/	Saline
.015" I.D.	(cc/sec)	Saline (cc/sec)	(cc/sec)

Approximate Values

Total Catheter Length (cm)	Dead Space (cc)	@ 100 psi (689 kPa)	@ 300 psi (2068 kPa)	@ 100 psi	@ 300 psi	@ 100 psi	@ 300 psi
175	0.34	0.03	0.1	0.1	0.3	0.2	0.6
155	0.32	0.04	0.1	0.1	0.3	0.2	0.6
75	0.23	0.1	0.1	0.2	0.6	0.4	0.9

Recommended Procedure

It is recommended that the Cerenovus Microcatheters be used with a guiding catheter, a compatible catheter sheath introducer, and a steerable guidewire.

1. The Cerenovus Microcatheters may be packaged in a protective coil dispenser fitted with a flushing luer connector or in a tray. The catheters have a hydrophilic coating and require hydration prior to removal from the coil dispenser and prior to use.

- Before removing the catheter from the coil dispenser, flush the dispenser with heparinized saline through the luer connector fitting on the end of the coil dispenser.

Note: The catheter may be hydrated in the tray prior to use.

- Remove catheter and inspect to verify that it is undamaged.

Note: Do not attempt to use microcatheters without flushing first to hydrate the coating. Failure to do so may compromise the coating and lubricity of the catheter.

WARNING: Do not use a catheter that has been damaged in any way. Damaged catheters may rupture causing vessel damage or tip detachment during the procedure.

2. Microcatheters may be packaged with a shaping mandrel. In order to maintain catheter burst integrity and dimensional stability, it is strongly recommended that the user follow these instructions when shaping a catheter.
 - Remove shaping mandrel from mounting card and insert into distal tip of catheter.
 - Bend catheter tip and shaping mandrel into desired shape. Over exaggeration of the desired shape is recommended to accommodate for slight catheter relaxation.
 - Hold shaping mandrel/catheter assembly directly over the steam source for approximately 30 seconds to set shape.
 - Remove shaped catheter assembly from heat source and allow to cool in either air or liquid prior to removing the mandrel.
 - Remove shaping mandrel from catheter and discard.

WARNING: Do not use shaping mandrel *in-vivo*.
3. Prior to use, flush the catheter lumen with heparinized saline solution by attaching a saline filled syringe to the catheter hub.
4. Remove the appropriate steerable guidewire from its package and inspect for damage.

5. Carefully insert guidewire into the funneled hub of the catheter and advance into the catheter lumen.

6. Place the appropriate guiding catheter using a percutaneous entry technique of choice. Connect a hemostatic valve to the guiding catheter hub and maintain a continuous heparinized saline flush.

Note: Microcatheters require a continuous flush during the procedure in order to maintain the lubricity of the coating.

7. Introduce the guidewire and microcatheter as a unit through the hemostatic valve into the lumen of the guiding catheter. Advance guidewire/microcatheter assembly to the distal tip of the guiding catheter.
8. Alternatively advance the guidewire and microcatheter until the desired site has been accessed.

Note: To facilitate catheter handling, the proximal portion of the microcatheters are uncoated to ensure a non-slip grip.

CAUTION: If strong resistance is met during manipulation, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, withdraw the catheter and guidewire as a system.

9. When ready to infuse, withdraw the guidewire completely from the catheter. Connect a syringe containing infusate to the microcatheter hub and infuse according to the manufacturer's instructions and precautions.

WARNING: If flow through the catheter becomes restricted, do not attempt to clear the catheter lumen by infusion. Determine and remedy the cause of blockage or replace the blocked catheter with a new catheter before resuming infusion.

Note: Infusion pressure should not exceed the maximum pressure rating.

10. After completing the procedure, withdraw the catheter and discard.

Disposal

Dispose of all used devices in accordance to the hospital policy. Do not discard this product as unsorted municipal waste. Contact your local authorities for disposal instructions or return it to Cerenovus for disposal.

Warranty

Cerenovus warrants that this medical device is free from defects in both materials and workmanship.

Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.

™ PROWLER and RAPIDTRANSIT are trademarks of Cerenovus or affiliated companies.

Visit www.e-ifu.com for the latest version of this document. For technical assistance, contact your local sales representative.

Interpretation of symbols

Do not re-sterilize

Do not use if package is open or damaged

Manufacturer

Consult instructions for use


Caution

Nonpyrogenic

Medical Device

**Temperature limits 15-30°C
(59-86°F)**

Packaging Unit

**Authorized European
Representative**

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© Cerenovus 2021-2022. All rights reserved.