EMBOTRAP™ III

EN

INSTRUCTIONS FOR USE


INDICATIONS FOR USE:

The EMBOTRAP™ III Revascularization Device (the Device) is intended to be used to restore blood flow in patients experiencing an acute ischemic stroke due to a large vessel neurovascular occlusion. The Device is designed for use in the anterior and posterior neurovasculature in vessels, such as the internal carotid artery, the M1 and M2 segments of the middle cerebral artery, the A1 and A2 segments of the anterior cerebral artery, the basilar, the posterior cerebral and the vertebral arteries.

The Device should only be used by physicians trained in neurointerventional catheterization and the treatment of ischemic stroke.

DEVICE DESCRIPTION:

The Device is composed of a three dimensional nitinol stent-like assembly at the distal end of a tapered nitinol shaft as shown in Figure 1 below. The Device is supplied pre-loaded in an Insertion Tool.

**Figure 1: Device Illustration**

1. Inner Channel 2. Proximal Markers 3. Outer Cage 4. Body Markers 5. Distal Markers 6. Distal Mesh

Table 1: Device Length Details
(From Figure 1)

Table 1: Device Length Details (From Figure 1)					Radiopaque Features							
Catalogue Number	A	B	C	D	E	F	Number of Radiopaque Markers			Indicated Vessel Range		
	Outer Cage Diameter	Working Length	Outer Cage Length	Overall Length	Radiopaque Tip Length	Proximal Radiopaque Coil Length	Proximal	Body	Distal			
	ET307522	5mm	22mm	34mm	194cm	4mm	20mm	2	8		3	1.5 – 5mm
	ET307537	5mm	37mm	49mm	195cm	4mm	20mm	2	16		3	1.5 – 5mm
ET307645	6.5mm	45mm	57mm	196cm	4mm	20mm	2	16	3	1.5 – 6.5mm		

COMPATIBILITY:

Microcatheter: The Device should be introduced through a microcatheter with an internal diameter (ID) of 0.021"- 0.027". See *Positioning and Deployment* for further information.

Balloon Guide Catheter / Guide Catheter: The Device should be used in conjunction with a Neurovascular Guide or Sheath. The guide/sheath chosen is to be compatible with the selected microcatheter.

Additional Accessories: The following devices may be selected for use during the procedure. Only accessory devices which are compatible with the selected microcatheter and guide catheter combination should be used.

- Intermediate Catheter / Aspiration Catheter
- Rotating Hemostasis Valve (RHV)

Note: RHV chosen should be appropriate for use with the EMBOTRAP III Insertion Tool, i.e. the insertion tool is required to seat fully into the microcatheter hub and be secured firmly by the RHV. See *Preparation and Delivery* for further information.

CONTRAINDICATIONS:

- Allergy or hypersensitivity to Nickel-Titanium.
- Excessive vessel tortuosity that may prevent device delivery or withdrawal.
- Patients with angiographic evidence of carotid dissection.

POTENTIAL COMPLICATIONS:

Possible complications include but are not limited to: air embolism, hematoma or hemorrhage at puncture site, infection, distal embolization, new territory embolization, vascular spasm, neurological deterioration including stroke and death, thrombosis, dissection or perforation, device malfunction or fracture, vascular occlusion, ischemia, intracranial hemorrhage, pseudoaneurysm formation, or other minor or major complications related to cerebral angiography.

WARNINGS & PRECAUTIONS: ⚠

- Store in a cool, dry place.
- Do not use the Device if the product sterilisation barrier or its packaging is compromised.
- Do not use after the expiration date printed on the product label.
- Before use inspect the Device carefully for damage. Do not use a damaged or kinked device.
- Do not deploy in vessels of less than 1.5mm diameter.
- Do not deploy in calcified lesions.
- The use of stent-retrievers in the removal of thrombus due to atherosclerotic lesions has not been evaluated and therefore the use of EMBOTRAP III in such lesions is not recommended.
- Do not advance the Device through the microcatheter against significant resistance.
- Do not torque or rotate the Device.
- Do not advance the Device after it has been deployed or partially deployed from the microcatheter.
- Do not withdraw the Device against significant resistance. Assess the cause of resistance using fluoroscopy and if necessary, advance the microcatheter over the Device to resheath or partially resheath to aid withdrawal. Resheathing a Device with captured clot may result in loss of clot and distal embolization.
- This Device is for use in one patient only. Do not resterilize or reuse in more than one patient. Reusing the Device in another patient could result in compromised Device performance, cross infection and other safety related hazards.
- Do not use any EMBOTRAP™ III Revascularization Device more than 3 times.
- Do not perform more than 3 retrieval attempts in any one vessel.
- Do not withdraw the deployed Device through a previously stented vessel. If it becomes necessary to cross a deployed stent to access a clot, first ensure that the stent can be crossed with the Guide or Intermediate Catheter. The Device should then be fully pulled back into Guide or Intermediate Catheter prior to retrieval through the stent.
- Users shall take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible. The risks of angiographic and fluoroscopic X-ray

radiation include alopecia, burns ranging in severity from skin reddening to ulcers, cataracts and delayed neoplasia that increase in probability as procedure time and number of procedures increase.

- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

PROCEDURAL INFORMATION:

The Device is delivered endovascularly under fluoroscopic guidance in a similar manner to that of other neurovascular clot-retrieval systems. Physicians are advised to select devices which are appropriately sized for the estimated occlusion length and location in the vasculature in accordance with the information presented in Table 1.

The use of antiplatelet and anticoagulation agents is recommended at the discretion of the treating physician.

Preparation & Delivery:

1. Using standard interventional techniques, access the arterial system and using angiography, determine the location of the occluded vessel.
2. Advance an appropriate Guide Catheter, Sheath or Balloon Guide Catheter as close to the occlusion site as possible. Connect a rotating hemostasis valve (RHV – see section “Compatibility”) to the proximal end of this catheter and connect to a continuous flush system.
Note: Depending on tortuosity and/or occlusion location the use of an Intermediate Catheter may be considered to provide additional support for ease of retrieval.
3. Select an appropriate microcatheter (see section “Compatibility”). Connect an RHV to the proximal end of the microcatheter and connect to a continuous flush system.
4. With the aid of a suitable guidewire, and using standard catheterization techniques and fluoroscopic guidance, advance the microcatheter up to and across the occlusion so that the distal end of the microcatheter is positioned distal to the occlusion.
5. Remove the guidewire from the microcatheter and, if desired, gently infuse contrast media through the microcatheter to visualize the distal end of the occlusion.
6. Remove the Insertion Tool with the preloaded Device from the packaging hoop, taking care not to accidentally retract or advance the Device from the Insertion Tool while doing so.
7. Insert the distal end of the Insertion Tool through the RHV of the microcatheter and wait until fluid is seen exiting the proximal end of the Insertion Tool, confirming that the Device is flushed. Advance the Insertion Tool until it is fully seated in the hub of the microcatheter. Fully tighten the RHV to hold the Insertion Tool securely in position.
8. Confirm that the Insertion Tool is fully seated in the hub of the RHV before proceeding to advance the Device. Advance the Device until at least half of the shaft length has been inserted into the microcatheter, at which point the Insertion Tool may be removed.
9. Do not discard the Insertion Tool as it will be needed if additional device passes are required during the procedure.

Positioning and Deployment:

1. Continue to advance the Device until the distal radiopaque tip of the Device approaches the distal region of the microcatheter. Care should be taken not to advance the distal tip of the Device out of the microcatheter.
Note: There are two silver bands positioned on the proximal portion of the Device shaft to assist in minimizing the amount of fluoroscopic exposure required during Device insertion. For example, if using a standard microcatheter (working length of 155cm and a 7cm RHV), then when the first silver band on the shaft approaches the RHV, the Device tip will be approximately 5cm to 8cm from the distal end of the microcatheter (when using the longest or shortest models, respectively). When the second silver

band on the Device shaft approaches the RHV the Device tip will be nearing the distal end of the microcatheter.

2. Advance the Device in the microcatheter and position the EMBOTRAP™ III such that the clot lies within the working length of the device, i.e. between the most proximal and most distal radiopaque markers.

Where possible, and if appropriate, position the Device proximally such that the proximal radiopaque markers align with the proximal margin of the clot (Figure 2).


Figure 2: Proximal positioning of EMBOTRAP™ III relative to clot
(5x22mm model shown as an example)

3. While stabilizing the Device shaft, retract the microcatheter to deploy the Device within the clot. The Device will be fully deployed once the distal tip of the microcatheter has been retracted to a position over the proximal radiopaque coil of the Device.
4. The Device may be left to embed within the clot for 3-5 minutes prior to withdrawal.

Retrieval:

1. Optimum control of the Device and microcatheter relative positions during retrieval may be achieved by advancing a torque device over the Device shaft and tightening it against the microcatheter RHV.
2. If using a balloon guide catheter inflate the balloon to occlude the vessel as per the balloon guide catheter IFU.
3. Withdraw the Device and microcatheter slowly and carefully as a unit to the guide catheter tip while aspirating through the guide catheter with a syringe.
4. As the device reaches the guide catheter apply vigorous aspiration, withdraw the Device and microcatheter into the guide catheter and continue to aspirate until the Device reaches the proximal RHV.

Note: If withdrawal into the guide catheter is difficult (as may be the case with a large clot burden) then deflate the balloon (if applicable) and withdraw the guide, microcatheter and Device together through the introducer sheath.

5. Disconnect the RHV from the guide catheter and remove the Device, microcatheter and RHV together from the guide catheter
6. Use a syringe to further aspirate the guide catheter to ensure it is clear of any thrombus material.

Cleaning and Re-use:

- The Device may be used for up to three retrieval attempts.
- If an additional pass is to be made with the Device then carefully remove any captured thrombus from the Device and clean the Device in heparinized saline, rubbing gently from proximal to distal to remove any residual thrombus material. Inspect the Device carefully and if any damage or deformity is observed do not use the Device – instead use a new EMBOTRAP™ III Revascularization Device for any subsequent passes. If using the same Device, replace the Insertion Tool onto the proximal shaft and retract the Device into the Insertion Tool until it is fully resheathed. Do not retract the device into the Insertion Tool against significant resistance –ascertain the cause of the resistance and if required use a new device for any subsequent passes.
- Do not attempt more than three retrieval attempts in one vessel.


Serious Incident Reporting:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- a) the death of a patient, user or other person,
- b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- c) a serious public health threat;

INFORMATIVE SYMBOLS FOR EMBOTRAP™ III:

	Warnings – see instructions for use		Lot Number
	Use only once – Do not re-use		Use By
	Sterilised using ethylene oxide		Contents
	Catalogue Number		Manufacturer, Date of Manufacture
	Minimum microcatheter inner diameter (ID)		Do not use if package is damaged
	Keep Dry		Consult Instructions for Use
	Keep away from sunlight/heat		Single Sterile Barrier System
	Unique Device Identifier		Do not re-sterilize
	Medical Device		

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Neuravi LTD