AGILITY® Steerable Guidewires

Medos International SÄRL
Chemin-Blanc 38, 2400
Le Locle, Switzerland

CE
2797

LCN 206431-001/D
© 2009–2018 Cerenovus
All rights reserved
Recycled symbol Revised 12/18

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

AGILITY® Steerable Guidewires

STERILE	EO
---------	----

Device Description

The hydrophilically coated **AGILITY®** Steerable Guidewires consist of a stainless steel wire core and a radiopaque platinum/tungsten coil on the distal tip. Guidewire length, diameter, and distal tip configuration are indicated on the product label. A steering/torquing device and a guidewire introducer are packaged with the **AGILITY** Guidewire.

Indications

The Cerenovus **AGILITY** Steerable Guidewires are intended for selective placement of microcatheters and other devices in the neuro and peripheral vasculature.

Contraindications

AGILITY Steerable Guidewires are contraindicated for use in chronic total occlusions in the peripheral vasculature.

Warnings

Do not reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Guidewires are delicate instruments and should be handled carefully. Special care should be taken when shaping the guidewire tip in order to prevent damage. Prior to use and when possible during the procedure, inspect the guidewire carefully for coil separation, bends, or kinks. Do not use a guidewire that shows signs of damage. Damage will prevent the guidewire from performing with accurate torque response and control.

- ***Guidewire manipulation/torquing should always be performed under fluoroscopic guidance.***

- ***Never advance, withdraw or auger the guidewire against resistance without first determining the cause of resistance under fluoroscopy. Torquing the guidewire against resistance may cause damage and/or fracture which may result in separation of the distal tip.***
- ***Should the guidewire tip become entrapped within the vasculature (i.e., small side branch), DO NOT TORQUE THE GUIDEWIRE. Advance the microcatheter distally, gently pull the guidewire back into the microcatheter, and remove the microcatheter/guidewire system as a unit.***
- ***Should torque control/tip response be compromised during use, confirm tip integrity using fluoroscopy. LOSS OF TORQUE CONTROL MAY BE DUE TO CORE WIRE FRACTURE. Under fluoroscopic guidance, advance the microcatheter to the distal end of the guidewire and remove the microcatheter/guidewire system as a unit.***
- ***The steerable guidewire should not be used during magnetic resonance procedures.***

Precautions

- Store in a cool, dark, dry place.
- Do not use open or damaged packages.
- Use prior to the "Use By" date.
- Do not autoclave. Exposure to temperatures above 54°C (130°F) may damage the guidewire.
- Do not expose to organic solvents.
- Movement of torque device or metal insertion tool on a guidewire's coating may compromise the integrity of the coating.

Complications

Procedures requiring percutaneous guidewire introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

Possible complications include, but are not limited to:

- infection
- hemorrhage
- emboli
- vasospasm
- perforation of the vessel wall
- ischemia and related neurological dysfunction, stroke, and death
- dissection of the vessel wall

How Supplied

The **AGILITY** Steerable Guidewires are intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Cerenovus will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

The device components have been tested and were determined to be nonpyrogenic, except for the torquing device, which is not tested.

Recommended Procedure

1. Use appropriately sized steerable guidewires with the selected microcatheter (refer to the microcatheter for compatible guidewire size).
2. Open the sterile package slowly. To activate the hydrophilic coating of the **AGILITY** Guidewire, and to aid in its removal from the dispensing tube, it is essential to thoroughly flush the dispensing tube by infusing saline into the Luer end of the dispensing tube until saline is flushed through at the other end. To prevent damage to the distal tip of the guidewire and to prevent the guidewire from springing onto a non-sterile field, carefully remove the guidewire from the dispensing tube. Do not allow product to dry and do not reinsert product into dispensing tube after it is hydrated.
3. Insert the distal tip of the guidewire into the microcatheter hub (or through the rotating hemostasis valve which is attached to the microcatheter hub) and advance the guidewire until the tip exits the microcatheter tip. The guidewire introducer may be used to aid introduction of the guidewire into the microcatheter hub.
Caution: Gently introduce and advance the guidewire to prevent damaging the distal tip.
Note: To aid in the navigation within and selection of arteries, the guidewire may be pre-shaped. If desired, gently shape the guidewire's distal, flexible tip according to standard practices. The guidewire introducer (packaged with the device) may be used for this purpose.

4. To aid in rotating or steering the guidewire, secure a torquing device to the proximal end of the guidewire.
Caution: To prevent vessel injury or *tip entrapment*, use fluoroscopy when advancing/torquing the guidewire.
5. Use accepted angiographic technique for guidewire positioning. Then continue the procedure.
Caution: If strong resistance is met during manipulation, discontinue the procedure and determine cause of resistance before proceeding. If the cause of resistance cannot be determined, withdraw the entire system. ***If the tip becomes fixed within the vasculature, advance the microcatheter as distally as possible, gently pull the guidewire back into the microcatheter and withdraw the microcatheter/guidewire system as a unit - NEVER TORQUE THE GUIDEWIRE IF THE TIP BECOMES FIXED.***
6. When desired procedure results are achieved, withdraw the guidewire and microcatheter slowly.
Note: When the guidewire is in the microcatheter's inner lumen, pressure measurement and contrast material injection through the inner lumen may be impaired. Consult the microcatheter's instructions for the recommended procedure in these circumstances.

Warranty

Cerenovus warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

® AGILITY is a registered trademark of Cerenovus or affiliated companies

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