CERENOVUS

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

CEREBASE DA Guide Sheath

Medos International SÀRL

Chemin-Blanc 38, 2400

Le Locle, Switzerland

CE
2797

LCN 208918-001/A

© 2020 Cerenovus.

All rights reserved.

 Printed 01/20

IMPORTANT INFORMATION

Please Read Before Use

CEREBASE DA Guide Sheath

STERILE

EO

Description

The CEREBASE DA Guide Sheath package contains the guide sheath, a hemostasis valve, and a dilator.

The CEREBASE DA Guide Sheath is a single lumen, variable stiffness guide sheath. It has a radiopaque marker band on the distal end and a luer hub at the proximal end. The guide sheath has an outer hydrophilic coating on the distal 20 cm that reduces friction during use. The lubricious PTFE lined inner lumen is designed to facilitate delivery of guide wires, catheters, and other devices. Please refer to the label for short sheath compatibility.

The dilator is a single lumen radiopaque catheter with a tapered distal end and a luer hub at the proximal end. The dilator is compatible with the CEREBASE DA Guide Sheath and up to 0.038 inch diagnostic guidewires. The dilator facilitates the percutaneous transition from the skin through the subcutaneous tissue to the vessel.

Indications

The CEREBASE DA Guide Sheath is indicated for the introduction of interventional devices into the neuro vasculature.

Contraindications

None known.

WARNINGS

- The CEREBASE DA Guide Sheath should only be used by physicians trained in neuro interventional endovascular procedures at medical facilities with the proper imaging equipment.
- Do not use with Ethiodol or Lipiodol contrast media, or other such contrast media which incorporates the components of these agents.
- Do not advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the device against resistance may cause vessel injury or damage to the device.
- Discard the sheath after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. This device is extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- Dispose of all used devices in accordance with hospital policy for biohazardous materials.
- To prevent thrombus formation and contrast media crystal formation, maintain a constant infusion of an appropriate flush solution through the catheter lumen.
- This device is coated with a hydrophilic coating at the distal end of the device for a length of 20 cm. Please refer to the *Instructions for Use* Section for further information on how to prepare and use this device to ensure it performs as intended. Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

Precautions

Inspect the sterile package carefully. Do not use if:

- o The package or seal appears damaged,
- o Contents appear damaged, or
- o The expiry date has passed.
- Carefully inspect all devices prior to use. Verify size, length, and condition are suitable for the specific procedure. Do not use a CEREBASE DA Guide Sheath, dilator, or hemostasis valve that has been damaged in any way; replace with another guide sheath, dilator, or hemostasis valve. A damaged device may cause complications.
- Exercise care in handling the CEREBASE DA Guide Sheath to reduce the chance of accidental damage.
- To control the proper introduction, movement, positioning and removal of the guide sheath within the vascular system, users should employ standard clinical angiographic and fluoroscopic practices and techniques throughout the interventional procedure.
- If removed from the patient, the hydrophilic coating on the guide sheath must remain hydrated with heparinized saline. Do not allow the coating to dry.
- Avoid wiping the device with dry gauze as this may damage the device coating.
- Avoid excessive wiping of the coated device.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the coating which could affect the device safety and performance.
- Limit the exposure to X-ray radiation doses to patients and physicians by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors when possible.

Adverse Events

Potential adverse events associated with the use of catheters or with endovascular procedures include, but are not limited to:

- Access site complications
- Allergic or anaphylactic reaction
- Arrhythmia
- Cardiac complications
- Cerebral edema
- Cerebral Infarct
- Death

- Embolism (air, foreign body, plaque, thrombus)
- Hematoma
- Hemorrhage
- Hypo/hypertension
- Infection
- Inflammation or granuloma at access site
- Ischemia
- Neurological deficits
- Pseudoaneurysm
- Renal insufficiency
- Respiratory complications
- Seizure
- Stroke
- Tissue necrosis
- Transient Ischemic Attack
- Vasospasm
- Vessel dissection
- Vessel occlusion
- Vessel perforation
- Vessel rupture
- Vessel thrombosis

Use of device requires fluoroscopy which presents potential risks to physicians and patients associated with x-ray exposure. Possible risks include, but are not limited to, the following:

- Alopecia
- Burns ranging in severity from skin reddening to ulcers
- Cataracts
- Delayed neoplasia

How Supplied

The CEREBASE DA Guide Sheath is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Cerenovus will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner unit is not opened or damaged, the product is sterile and nonpyrogenic.

Instructions for Use

1. Select an appropriately sized CEREBASE DA Guide Sheath based on the patient anatomy and length.
2. Gently remove the device and accessories from the packaging.
3. Inspect the guide sheath upon removal from packaging to verify that it is undamaged.
CAUTION: Do not use a guide sheath, dilator, or hemostasis valve that has been damaged in any way. If damage is detected, replace with another guide sheath, dilator, or hemostasis valve that is not damaged.
4. Prior to use, flush the lumen of the guide sheath and wet the outer body of the catheter with a heparinized saline solution. Connect a hemostasis valve to the hub of the guide sheath. Set up a continuous heparinized saline flush through the hemostasis valve.
5. If using as a primary access sheath with the dilator, flush and wet the dilator with saline. Insert the dilator and advance until the tapered distal tip is beyond the distal end of the CEREBASE DA Guide Sheath. Close the hemostasis valve around the dilator.
6. Gain primary access to the vasculature using a standard technique of choice.
7. Under fluoroscopy, advance the guide sheath and/or guide sheath/dilator assembly over the guidewire and into vasculature. If using the dilator, loosen the hemostasis valve and remove the dilator.
8. Insert appropriately sized catheter and/or guidewire as needed.
9. Under fluoroscopic guidance, advance the CEREBASE DA Guide Sheath over the catheter and/or guidewire through the vasculature to the desired location.
WARNING: Do not advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy.
10. Remove the catheter and/or guidewire.
11. When procedure is complete, remove the CEREBASE DA Guide Sheath using standard technique.
12. If re-access to the vasculatures with the same device is desired, flush and clean the inner lumen of the device by infusion. Inspect the device for any damage.
Caution: Do not use the device if any damage or irregularities are observed.

Storage

Store in a cool, dark, dry place.

Warranty

Cerenovus warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

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4. Prior to use, flush the catheter with saline solution and prime the catheter with heparinized saline solution. Connect the catheter to the saline solution. Flush the catheter with saline solution. Flush the catheter with saline solution.
 5. If you use the catheter as a primary catheter for dialysis, flush the catheter with saline solution and prime the catheter with heparinized saline solution. Connect the catheter to the saline solution. Flush the catheter with saline solution. Flush the catheter with saline solution.
 6. If you use the catheter as a primary catheter for dialysis, flush the catheter with saline solution and prime the catheter with heparinized saline solution. Connect the catheter to the saline solution. Flush the catheter with saline solution. Flush the catheter with saline solution.
 7. If you use the catheter as a primary catheter for dialysis, flush the catheter with saline solution and prime the catheter with heparinized saline solution. Connect the catheter to the saline solution. Flush the catheter with saline solution. Flush the catheter with saline solution.
 8. If you use the catheter as a primary catheter for dialysis, flush the catheter with saline solution and prime the catheter with heparinized saline solution. Connect the catheter to the saline solution. Flush the catheter with saline solution. Flush the catheter with saline solution.
 9. If you use the catheter as a primary catheter for dialysis, flush the catheter with saline solution and prime the catheter with heparinized saline solution. Connect the catheter to the saline solution. Flush the catheter with saline solution. Flush the catheter with saline solution.
- PROEIDOPHOSH:** Mην προωθήσετε ή αποσύρετε ενδοαγγειακό προϊόν αν συναντήσετε αντίσταση μέχρι να προσδιοριστεί ακτινοσκοπικά η αιτία της αντίστασης.
10. Αφαιρέστε τον καθετήρα ή/και το οδηγό σύρμα.
 11. Όταν η διαδικασία ολοκληρωθεί, αφαιρέστε το θηκάρι οδηγό CEREBASE DA χρησιμοποιώντας την καθιερωμένη τεχνική.
 12. Εάν επιθυμείτε εκ νέου πρόσβαση στα αγγεία με το ίδιο προϊόν, ξεπλύνετε και καθαρίστε τον εσωτερικό αυλό του προϊόντος με έγχυση. Επιθεωρήστε το προϊόν για τυχόν ζημιές.
- Προσοχή:** Μην χρησιμοποιείτε το προϊόν εάν παρατηρηθούν ζημιές ή ανωμαλίες.

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