MICRUSFRAME™
DELTA FILL™
DELTA XSFT™
MICROCOIL DELIVERY SYSTEMS

Medos International SÄRL
Chemin-Blanc 38, 2400
Le Locle, Switzerland

CE
2797

LCN 208535-001/E
© 2015-2019 Cerenovus
All rights reserved.
 Revised 11/19

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

MICRUSFRAME™

DELTAfill™

DELTAxSFT™

Microcoil Delivery Systems


STERILE R

Rx Only

System Description

The **MICRUSFRAME™**, **DELTAfill™**, and **DELTAxSFT™** Microcoil Delivery Systems consist of three components, a Microcoil System, a connecting cable, and an EnPOWER® Detachment Control Box (DCB). Each component is sold separately.

- As shown in Figure 1, the Microcoil System consists of an embolic microcoil attached to a device positioning unit (DPU) wire. It is covered by an introducer sheath system. The microcoil's shape can be spherical, complex, or helical. The DPU wire is a variable stiffness pushing system and has a radiopaque marker band located three (3) cm from its distal end. The introducer sheath system has three main components: an introducer tip, a translucent introducer body, and a re-sheathing tool.

- The connecting cable delivers the energy needed to detach the coil from the Microcoil System's detachment zone. The connecting cable is connected between the Microcoil System's hub connector on the DPU Wire and the output connector on the DCB. Its length is approximately 5–6 ft (1.5–1.8 m). The connecting cable may be one of two types: one with a remote detach button (the EnPOWER® Control Cable) catalog no. ECB000182-00, or one without a detach button (the standard connecting cable) catalog no. CCB000157-00.
- The EnPOWER Detachment Control Box (DCB), catalog no. DCB2000500, provides the energy necessary to allow for a thermo-mechanical detachment of the microcoil from the DPU wire. The EnPOWER Detachment Control Box works with the EnPOWER Control Cable or with the standard connecting cable. Both cables and the EnPOWER Detachment Control Box are available separately from the microcoil delivery system.

Indications for Use

MICRUSFRAME, DELTAfill, and DELTAxSFT Microcoil Delivery Systems are intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and are also intended for arterial and venous embolizations in the peripheral vasculature.

WARNINGS

- The Microcoil Delivery System should only be used as a complete system. A complete system requires a Detachment Control Box, a connecting cable, and Microcoil System. DO NOT SUBSTITUTE any components or devices from other manufacturers with the Cerenovus Detachment Control Box or injury to the patient or user could result.
- Explosion Hazard: The Detachment Control Box is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide or oxygen enriched atmospheres.**

PRECAUTIONS

- The Microcoil System shall be used only by physicians trained in interventional neuroradiology and in all aspects of the system.
- Do not attempt to sterilize any component of the Microcoil Delivery System.
- Verify the functionality of the Microcoil Delivery System before proceeding with microcoil placement. Refer to the section *Microcoil System Preparation* on how to verify the functionality.
- Store in a cool, dark, dry place.
- Do not use the device if the sterile barrier is compromised or the device is damaged.
- Do not use device that is past the "Use By" date.

Adverse Events

Possible adverse events include, but are not limited to, the following: hematoma at site of entry, vessel perforation, infection, emboli, hemorrhage, ischemia or vasospasm, neurological deficits including stroke, and possibly death.

How Supplied

2

The **MICRUSFRAME**, **DELTAfill**, and **DELTAxSFT** Microcoil System and cables are intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Discard the cables after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Cables are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Microcoils are permanent implants and are NOT re-usable. Discard the Microcoil System if it is retrieved from the patient without being implanted. Structural integrity and/or function may be impaired through reuse or cleaning. Microcoils are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Cerenovus will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner unit is not opened or damaged, the Microcoil System and the cables are sterile. The Microcoil System is nonpyrogenic.

Required Additional Materials

The following items are either supplied with the system or sold separately and must be on-hand prior to commencing a procedure.

- Femoral sheath
- Guiding catheter, 5 to 7 Fr
- Infusion microcatheter with 2 tip markers located 3 cm apart
- Guidewire compatible with the selected infusion microcatheter
- Three (3) continuous saline (or heparinized saline) flush set ups with pressure bags: one (1) flush system for the femoral sheath, one (1) for the guiding catheter, and one (1) for the microcatheter
- Two (2) rotating hemostatic valves (RHV)
- Three-way stopcock
- One-way stopcock
- IV Pole (not required if DCB will be placed in sterile sleeve and placed in sterile field)

Note: A backup EnPOWER Detachment Control Box (DCB2000500) is suggested for all procedures.

Instructions

Microcoil Size Selection

Microcoil selection is at the discretion of the physician. The appropriate microcoil size should be chosen based upon pre-embolization angiographic assessment of the diameter, height, and width of the aneurysm, as well as, the width of the aneurysm ostium (neck). In most cases, the initial microcoil implanted should be a three-dimensional spherical or complex shape. To minimize the potential of microcoil migration away from the aneurysm, the diameter of the first microcoil selected should not be less than the width of the aneurysm ostium (neck). Subsequent implanted microcoils may either be spherical, complex, or helical in shape. The selected coils typically will be of decreasing size and the physician may continue to implant microcoils until he or she determines that the aneurysm has been successfully treated.

Microcatheter Selection

Proper selection of the appropriately sized microcatheter is required to avoid damage to the Microcoil System and to minimize potential complications. Microcatheter selection is also determined by the physician and is predicated by the location of the aneurysm, patient safety, and physician preference. To ensure proper placement and detachment of the Microcoil, the selected microcatheter must have two (2) radiopaque tip markers three (3) cm apart.

The MICRUSFRAME, DELTAFILL, and DELTAXSFT Microcoil 10 Systems are compatible with microcatheters with inner lumen diameters ranging from 0.0165 to 0.017 in (0.419 to 0.432 mm).

The MICRUSFRAME and DELTAFILL Microcoil 14 and 18 Systems are compatible with microcatheters with inner lumen diameters ranging from 0.0165 to 0.021 in (0.419 to 0.533 mm).

Note: Fluoro Saver Markers are compatible with a 150 cm length microcatheter.

Continuous Saline Flush Set-Up

To achieve optimal performance of the Microcoil System, it is important that a continuous infusion of an appropriate flush solution be maintained. Figure 2 illustrates the connections necessary for the Microcoil Delivery System including a typical continuous saline flush setup with pressure bag for the catheter systems.

1. Attach a rotating hemostatic valve (RHV) to the hub of the guiding catheter.
2. Connect a three-way stopcock to the RHV sidearm, and attach a line to the stopcock for the continuous infusion of solution.
3. Carefully select an appropriately sized infusion microcatheter based on the size of the selected Microcoil System.

4. Insert the infusion microcatheter into the RHV connected to the guiding catheter hub.
5. Attach a second RHV to the hub of the infusion microcatheter.
6. Connect a one-way valve stopcock to the RHV's sidearm, and attach a flush line.
7. Adjust the hydrostatic pressure to 300 mm of Hg and maintain an open flush throughout the procedure.
8. Check to ensure that all fittings are secure and that no air is being introduced into the system during active flushing.


Fig. 2: Delivery System Set Up with Continuous Flush

Microcoil System Preparation

Verify the functionality of the Microcoil Delivery System before proceeding with microcoil placement. This needs to be done with the microcoil still in the hoop. To verify proper EnPOWER DCB and microcoil functionality a connecting cable and microcoil must be connected to the DCB unit. After verification of the DCB and connecting cable, turn off power to the DCB and disconnect the connecting cable from the microcoil until the microcoil is ready to be detached. Please refer to the Instructions for Use packaged with your EnPOWER DCB and your cable before proceeding.

Removing Microcoil System from Packaging Hoop


1. Grasp the base of the DPU's hub connector, and gently slide it completely out of the retaining clip. Be sure to keep the hub connector in line with the retaining clip until the entire connector is out.
2. Gently grasp the DPU wire hub and slowly pull the entire Microcoil System from the packing hoop. Do not bend the DPU wire as the device is pulled from the packing hoop as this may result in damage to the device.

Microcoil Placement

Introducing the Microcoil System


1. Loosen the main valve of the RHV attached to the infusion microcatheter hub.
2. Gently insert the introducer tip into the RHV until it can go no further and is in close alignment with the hub of the infusion microcatheter. (There may be a small space, as shown in Figure 3, between the tip and hub depending on the type of microcatheter used.) Gently tighten the RHV main valve around the introducer sheath to prevent back-flow of blood.

Visually inspect the alignment of the introducer tip and the microcatheter hub to ensure they have not slipped apart.

CAUTION: Do not fasten the RHV valve too tightly around the introducer sheath since excessive pressure may cause damage to the introducer sheath and/or the microcoil as it is advanced into the infusion microcatheter. Additionally, if the introducer tip and microcatheter hub are misaligned, damage may occur to the microcoil as it passes through this transition.


3. Grasp the distal end of the re-sheathing tool between your left thumb and forefinger. Grasp the clear tab near the proximal end of the introducer sheath body with the thumb and forefinger of your other hand. (See Figure 4). Gently pull the clear tab of the introducer sheath out and away from the re-sheathing tool at a 45-degree angle to unlock the microcoil. Continue to pull the tab until an additional 0.5 to 1.0 in (1.3 to 2.5 cm) of the translucent material is exposed.
4. Gently fold the translucent tab towards the distal end, and firmly grasp the distal end of the re-sheathing tool and the translucent tab between your thumb and forefinger, as shown in Figure 4.

5. While holding the distal end of the re-sheathing tool and translucent introducer sheath together between your thumb and forefinger, advance the DPU wire through the introducer sheath into the infusion microcatheter. **Note:** If unusual friction is noticed during advancement or retraction of the Microcoil System through the introducer, open the RHV main valve, and partially withdraw the distal end of the introducer to expose its tip within the RHV. Tighten the RHV main valve, and flush the Y-connector of the RHV with sterile saline and verify that fluid exits the slit in the clear portion of the introducer. Loosen the RHV main valve and fully re-insert the introducer tip into the infusion microcatheter hub. Gently tighten the RHV main valve around the introducer sheath to prevent back-flow of blood. Visually inspect the alignment of the introducer tip and the microcatheter hub to ensure they have not slipped apart.
6. As the microcoil passes through the introducer tip into the microcatheter hub, continuously verify that the introducer tip and microcatheter hub remain aligned. Continue to advance the DPU wire until its hub connector reaches the proximal end of the re-sheathing tool.
7. Return to the infusion microcatheter's RHV. Loosen the RHV and gently slide the introducer tip out from the RHV, over the DPU wire. Once a small section of the exposed DPU wire is visible, tightly grasp it with the thumb and forefinger of the same hand that is holding the RHV. Using the thumb and forefinger of your other hand, grasp the introducer tip, slowly slide

it away from the RHV over the DPU wire. Continue sliding the introducer tip until just before the tip reaches the distal end of the re-sheathing tool, leaving approximately 1 in (2-3 cm) of the unsheathed introducer sheath still visible.

8. Advance the microcoil system through the infusion microcatheter until you see the fluoro saver marker on the DPU wire approaching the RHV. The microcoil delivery system includes a group of five (5) fluoro saver markers on the shaft of the device. The markers are intended to indicate when the tip of the embolic coil is approaching the tip of the microcatheter when used with a 150 cm long microcatheter in conjunction with an RHV. When the distal-most marker reaches the proximal end of the RHV on the microcatheter, the tip of the coil is approaching the tip of the microcatheter, and fluoroscopy guidance should be used to guide further coil insertion. (See Figure 5). Continually verify the Microcoil System's progress and final position in the aneurysm using fluoroscopic guidance.

CAUTION: If unusual friction is noticed during advancement or retraction of the Microcoil System, verify the clear tab is unlocked and pulled out from the re-sheathing tool approximately 1 in (2-3 cm).

CAUTION: If unusual friction is noticed during advancement or retraction of the Microcoil System through the introducer, open the RHV main valve, and partially withdraw the distal end of the introducer to expose its tip within the RHV. Tighten the RHV main valve, and flush the Y-connector of the RHV with sterile saline and verify that fluid exits the slit in the clear portion of the introducer. After flushing, reinsert the introducer into the infusion catheter hub as described in "Microcoil Placement" section above.

CAUTION: If unusual friction is still noticed during advancement or retraction of the Microcoil System, verify flush lines are open and properly pressurized. Then slowly withdraw the entire Microcoil System and examine for damage. Replace it with a new Microcoil System. If friction still exists, withdraw and examine the delivery catheter system.

CAUTION: If the Microcoil System becomes immobile in the infusion microcatheter, apply a gentle push-pull motion to free it. If unsuccessful, remove both microcatheter and Microcoil System together as a unit and replace with new devices.

CAUTION: Do not attempt to use the Microcoil System as a guidewire if positioning of the microcatheter is lost during microcoil deployment.

CAUTION: If repositioning of the microcoil is necessary, carefully observe the motion of the microcoil in respect to the DPU wire while retracting the microcoil under fluoroscopy. If the microcoil movement is not one-to-one with the DPU wire, or if repositioning is difficult, the microcoil may have become stretched and could possibly break. Gently remove and discard the Microcoil System.

CAUTION: If the microcoil is positioned at a relative sharp angle to the microcatheter, a microcoil may stretch or break as it is being withdrawn. By repositioning the distal tip of the microcatheter at or slightly inside the ostium (neck) of the aneurysm, the microcoil may be more easily funneled back into the microcatheter.

Positioning the Microcoil

1. The proper alignment of the marker bands is illustrated in Figure 6. The infusion microcatheter has two distal tip marker bands located three (3) cm apart. The Microcoil System has a distal marker band three (3) cm from the detachment mechanism. To obtain proper positioning of the

Microcoil System for detachment under fluoroscopic guidance, align the radiopaque marker on the DPU wire just past the proximal marker band on the tip of the microcatheter.

CAUTION: Do not advance the DPU wire outside of the microcatheter as this may prevent detachment.

2. Once the desired microcoil placement has been achieved, gently tighten the RHV around the DPU wire to maintain its position.


Fig. 6: Proper Marker Band Alignment

Detaching the Microcoil

1. Re-verify the position of the microcoil within the aneurysm and the detachment zone position per Figure 6 via fluoroscopy. The Microcoil System is now ready for detachment.
2. Press the power button on the EnPOWER Detachment Control Box (DCB), if not already powered on.
3. Attach the connecting cable to the connector end of the DPU wire, ensuring that it is fully seated. If not already connected, attach the other end of the connecting cable to the output connector of the DCB.
4. Verify that the Microcoil Delivery System is fully connected and no faults are indicated on the DCB. If a fault exists, reseal all connections between the DPU, the DCB, and the connecting cable. If a fault still persists, replace the connecting cable. If this does not correct the error, replace the DCB. If the Microcoil Delivery System still continues to have a fault, retrieve the microcoil as described in the following section, re-sheathing the Microcoil System, and replace with a new Microcoil System.
5. If the system is free of faults, depress the **Detach** button on the EnPOWER DCB or the EnPOWER Control Cable. The **Detach Cycle** light next to the button will illuminate and an intermittent tone will sound for the duration of the detachment cycle. If the light and audible tone do not activate, replace the DCB.
6. After the light goes out and the tone stops, detachment of the microcoil from the DPU wire must be fluoroscopically verified by gently withdrawing the DPU wire back approximately one (1) mm. Observe if the microcoil has detached.
 - If the microcoil did not detach and no fault light is illuminated or flashing, repeat the steps above.
 - If a fault light is detected, the **System Ready** light is not illuminated, or there is no detachment after two detachment attempts, replace the DCB unit.
 - If detachment still does not occur, carefully withdraw the entire Microcoil System and redeploy a new Microcoil System.
7. After detaching the coil, remove the DPU wire from the microcatheter and discard.

Note: Do not disconnect the cable from either the DCB or DPU wire during a detachment cycle.

CAUTION: Always perform fluoroscopic verification of detachment prior to withdrawing the DPU wire fully. Failure to do so may lead to an embolic complication.
8. Repeat the above sequence for all additional microcoils until the procedure is complete.

Note: If it is determined that the microcoil will not be detached, simply disconnect it from the connecting cable and remove it from the microcatheter as outlined in the following section.

Re-sheathing The Microcoil System

When necessary, the Microcoil System can be reloaded into the introducer sheath using the re-sheathing tool.

1. Loosen the RHV. Using the fluoroscopic guidance, retract the microcoil from the aneurysm back into the microcatheter.
2. Locate the introducer tip. Replace the introducer tip back inside the RHV. Tighten the RHV. Grasp the introducer tip with your left hand and the re-sheathing tool with your right hand. Pull with your right hand while holding on the re-sheathing tool towards the connector. This will start the re-sheathing of the coil and the DPU wire.
3. With your right hand, grasp the connector and continue to pull the DPU wire out of the sheath until coil is completely inside the distal end of the introducer tip.


Fig. 7: Pulling Back the DPU Wire Hub Connector

CAUTION: Pulling the exposed microcoil through the RHV grommet may damage the coil.

CAUTION: Do not pull the DPU wire back too far since it may cause a softer section of the DPU wire to become exposed. An arm's length pull should re-sheath the coil.

4. Once the device is fully retracted, slide the re-sheathing tool over the DPU/introducer sheath body until the system is relocked.
5. Loosen RHV and remove device. If needed, the Microcoil System is now ready to be re-inserted.

Magnetic Resonance Imaging (MRI) Safety Information

Non-clinical testing demonstrated that the detachable coil is MR Conditional. A patient with this device can be scanned safely, immediately after implant, under the following conditions:

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 1000 Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the coil produced a temperature rise of +1.4 degrees centigrade during MRI performed for 15 minutes of scanning in a 3 Tesla MR system. Therefore, the MRI-related heating experiments for the detachable coil at 3 Tesla using a transmit/receive RF body coil at an MR system reported whole-body-averaged specific absorption rate (SAR) of 3.6-W/kg indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.4 degrees centigrade. MR image quality can be compromised if the area of interest is in the same area or relatively close to the position of the detachable coil. Optimization of MR imaging parameters is recommended.

Warranty

Cerenovus warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

® EnPOWER is a registered trademark of Cerenovus or affiliated companies.
™ MICRUSFRAME, DELTAFILL, and DELTAXSFT are trademarks of Cerenovus or affiliated companies.

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