EN Do not reuse!
ES No reutilizar!
PT Não reutilizar!
IT Non riutilizzare!
DE Nur zum Einmalgebrauch!
NL Voor éénmalig gebruik!
TR Tek Kullanmli!
EL Δεν να χρησιμοποιηθεί повторно!
PL Nie nadaje się do powtórnego stosowania!
DA Må ikke genbruges!
FR Ne pas réutiliser!
SR Samo za jednokratno korišćenje!
SK Nepoužívať opakovaně!
CZ Nepoužívat opakovaně!
ZH 请勿再次使用!

EN Date of Manufacture, Year/Month.
ES Fecha de fabricación, año/mes.
PT Data de fabricação, ano/mês.
IT Data di fabbricazione, Anno/Mese.
DE Herstellungsdatum (Jahr/Monat).
NL Productiedatum, jaar/maand.
TR Üretim Tarihi, Yıl/Ay.
EL Ημερομηνία κατασκευής - Έτος/Μήνας.
PL Data Produkcji, Rok/Miesiąc.
DA Produktionsdato, år/måned.
FR Date de fabrication, année/mois.
SR Datum proizvodnje, godina/mesec.
SK Dátum výroby, rok/mesiac.
CZ Datum výroby, rok/měsíc.
ZH 生产日期, 年/月

EN Expiry Date, Year/Month.
ES Fecha de caducidad, año/mes.
PT Data de expiração, ano/mês.
IT Data di Scadenza, Anno/Mese.
DE Haltbarkeitsdatum (Jahr/Monat).
NL Houdbaarheidsdatum (Jaar/Maand).
TR Son Kullanma Tarihi, Yıl/Ay.
EL Ημερομηνία λήξης, Έτος/Μήνας.
PL Data Ważności, Rok/Miesiąc.
DA Udløbsdato, år/måned.
FR Date d'expiration, année/mois.
SR Datum isteka roka trajanja, godina/mesec.
SK Dátum expirácie, rok/mesiac.
CZ Datum vypršení, rok/měsíc.
ZH 有效期, 年/月

EN Batch Number.
ES Código de lote.
PT Lote nº.
IT Numero di Lotto.
DE Chargennummer.
NL Lotnummer.
TR Seri No.
EL Αριθμός Παρτίδας.
PL Numer Partii.
DA Batchnummer.
FR Numéro de lot.
SR Broj šarže.
SK Kód dávky.
CZ Číslo šarže.
ZH 批号

EN Store between 5°C and 25°C.
ES Conservar a una temperatura de 5 °C a 25 °C.
PT Armazenar entre 5°C e 25°C.
IT Conservare ad una temperatura fra 5°C e 25°C.
DE Bei 5-25 °C lagern.
NL Tussen 5°C en 25°C bewaren.
TR 5°C ile 25°C arasında saklayınız.
EL Φυλάξτε σε θερμοκρασία μεταξύ 5°C και 25°C.
PL Przechowywać w temperaturze 5°C - 25°C.
DA Opbevares mellem 5 °C og 25 °C.
FR Conserver à une température comprise entre 5 °C et 25 °C.
SR Skladištili između 5°C i 25°C.
SK Skladujte medzi 5°C a 25°C.
CZ Uchovávejte při teplotě 5°C až 25°C.
ZH 置于 5° - 25°C 保存

EN Keep away from sunlight.
ES Mantener alejado de la luz del sol.
PT Manter longe da luz do sol.
IT Tenere lontano dalla luce diretta.
DE Vor Sonneneinstrahlung schützen.
NL Uit de zon bewaren.
TR Güneşten uzak tutunuz.
EL Να μην εκτίθεται στον ήλιο.
PL Chronić przed działaniem promieni słonecznych.
DA Må ikke udsættes for direkte sollys.
FR Conserver à l'abri de la lumière du soleil.
SR Čuvati zaštićeno od sunčeve svetlosti.
SK Chránite pred slnkom.
CZ Chránit před sluncem.
ZH 避免阳光直射

EN Protect from humidity.
ES Proteger de la humedad.
PT Proteger da umidade.
IT Proteggere dall'umidità.
DE Vor Feuchtigkeit schützen.
NL Droog bewaren.
TR Kuru yerde muhafaza ediniz.
EL Να προφυλάσσεται από υγρασία.
PL Chronić przed wilgocią.
DA Beskyttes mod fugt.
FR Conserver à l'abri de l'humidité.
SR Zaštiti od vlage.
SK Chránite pred vlhkom.
CZ Chránit před vlhkem.
ZH 避免受潮

EN Sterile R: Gamma Irradiation.
ES Esterilizado mediante irradiación: Irradiación gamma.
PT Steril R: irradiação gama.
IT Sterile R: Raggi Gamma.
DE Steril R: Gammastrahlen.
NL Steriel R: Gammastralen.
TR Steril R: Radyasyon.
EL Αποστειρωμένο.
PL Produkt jalowy: sterylizowany promieniowaniem gamma.
DA Sterile R: Gamma Irradiation.
FR Stérile R : Irradiation aux rayons gamma.
SR Sterilno R: gama zračenje.
SK Sterilizované žiarením.
CZ Sterilováno gama zářením.
ZH 无菌 R: 伽马辐照

CURACEL®

Absorbable Haemostat Oxidised Regenerated Cellulose

CURACEL^{HD}
(High Density)

Absorbable Haemostat Oxidised Regenerated Cellulose

CURACEL® FIBRILLAR

Absorbable Haemostat Oxidised Regenerated Cellulose

- EN Instructions for use
- ES Instrucciones de uso
- PT Instruções de utilização
- IT Istruzioni per l'uso
- DE Gebrauchsanleitung
- NL Gebruiksaanwijzing
- TR Kullanma Talimatı
- EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
- PL Wskazówki dotyczące stosowania opatrunku
- DA Brugsanvisning
- FR Instructions d'utilisation
- SR Uputstvo za korišćenje
- SK Návod na použitie
- CZ Návod k použití
- ZH 使用说明书

CE 1014

CuraMedical B.V.
Industrieweg 6B, 1566 JP ASSENDELFT, The Netherlands
Phone: +31 75 657 18 85 Fax: +31 75 657 11 84
www.curamedical.com

IMJ-CU-REV-05-21-08-2017

CURACEL[®]

Absorbable Haemostat Oxidised Regenerated Cellulose

Labelling:
Tyvek Cover

Labelling :
Individual (Sachet)
Pack / Sterile Pack

Labelling :
Box (Protective Cover) - Front

CURACEL^{HD} (High Density)

Absorbable Haemostat Oxidised Regenerated Cellulose

Labelling:
Tyvek Cover

Labelling :
Individual (Sachet)
Pack / Sterile Pack

Labelling :
Box (Protective Cover) - Front

CURACEL[®] FIBRILLAR

Absorbable Haemostat Oxidised Regenerated Cellulose

Labelling:
Tyvek Cover

Labelling :
Individual (Sachet)
Pack / Sterile Pack

Labelling :
Box (Protective Cover) - Front

Index

EN	
Instructions for use	1
ES	2
Instrucciones de uso	3
PT	4
Instruções de utilização	5
IT	6
Istruzioni per l'uso	7
DE	8
Gebrauchsanleitung	9
NL	10
Gebruiksaanwijzing	11
TR	12
Kullanma Talimatı	13
EL	14
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	15
PL	16
Wskazówki dotyczące stosowania opatrunku	17
DA	18
Brugsanvisning	19
FR	20
Instructions d'utilisation	21
SR	22
Uputstvo za korišćenje	23
SK	24
Návod na použitie	25
CZ	26
Návod k použití	27
ZH	28
使用说明书	29
	30

CURACEL®

CURACEL® HD®
(High Density)

CURACEL® FIBRILLAR

100 % ABSORBABLE HAEMOSTAT
(oxidized regenerated cellulose) U.S.P

EN

Instructions for use

DESCRIPTION

Curacel® Absorbable Haemostat is a sterile product made of oxidized regenerated cellulose (polyoxan hydro glucuronic acid) available in the flat form with textile structure or in fibrillose form. Curacel® Haemostat can be cut into the desired dimension. Its structure is stable and should be stored at controlled room temperature. The fabric is pale yellow and has a faint, caramel-like aroma. A slight discoloration may occur with time, but this does not affect its performance.

ACTIONS

Curacel® Haemostat is applied dry. It can be laid on, held against or wrapped around a bleeding surface. As special features, Curacel High Density can be sutured, and Curacel Fibrillar can be peeled off in layers for optimal dosage. Action mechanism of Curacel® is independent from blood coagulation mechanism of the body. When it contacts with blood, Curacel® Haemostat forms a brownish or black gelatinous mass which aids clot formation in 3-4 minutes. This gelatinous mass acts as a physical matrix to which platelets can adhere. With platelet aggregation and formation of platelet-fibrin plug, haemostasis occurs. When Curacel® Haemostat is used properly in minimal amounts, it is absorbed from the sites of implantation without tissue reaction. Absorption depends on; Haemostat amount used, tissue type and blood saturation degree. It is completely absorbed within 7-14 days. Curacel® Haemostat has been shown to be bactericidal in vitro against at least 40 types of (pathogenic) Gr(+) (Gram-positive) and Gr (-) (Gram-negative) microorganisms. These microorganisms include;

Staphylococcus aureus
Methicillin Resistant Staphylococcus epidermidis (MRSE)
Methicillin Resistant Staphylococcus aureus (MRSA)
Staphylococcus epidermidis
Micrococcus luteus
Vancomycin Resistant Enterococcus sp (VRE)
Penicillin Resistant Streptococcus pneumoniae (PRSP)
Enterococcus faecalis
Enterococcus faecium
Enterococcus sp
Streptococcus pyogenes Group A
Streptococcus pyogenes Group B
Streptococcus salivarius
Corynebacterium xerosis
Escherichia coli
Salmonella typhimurium

Proteus mirabilis
Proteus vulgaris
Shigella boydii
Shigella flexnerii
Shigella dysenteriae
Salmonella enteritidis
Bacteroides fragilis
Lactobacillus casei
Lactobacillus sp.
Pseudomonas aeruginosa
Pseudomonas stutzeri
Enterobacter cloacae

Klebsiella pneumoniae
Klebsiella aerogenes
Bacillus subtilis
Branhamella catarrhalis
Serratia marcescens
Clostridium perfringens
Clostridium tetani
Listeria monocytogenes
Bordetella pertussis
MOT (Mycobacterium Other than Tuberculosis)
Mycobacterium phlei
Streptococcus faecalis

INDICATIONS

Curacel® Haemostat is designed to arrest capillary bleeding and bleeding from parenchymatous organs and resection areas at surgical interventions. It shall be used to stop bleeding when other haemostatic procedures are ineffective or impractical. It is suitable for use in general surgery and digestive surgery, neurosurgery (**especially cerebral operations**), plastic surgery, orthopaedic, gynaecology, urology, stomatology, traumatology, and many other branches of surgery. Curacel® Haemostat can be applied into cavities (after extirpation of tumours) as well as endoscopic interventions or dental praxis.

CONTRAINDICATIONS

Curacel® Haemostat should not be used for packing. If you have to use in this manner, it must be removed after haemostasis is achieved. Curacel® Haemostat should not be used for implantation in bone defects, such as fractures, since it may interfere with callus formation and may cause cyst formation. Curacel® Haemostat should not be used to control haemorrhage from large arteries. Since body fluids, except whole blood, such as serum, don't react with Curacel® Haemostat, it should not be used to produce haemostatic effect on non-haemorrhagic serous oozing surface.

WARNINGS

Curacel® Haemostat should not be re-sterilized, not even by gamma-irradiation, because product conformity cannot be guaranteed. Closing Curacel® Haemostat in a contaminated wound without drainage may lead to complications and should be avoided. The haemostatic effects of Curacel® Haemostat are greater when it is applied dry; it should not be moistened with water or saline prior to use. Impregnation of Curacel® Haemostat with anti-infective agents, buffering materials or haemostatic substances should not be done. Since thrombin activity is destroyed due to the low pH of the product, haemostatic effect of Curacel® Haemostat is not enhanced by addition of thrombin. Curacel® Haemostat may be left in tissue when necessary, but it is advisable to remove it once haemostasis is achieved. It must always be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm regardless of the type of surgical procedure, because Curacel® Haemostat may exert pressure resulting in paralysis and/or nerve damage by swelling. Re-motion of Curacel® Haemostat could possibly occur by means such as repacking, further intraoperative

manipulation, lavage, etc. There have been reports that Oxidized (Regenerated) Cellulose Haemostat, when left in the patient after closure during procedures such as lobectomy, laminectomy and repair of a frontal skull fracture and lacerated lobe, migrated from the site of application into foramina in bone around the spinal cord resulting in paralysis and, in another case, migrated to the left orbit of the eye causing blindness. While these reports cannot be confirmed, physicians must act cautiously, regardless of the type of surgical procedure, and consider the advisability of removing Curacel® Haemostat once haemostasis is achieved. Although Curacel® Haemostat is bactericidal against a wide range of pathogenic microorganisms, it should not be used as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections. Curacel® Haemostat should not be used as an adhesion prevention product.

PRECAUTIONS

Use only minimal amounts of Curacel® Haemostat in appropriate size needed to obtain haemostasis, holding it firmly in place until bleeding stops. Remove any excess before closure in order to facilitate absorption and minimize the possibility of foreign body reaction. In urological procedures, minimal amounts of Curacel® Haemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product. Since absorption of Curacel® Haemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals. If Curacel® Haemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped. Precautions should be taken in otorhinolaryngology surgery to assure that none of the material is aspirated by the patient. (Examples: controlling haemorrhage after tonsillectomy and controlling epistaxis). Care should be taken not to apply Curacel® Haemostat tightly when it is used as a wrap during vascular surgery.

ADVERSE REACTIONS

"Encapsulation" of fluid and foreign body reactions has been reported for Oxidized (Regenerated) Cellulose Haemostat. There have been reports of stenotic effect when Oxidized (Regenerated) Cellulose Haemostat has been applied as a wrap during vascular surgery. Although it has not been established that the stenosis was directly related to the use of Oxidized (Regenerated) Cellulose Haemostat, it is important to be cautious and avoid applying the material tightly as a wrapping. Paralysis and nerve damage have been reported when Oxidized (Regenerated) Cellulose Haemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures. Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when Oxidized (Regenerated) Cellulose Haemostat was placed in the anterior cranial fossa. Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy has been reported. There has been one report of a blocked ureter after kidney resection, in which post-operative catheterization was required. Occasional reports of "burning" and "stinging" sensation and sneezing when Oxidized (Regenerated) Cellulose Haemostat has been used as packing in epistaxis, are believed to be due to the low pH of the product. Burning has been reported when Oxidized (Regenerated) Cellulose Haemostat was applied after nasal polyp removal and after haemorrhoidectomy. Headache, burning, stinging, and sneezing in epistaxis control and other rhinological procedures, and stinging when Oxidized (Regenerated) Cellulose Haemostat was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) have been reported as well.

DOSAGE AND ADMINISTRATION

Sterile technique should be observed in removing Curacel® Haemostat from its sterile container. Minimal amount of Curacel® Haemostat in appropriate size is laid on the bleeding site or held firmly against the tissues until haemostasis is achieved. Required amount of Curacel® Haemostat depends on the nature and intensity of the haemorrhage to be stopped. The haemostatic effect of Curacel® Haemostat is optimal when used dry. Moistening the material with water or physiological saline solution prior to use is not recommended.

CAUTION

Since Curacel® Haemostat cannot be re-sterilized, opened and unused Curacel® Haemostat should be discarded.

NOTES AND INFORMATION ON SHELF-LIFE

Curacel® Haemostat should be stored dry under controlled conditions (5°C - 25°C) in the original packaging, and protected from direct sunlight. Shelf life of the product is 3 years. The expiry date of Curacel® Haemostat is printed on the pack. Do not use Curacel® Haemostat after this date.

PRESENTATION AND CONTENTS

Curacel® Haemostat is available sterile, as an original pack in different sizes.

Curacel® Haemostat strips are sterile-packed individually.

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