SUTMEDICAL NON ABSORBABLE SURGICAL SUTURES®

PRESENTATION

12 (twelve) Unit box in an individual sterile package. **SUTMEDICAL®** sutures are available in a variety of sizes, lengths and needles.

They are available in the following types of needles:

- **Sharp:** those needles have unless two opposite sharp edges. They are sharpened to cut through hard tissues which are difficult to penetrate. Due to its sharp edge, it is necessary to take care with some tissues to avoid cutting more than the necessary limit.

- **Round:** these needles do not have a cutting edge, they are the most used to suture tissues which are easy to penetrate as peritoneum, abdominal viscera, cardiovascular tissue and neurological tissues.

- **Cutting Reversal:** these needles were created specifically for tissues which are difficult to penetrate as the skin or the oral mucus. This needle is so sharp as the conventional cutting needle but its design is different, since the third cutting edge is located in the external convex curve of the needle. This feature provides more resistance to this needle.

- **Spatulate:** these needles does not have a Sharp point, since due to their structure they are used in proceedings in which it is necessary to generate the minor trauma to the tissue, as in ophthalmological proceedings.

DESCRIPTION

- **Nylon**
Non absorbable monofilament surgical suture, synthetic, sterile which is resistant to abrasion. The nylon fibers are characterized for its high tension force and minimum tissue reaction, which is followed by a gradual encapsulation of the suture in the connective fibrous tissue. They are resistant to organic enzymes and they avoid bacterial growing. They do not present absorption by tissue enzymes and they are degraded. In vivo through hydrolysis at an approximate rate of 15-20% per year. The suture is tied with black color to remark the tissue visibility.

- **Silk**
Non absorbable multifilament surgical suture, coated. It is a natural material of animal origin obtained from silk worm. It does not present degrading or absorption by tissue enzymes, it is black color in order to remark the visibility in the tissue, it loses the greatest part of its tensile force per year and it has a high limit of elasticity. Due to it, one of the most notorious characteristics is the easiness to manage it.

- **Polypropylene**
Surgical suture sterile monofilament. It is a synthetic fiber, its chemical structure permits a greater flexibility and management. Polypropylene is an inert material in the tissue, it manages a high tensile force which is retained by periods up to two years. It has a low density, high hardness and resistance to abrasion, high rigidity, good heat resistance, excellent chemical resistance and versatility. It does not present degrading or absorption by tissue enzymes. Besides, it is not adhered to tissue. Therefore, it is efficient as an easy removal suture.

- **Polyester**
Multifilament braided surgical suture, sterile, recovered, White and green color to remark the visibility in the tissue and it has a high limit of elasticity. It has a high resistance to weakening and corrosion; it has a good chemical and thermal resistance; it presents a high rigidity, hardness, high tensile force and low humidity absorption. Its uniform surface minimizes the trauma. It does not present degrading or absorption by tissue enzymes.

- **PVDF (Polivinil denfluoruro)**
Monofilament sterile surgical suture, with a good resistance to abrasion. It has a high flexibility and a greater resistance to pressure than other materials. It also has a low friction coefficient and a high dimensional stability. PVDF fibers are characterized by their great biostability, high tension force and the minimum tissue or hystic reaction, which is followed by a gradual encapsulation of the suture in the connective fibrous tissue. They do not present degrading or absorption by tissue enzymes. Besides, they are easy to manage, which provides a joint process which is tighter with less knots. The suture is blue color to remark visibility and the same as polypropylene, it maintains knots in a better way in comparison to the majority of other synthetic materials of monofilament.

SUTMEDICAL® sutures meet all requirements established by the United States Pharmacopeia (U.S.P) and NTC 2267 standards for non absorbable surgical sutures which are described as follows:

USP Designation	Metric size (caliber No.)	Diameter in mm		Resistance to tension in knots N (kg-f)	
		Average values		Average minimum values	Individual values in min.
		Minim	Maxim		
9-0	0,4	0,04	0,049		
8-0	0,5	0,05	0,069	0,441 (0,045)	0,245 (0,025)
7-0	0,7	0,07	0,099	0,68 (0,07)	0,54 (0,055)
6-0	1	0,1	0,149	1,77 (0,18)	0,98 (0,10)
5-0	1,5	0,15	0,199	3,73 (0,38)	1,96 (0,20)
4-0	2	0,2	0,249	7,55 (0,77)	3,92 (0,40)
3-0	3	0,3	0,349	12,25 (1,25)	6,67 (0,68)
2-0	3,5	0,35	0,399	19,61 (2,00)	10,2 (1,04)
0	4	0,4	0,499	27,16 (2,77)	14,22 (1,45)
1	5	0,5	0,599	37,27 (3,80)	19,12 (1,95)
2	6	0,6	0,699	44,23 (4,51)	23,6 (2,40)
3	7	0,7	0,799	57,89 (5,90)	29,32 (2,99)
4	8	0,8	0,899	68,65 (7,00)	34,23 (3,49)

INSTRUCTIONS

Sutures are sterilized with Ethylene Oxide.

Polypropylene, PVDF and Nylon. Being monofilament sutures, these are made of only one material thread and due to their simplified structure, these sutures provide a greater easiness in the passing through the tissue.

These sutures are indicated to approximation of soft tissue in general and/or ties, including the use of general cardiovascular, ophthalmological, orthopedic and neurological proceedings. They are also used when a minimum tissue reaction is desired, as in contaminated or infected wounds, in which fistula formation and suture extrusion can be minimized.

Polyester and silk sutures, due to those are multifilament, are formed by several filaments, threads or strands; both plaining. This feature provides a greater force of tension and flexibility.

These sutures are indicated to approximation of soft tissue in general and/or tie, including the use in general, cardiovascular, gastrointestinal, odontologic, ophthalmologic, neurologic and plastic surgery proceedings. These sutures are commonly used in mucus and Cutaneous regions with folds or constant friction.

INSTRUCTIONS OF USE

Sutures should be used as it is required in the surgical proceeding. Size and diameter of suture and the needle used for each type of surgical operation, are defined by medical criteria.

Recommendations for the management of sutures:

1. Read the expiration date with care.
2. Only open the sutures which are necessary for the proceeding.
3. Straight sutures with a soft traction with the purpose of removing the package memory.
4. Do not press or tight the suture threads with surgical instruments, so friction which weakens suture integrity should be avoided, exerting an excessive tension on the suture which could break and/or cut the tissue, should be also avoided.
5. Do not pull the needle.
6. To avoid damaging needles point, the needle should be tightened in an area of one third (1/3) or one half (1/2) of the distance, from the end to the tip proximal medium point. Forcing needles should cause they lose their force and they are less resistant to folding or breaking.
7. Maintaining dryness of silk.
8. The finished knot should be strong to eliminate sliding; it is recommended to make a knot which is the smallest as possible and ends should be the shorter as possible in order to minimize the tissue reaction of a foreign matter; traction in a thread end should be maintained after making the first knot in order to avoid that it is loosened and making the final knot the most horizontal as possible.
9. Taking care of managing surgical needles in order to avoid a puncture.
10. Discard used needles in a sharp container

CONTRAINDICATIONS

No known contraindications.

WARNINGS

The user should be familiarized with the surgical proceedings and the techniques related to the use of non absorbable sutures, before using SUTMEDICAL for wounds closing. Product should not be re-sterilized. Discard opened packages and non used sutures. Transporting and storing the product at a maximum temperature of 25C grades, in a fresh and dry place protected from humidity, sun light and other heat sources. Avoid use after the expiration date.

PRECAUTIONS

Whenever tissue integrity is broken due to an accident or dissection, the patient is vulnerable to an infection and its complications. Although the surgical equipment follows the adequate proceeding rigorously, complications in some patients could happen and they delay recovery. The two greater problems which could be faced are infection and separation of the wound edges.

In the moment of using SUTMEDICAL as a surgical suture, there are factors which have an impact on the proceeding dealing and the selection of materials to close the wound. The necessary factors to be taken into account in the patient and which could affect the surgical suture purpose are age, weight, nutritional conditions, dehydrating, inadequate blood supply to the wound place, immunologic answer, the presence of chronic diseases, neoplasias, weakening injuries or located infection.

In order to avoid the risk of dehiscence, a full knowledge of suture and material techniques should exist, since depending on the type of suture and the management of the same, this complication can be avoided.

ADVERSE REACTIONS

As it happens with the use of any medical device, surgical sutures can have adverse effects related to their use, which is not only and proper of SUTMEDICAL, but it is intrinsic of all materials used in surgical sutures around the world. Adverse reactions could include wound dehiscence, formation of bilestone, or urinary stones, infections, minimum acute inflammatory tissue reaction and local transitory irritation in the wound place. The broken needles could cause a prolonged or additional surgery, or also residual foreign bodies.

SYMBOLS USED

-  Sharp needle
-  Reverse cutting needle
-  Spatulate needle
-  Round needle
-  Lot number
-  Product manufacturing date
-  Expiration date
-  Read use instructions
-  Not reuse
-  Care of environment conditions
-  Sterilization with Ethylene Oxide

BIBLIOGRAPHY:

- Moredó K, Armas B, Segura L, *et al.* Materiales de sutura quirúrgica. *AMC* 2009; 13:5 (Surgical suture materials)
- Kronenthal RL. Intraocular degradation of nonabsorbable sutures. *Amer Intra-Ocular Implant Soc J* 1977; 3:222-28
- Cagici CA, Cakmak O, Bal N, *et al.* Effects of Different Suture Materials on Cartilage Reshaping. *Arch Facial Plast Surg* 2008; 10(2):124-29
- Silva RF, Eraso MF, Villalobos ME, *et al.* Comparison between polypropylene and polyamide as suture elements for the canines skin. *Vet Zootec* 2007; 1(1):36-40
- Israelsson L, Jonsson T. Suture length to wound length ratio and healing of midline laparotomy incisions. *The British Journal of Surgery* 1993; 80(10):1284-86
- Rath A, Attali P, Dumas J, *et al.* The abdominal linea alba: an anatomologic and biomechanical study. *Surgical and Radiologic Anatomy* 1996; 18(4):281-88
- Postleimiarr RW. Long-Term Comparative Study of Nonabsorbable Sutures. *Annals of Surgery* 1970; 171(6):892-97
- Bellón JM, Rodríguez M, Serrano N, *et al.* Polipropileno frente a polidioxanona: igual eficacia biomecánica en cierres de la línea media. *Cir Esp* 2005;78(6):377-81 (Polypropylene in comparison to polidioxanone: equal biomechanical efficiency in medium line closings)
- Dineen P. The effect of suture material in the development of vascular infection. *Vasc Surg* 1977; 11(1):29-32,
- BHerrman JB. Tensile strength and knot security of surgical suture materials. *Am Surg* 1971; 37(4):209-17
- Rodríguez JI, Cazador AC, Ruiz B, *et al.* Factores de riesgo de dehiscencia aguda de la pared abdominal tras laparotomía en adultos. *Cir Esp* 2005;77(5):280-6 (Acute dehiscence risk factors of the abdominal wall after laparotomy in adults)

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