PROOFREADING APPROVAL/RELEASED ARTWORK		Check as applicable	PROOFREADING APPROVAL/RELEASED ARTWORK		Check as applicable
Printed Name		☐ Excluding dimensions ☐ Excluding bar code dimensions ☐ U.S. English ☐ Other Language:	Printed Name		☐ Excluding dimensions ☐ Excluding bar code dimensions ☐ U.S. English ☐ Other Language:
Signature			Signature		
Date			Date		

7N8399

Baxter

ONE-LINK
Needle-free IV Connector
with Neutral Fluid Displacement
Power Injectable (325 psi / 2241 kPa)

0.08 mL
Sterile, Nonpyrogenic.

Do not use if package has been opened or damaged.

Indications for Use:

The Baxter ONE-LINK Needle-free IV connector is intended for single patient use with a vascular access device for the administration of drugs and solutions without needles, thus eliminating the potential for needle-stick injuries during use. This device is an in-line injection site which can be connected to standard male Luer adapters (e.g., syringe or sets) for continuous or intermittent fluid administration or withdrawal of fluids. This device may be used with low pressure power injectors.

Directions: Use aseptic technique.

- When using ONE-LINK connector, grasp finger grip to stabilize.
- Swab** Luer activated surface with preferred antiseptic **prior to first use and before every subsequent connection.**
- Prime to purge air. Remove protector.
- For each access, firmly push male Luer of syringe or administration set directly against ONE-LINK connector's Luer activated surface and rotate until connection is secure. Engage Luer locking features, if applicable.
- Hold the ONE-LINK connector securely while detaching syringe or administration set.
- Connect to female Luer of IV catheter or extension set. Do not over-tighten.
- Clamp extension set or non-valved catheter when not in use.

Cautions:

- To avoid the risks associated with inadvertent disconnection, clamp the extension set or non-valved catheter when not in use.
- Prior to power injection, ensure all parts of the I.V. system are compatible. Replace if a pressure over 325 psi / 2241 kPa is applied to the connector.
- Do not allow air to be trapped in device.
- Use of Luer lock connection is recommended. If Luer slip connection is used, insert into ONE-LINK connector using a firm push and twist motion. Do not leave Luer slip connection unattended.**
- Do not use ONE-LINK connector with needles or cannulas. Attempting such access will render the product damaged; replace immediately.**
- If ONE-LINK connector cannot be cleared of blood after blood infusion/sampling, replace immediately.
- If blood remains on ONE-LINK connector's Luer activated surface after sampling, swab surface.
- Trace lines before connection. Do not connect any compressed gas devices to vascular injection sites.
- Replace the ONE-LINK connector if exposed to gross contamination.
- Replace the ONE-LINK connector if damaged.
- USA: Rx Only.
- Single use only. Do not resterilize. Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.

Notes:

- This product is not made with natural rubber latex.
- This product does not contain PVC or DEHP.
- This product does not contain metal components.
- This product can be used to administer lipid-based solutions.
- Flush ONE-LINK connector after injection to prevent inadvertent mixing of incompatible medications/fluids.
- Flush ONE-LINK connector with 10 mL or more after blood infusion/sampling. If 10 mL flush cannot be performed after blood infusion/sampling, replace ONE-LINK connector immediately.
- Luer lock blood collection devices, syringes and accessories are recommended for blood sampling.
- Always protect the ONE-LINK connector with an impermeable cover during showering or bathing.
- If intermittently disconnecting the administration set from the ONE-LINK connector, immediately cover the male Luer of the administration set with a sterile replacement protector.
- Replace per institutional protocol. The ONE-LINK Needle-free IV connector can be used for up to 200 actuations and over a period of 7 days; replace device whichever comes first.
- Volume is approximate.
- This device may be used with low pressure power injectors having a maximum pressure of 325 psi/2241 kPa and a maximum flow rate of 10 mL/s.


For Product Information 1-800-933-0303

Baxter Healthcare Corporation
Deerfield, IL 60015 USA

Baxter and One-Link are trademarks of Baxter International Inc.
07-19-77-205 Rev 2016-07-01

071977205

Symbols Glossary

Symbol	Title of Symbol/Definition/Reference	Symbol	Title of Symbol/Definition/Reference
	Consult instructions for use Indicates the need for the user to consult the instructions for use. Ref 5.4.3 in ISO 15223-1:2012		Do not re-use Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. Ref 5.4.2 in ISO 15223-1:2012
	Do not use if package is damaged Indicates a medical device that should not be used if the package has been damaged or opened. Ref 5.2.8 in ISO 15223-1:2012		Sterile Indicates a medical device that has been subjected to a sterilization process. Ref 5.2.1 in ISO 15223-1:2012
	Do not resterilize Indicates a medical device that is not to be resterilized. Ref 5.2.6 in ISO 15223-1:2012		Non-pyrogenic Indicates a medical device that is non-pyrogenic. Ref 5.6.3 in ISO 15223-1:2012

Note: Symbols are contained in ISO 15223-1:2012 Medical devices — Symbols to be used with medical device labels, labeling and information to be supplied — Part 1: General requirements

PROOFREADING APPROVAL/RELEASED ARTWORK	Check as applicable	PROOFREADING APPROVAL/RELEASED ARTWORK	Check as applicable
Printed Name	☐ Excluding dimensions ☐ Excluding bar code dimensions ☐ U.S. English ☐ Other Language:	Printed Name	☐ Excluding dimensions ☐ Excluding bar code dimensions ☐ U.S. English ☐ Other Language:
Signature	Date	Signature	Date

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Baxter et One-Link sont des marques de commerce de Baxter International Inc.