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
Certificado de calidad Quality certification Bionova® BT91

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ES

Biological Indicator

For Hydrogen Peroxide sterilization

Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores deposited on the inner surface of the plastic tube base. It also has a growth indicator media of purple color contained in the glass ampoule.

Product description

Bionova® BT91 Biological Indicator is specifically designed for the monitoring of Plasma or Vaporized Hydrogen Peroxide sterilization processes. If sterilization process was not successful, the indicator media will change from purple to yellow after incubation at 60±2 °C, thus indicating the presence of live *Geobacillus stearothermophilus* spores. If the sterilization process was successful, the indicator media will remain purple after incubation. The final readout should be made after 24 hours of incubation at 60±2 °C.

Warning!

Do not use Bionova® BT91 Biological Indicator for monitoring Steam sterilization cycle, Ethylene Oxide, Dry Heat, Formaldehyde or other sterilization processes. Do not re-use the biological indicators.

Storage

Store it in a dark place at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Directions for use

1. Identify the Bionova® BT91 Indicator by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the indicator label.
2. Pack the biological indicator along with materials to be sterilized in an appropriated package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for sterilizing agent.
3. Sterilize as usual.
4. After the sterilization process has finished, remove the biological indicator from sterilization package for incubation.
5. Check the chemical indicator on the label of biological indicator. A color change to green confirms that the biological indicator has been exposed to Hydrogen Peroxide. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.
6. Crush the glass ampule contained in the biological indicator and incubate at 60±2 °C. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control each time a processed biological indicator is incubated. The positive control ensure that correct incubation conditions were met.
7. Incubate the processed biological indicator and the indicator used a positive control for a maximum of 24 hours at 60±2 °C.

A color change from purple to yellow of the growth indicator media means that sterilization process failure has occurred. If after 24 hours there is no color change in the processed indicators, a final negative result is made (the sterilization process was acceptable). The positive control indicator should show a purple to yellow color change for the result to be valid.

NOTE: To extend the incubation time to more than 24 hours, we recommend using a humidified environment to prevent the complete evaporation of the culture medium contained in the indicator.

Record the results and discard immediately as it is shown below.

WARNING! Do not use the sterilizer until the biological indicator test results are negative (process indicator remains purple).

Disposal

Dispose of used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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ATCC is a registered trademark of American Type Culture Collection