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Firma válida

Digitally signed by ELKIN
HERNAN OTALVARO
CIFUENTES
Date: 2016.04.26
10:35:11 -0500
Reason: Firma
Location: Bogota, CO

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

 **LG Life Sciences**
129, Seokam-ro, Iksan-si,
Jeollabuk-do, Korea

LOT  

【Intended use】 For use as a symptomatic treatment for osteoarthritis of weight bearing joints including knee, shoulder, hip and ankle.

【Composition】

Active ingredient: BDDE-crosslinked sodium hyaluronate gel - 3.0g

(BDDE-crosslinked sodium hyaluronate - 60mg)

【Contents】 3.0mL/syringe X 1 syringe

【Storage】 Store at 1~30°C free from light, in a hermetic container.

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Cross-linked hyaluronic acid

Description

Hyruan ONE is a milestone viscosupplementation treatment used to relieve pain associated with osteoarthritis of the knee, shoulder, hip and ankle. Hyruan ONE is a clear, colorless and viscous gel containing 60mg BDDE(1,4-butanediol diglycidyl ether)-crosslinked sodium hyaluronate in a 3mL pre-filled syringe with rubber stoppers on the both ends.

Contents

Each Hyruan ONE (volume of 3.0mL/syringe) contains:
Active ingredient: BDDE-crosslinked sodium hyaluronate gel - 3.0g
(BDDE-crosslinked sodium hyaluronate - 60mg)

Intended Use

For use as a symptomatic treatment for osteoarthritis of weight bearing joints including knee, shoulder, hip and ankle.

Dosage

Hyruan ONE is a single dose for single injection, and should only be injected once per treatment course. The recommended dose is 3mL for knee, shoulder or hip joint and 1-2mL for ankle joint.

Mode of Action

Naturally-occurring hyaluronic acid in the various part of the body acts in the joints as a shock absorber, protects cartilage, and relieves pain associated with osteoarthritis.

Contraindications

1) Hyruan ONE should not be injected:

- To patients who are known to be sensitive to the product or its ingredient(s)
- To patients with an infection or severe inflammation at joint cavity
- To patients with a disease or infection of skin near the injection site

2) Hyruan ONE should be injected with cautions:

- In patients with hypersensitivity to other substance
- In patients with liver disease or prior history of the liver disease

Adverse Events

1) Serious adverse events
Shock: since symptoms of shock (frequency unknown) may occur, careful observation should be made. In case abnormalities are noticed, administration should be discontinued and proper measures should be taken.
2) Reported adverse events in clinical trials
In a pivotal clinical trial of Hyruan ONE in patients with the knee osteoarthritis (a total of 285 subjects), the occurrence rate of local reaction at injection site after injection into articular cavity was 48.9% (68/139 subjects) in the study group and 49.3% (72/146 subjects) in the control group. The reported adverse events are shown in Table 1. The serious adverse events were reported in the decreasing order of pain (7.2% of the study group, 6.2% of the control group, redness (5.0% of the study group, 2.1% of the control group) and so on. Adverse events lasting over a week were pain 5.3%, redness 1.4%, swelling and warmth each accounting for 1.0%, yet all of them completely resolved in two weeks without any special treatment.
Table 1. Reported local reaction at injection area in the clinical trial

Adverse events	Hyruan ONE N=139(%)	Hyruan Plus*(Control) N=146(%)
Pain	60(43.2)	55(37.7)
Redness	27(19.4)	31(21.2)
Swelling	17(12.2)	21(14.4)
Warmth	33(23.7)	30(20.5)

* Hyruan Plus, requiring three weekly injections, is a viscosupplement for osteoarthritis developed by LG Life Sciences, Ltd.

3) Of the patients with the knee osteoarthritis (285 subjects), the occurrence rate of the adverse events excluding ones occurred at the injection site is 34.5% (48/139, 73 cases) for the study group and 28.8% (42/146, 63 cases) for the control group. Most of them were mild to moderate. Table 2 indicates the adverse events occurred in more than 1% of the study group.

Table 2. Adverse events reported in more than 1% of the subjects treated with Hyruan ONE in the clinical trial

Adverse events	Hyruan ONE N=139(%)	Hyruan Plus*(Control) N=146(%)
Nasopharyngitis	7(5)	3(2.1)
Pain	4(2.9)	1(0.7)
Warmth	3(2.2)	0(0)
Upper respiratory tract infection	3(2.2)	0(0)
Pain in extremity	3(2.2)	6(4.1)
Cystitis	2(1.4)	2(1.4)
Dyspepsia	2(1.4)	1(0.7)
Parathesia	2(1.4)	0(0)
Redness	2(1.4)	0(0)
Joint swelling	2(1.4)	0(0)
Musculoskeletal pain	2(1.4)	0(0)
Plantar fasciitis	2(1.4)	0(0)

125*190mm(back)

Precautions

1) General precaution:

- Administration of Hyruan ONE to severely inflamed joints caused by deformative osteoarthritis can lead to exacerbation of the local inflammation symptom. Thus, Hyruan ONE is desired to be given after the removal of existing inflammation symptoms.
- Administration of Hyruan ONE may occasionally cause local pain or swelling. Thus, patients should be informed to avoid strenuous exercise or action leading joint pain of knee up to 48 hours after administration, and actions such as local relaxation should be guided after injection.
- Leakage of Hyruan ONE other than articular cavity might cause pain. Thus, Hyruan ONE should be accurately injected into articular cavity.

2) Precautions in use:

- The product should be administered by trained doctors.
- Since Hyruan ONE is injected directly into the joints, administration should be performed under intact sterilization status.
- In case of retention of articular fluid, the fluid should be removed by puncturing prior to administration of Hyruan ONE.
- Intravascular injection, extra-articular injection or injection in the synovial tissues should be avoided.
- It is desirable to administer the product using enclosed needles.
- Care should be taken with a disinfectant fourth-grade ammonium salt such as benzalkonium chloride and chlorhexidine.
- Hyruan ONE is intended for single use only. No re-sterilization or reuse is allowed.
- The injection site must be sterilized either by alcohol or other disinfecting solution prior to administration.
- Used syringe, needle and unused materials all need to be discarded after administration.
- In case Hyruan ONE is administered into both sides of knee, a separate product should be applied for each administration.

3) Precautions in handling

- Hyruan ONE should be kept away from children.
- Hyruan ONE should be kept in original packaging; storing Hyruan ONE in a different packaging can cause misuse of product or decline the product quality.

Incompatibilities

Incompatibilities of Hyruan ONE with other drug substances have not been established.

Use in specific populations

1) Pregnancy
Safety and efficacy of Hyruan ONE in pregnant women have not been established so that Hyruan ONE should be cautiously administered to pregnant women or woman having possibility to be pregnant only when the benefits of treatments outweigh danger.
2) Nursing mothers
Safety and efficacy of Hyruan ONE in nursing mother have not been established, thus Hyruan ONE should be cautiously administered to nursing mother only when the benefits of treatments outweigh danger.
3) Pediatrics
The pediatric Safety and efficacy of Hyruan ONE have not been established thus Hyruan ONE should be administered to children with caution when administered unavoidably.
4) Geriatrics
Since physiological function of geriatrics has a declining tendency, administration should be made carefully. Difference in the incidence of reaction rate after administration of Hyruan ONE between geriatrics and non-geriatric group was not observed in the clinical trial.

Storage condition

Store at 1-30℃ free from light, in a hermetic container.

How supplied

3.0mL/syringe × 1syringe

Shelf life

24 months

※ Date of Issue : February 11, 2015

Symbols on packaging

 Do not re-sterilize

 Do not re-use

 Consult instructions for use

 Do not use if package is damaged

 Caution

 Manufacturer

 Batch code

 Temperature limitation

 Keep away from sunlight

 Sterilized stream or dry heat

 Keep away from rain

 Date of manufacture

 Use-by date

 **LG Life Sciences**
129, Seokam-ro, Iksan-si,
Jeollabuk-do, Korea

125*190mm(front)

Absorbable viscosupplement for osteoarthritis

IP M206472-2

Description

Hyruan ONE is a milestone viscosupplementation treatment used to relieve pain associated with osteoarthritis of the knee, shoulder, hip and ankle. Hyruan ONE is a clear, colorless and viscous gel containing 60mg BDDE(1,4-butanediol diglycidyl ether)-crosslinked sodium hyaluronate in a 3mL pre-filled syringe with rubber stoppers on the both ends.

Contents

Each Hyruan ONE (volume of 3.0mL/syringe) contains:
Active ingredient: BDDE-crosslinked sodium hyaluronate gel - 3.0g
(BDDE-crosslinked sodium hyaluronate - 60mg)

Intended Use

For use as a symptomatic treatment for osteoarthritis of weight bearing joints including knee, shoulder, hip and ankle.

Dosage

Hyruan ONE is a single dose for single injection, and should only be injected once per treatment course. The recommended dose is 3mL for knee, shoulder or hip joint and 1-2mL for ankle joint.

Mode of Action

Naturally-occurring hyaluronic acid in the various part of the body acts in the joints as a shock absorber, protects cartilage, and relieves pain associated with osteoarthritis.

Contraindications

- 1) Hyruan ONE should not be injected:
 - To patients who are known to be sensitive to the product or its ingredient(s)
 - To patients with an infection or severe inflammation at joint cavity
 - To patients with a disease or infection of skin near the injection site
- 2) Hyruan ONE should be injected with cautions:
 - In patients with hypersensitivity to other substance
 - In patients with liver disease or prior history of the liver disease

Adverse Events

- 1) Serious adverse events
Shock: since symptoms of shock (frequency unknown) may occur, careful observation should be made. In case abnormalities are noticed, administration should be discontinued and proper measures should be taken.
- 2) Reported adverse events in clinical trials
In a pivotal clinical trial of Hyruan ONE in patients with the knee osteoarthritis (a total of 285 subjects), the occurrence rate of local reaction at injection site after injection into articular cavity was 48.9% (68/139 subjects) in the study group and 49.3% (72/146 subjects) in the control group. The reported adverse events are shown in Table 1. The serious adverse events were reported in the decreasing order of pain (7.2% of the study group, 6.2% of the control group), redness (5.0% of the study group, 2.1% of the control group) and so on. Adverse events lasting over a week were pain 5.3%, redness 1.4%, swelling and warmth each accounting for 1.0%, yet all of them completely resolved in two weeks without any special treatment.

Table 1. Reported local reaction at injection area in the clinical trial

Adverse events	Hyruan ONE N=139(%)	Hyruan Plus*(Control) N=146(%)
Pain	60(43.2)	55(37.7)
Redness	27(19.4)	31(21.2)
Swelling	17(12.2)	21(14.4)
Warmth	33(23.7)	30(20.5)

* Hyruan Plus, requiring three weekly injections, is a viscosupplement for osteoarthritis developed by LG Life Sciences, Ltd.

- 3) Of the patients with the knee osteoarthritis (285 subjects), the occurrence rate of the adverse events excluding ones occurred at the injection site is 34.5% (48/139, 73 cases) for the study group and 28.8% (42/146, 63 cases) for the control group. Most of them were mild to moderate. Table 2 indicates the adverse events occurred in more than 1% of the study group.

Table 2. Adverse events reported in more than 1% of the subjects treated with Hyruan ONE in the clinical trial

Adverse events	Hyruan ONE N=139(%)	Hyruan Plus*(Control) N=146(%)
Nasopharyngitis	7 (5)	3 (2.1)
Pain	4 (2.9)	1 (0.7)
Warmth	3 (2.2)	0 (0)
Upper respiratory tract infection	3 (2.2)	0 (0)
Pain in extremity	3 (2.2)	6 (4.1)
Cystitis	2 (1.4)	2 (1.4)
Dyspepsia	2 (1.4)	1 (0.7)
Parosmia	2 (1.4)	0 (0)
Redness	2 (1.4)	0 (0)
Joint swelling	2 (1.4)	0 (0)
Musculoskeletal pain	2 (1.4)	0 (0)
Plantar fasciitis	2 (1.4)	0 (0)

125*190mm(back)

Precautions

- 1) General precaution:
 - Administration of Hyruan ONE to severely inflamed joints caused by deformative osteoarthritis can lead to exacerbation of the local inflammation symptom. Thus, Hyruan ONE is desired to be given after the removal of existing inflammation symptoms.
 - Administration of Hyruan ONE may occasionally cause local pain or swelling. Thus, patients should be informed to avoid strenuous exercise or action leading joint pain of knee up to 48 hours after administration, and actions such as local relaxation should be guided after injection.
 - Leakage of Hyruan ONE other than articular cavity might cause pain. Thus, Hyruan ONE should be accurately injected into articular cavity.
- 2) Precautions in use:
 - The product should be administered by trained doctors.
 - Since Hyruan ONE is injected directly into the joints, administration should be performed under intact sterilization status.
 - In case of retention of articular fluid, the fluid should be removed by puncturing prior to administration of Hyruan ONE.
 - Intravascular injection, extra-articular injection or injection in the synovial tissues should be avoided.
 - It is desirable to administer the product using enclosed needles.
 - Care should be taken with a disinfectant fourth-grade ammonium salt such as benzalkonium chloride and chlorohexidine.
 - Hyruan ONE is intended for single use only. No re-sterilization or reuse is allowed.
 - The injection site must be sterilized either by alcohol or other disinfecting solution prior to administration.
 - Used syringe, needle and unused materials all need to be discarded after administration.
 - In case Hyruan ONE is administered into both sides of knee, a separate product should be applied for each administration.
- 3) Precautions in handling
 - Hyruan ONE should be kept away from children.
 - Hyruan ONE should be kept in original packaging; storing Hyruan ONE in a different packaging can cause misuse of product or decline the product quality.

Incompatibilities

Incompatibilities of Hyruan ONE with other drug substances have not been established.

Use in specific populations

- 1) Pregnancy

Safety and efficacy of Hyruan ONE in pregnant women have not been established so that Hyruan ONE should be cautiously administered to pregnant women or woman having possibility to be pregnant only when the benefits of treatments outweigh danger.
- 2) Nursing mothers

Safety and efficacy of Hyruan ONE in nursing mother have not been established, thus Hyruan ONE should be cautiously administered to nursing mother only when the benefits of treatments outweigh danger.
- 3) Pediatrics

The pediatric Safety and efficacy of Hyruan ONE have not been established thus Hyruan ONE should be administered to children with caution when administered unavoidably.
- 4) Geriatrics

Since physiological function of geriatrics has a declining tendency, administration should be made carefully. Difference in the incidence of reaction rate after administration of Hyruan ONE between geriatrics and non-geriatric group was not observed in the clinical trial.

Storage condition

Store at 1~30°C free from light, in a hermetic container.

How supplied

3.0mL/syringe × 1syringe

Shelf life

24 months

※ Date of Issue : February 11, 2015

Symbols on packaging

	Do not re-sterilize		Temperature limitation
	Do not re-use		Keep away from sunlight
	Consult instructions for use		Sterilized stream or dry heat
	Do not use if package is damaged		Keep away from rain
	Caution		Date of manufacture
	Manufacturer		Use-by date
	Batch code		

 **LG Life Sciences**
129, Seokam-ro, Iksan-si,
Jeollabuk-do, Korea