ENTERAL FEEDING TUBES

GB Nutricare™ ISOSAF™ Long Term Polyurethane Nasogastric Feeding Tube with ENFit® ISO 80369-3 Connector

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GBUK Group Ltd.
Woodland House
Blackwood Hall Business Park
North Duffield, Selby
North Yorkshire
YO8 5DD, UK

EU: +44 (0)1757 282 945
US: 18665 289 422
www.gbukenteral.com
info@gbukenteral.com

EC REP QualRep Services BV
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

GBUK Enteral is a trading name of
GBUK Enteral Ltd.

www.gbukgroup.com

NC-ROW-03 08

GB

ENTERAL FEEDING TUBES

Instructions For Use

Please read these instructions before use.

FOR ENTERAL USE ONLY

Nutricare™-ISOSAF™ is a single-use medical device, sterilised by ethylene oxide. Nutricare™-ISOSAF™ fitted with ENFit® safety connector in accordance with ISO 80369-3. Fully radio-opaque tubing. Graduations printed every centimetre on the tube to aid insertion and checking for signs of migration. The device is intended for adult and paediatric patients who require enteral administration of medication, fluid and/or nutrition, or drainage of gastric secretions.

INTENDED USE

The Nutricare devices are intended for insertion into the alimentary canal through either the nasal or oral cavity for the purpose of supplying fluids, medication or nutrition. The device is also intended to be used for draining or aspiration of gastric secretions.

INTENDED USER

Clinical personnel and others with appropriate training and knowledge.

INDICATIONS FOR USE

The Nutricare devices are indicated for use with liquids and dispersible / soluble formulations.

CONTRAINDICATIONS

Severe facial trauma (cribriform plate disruption), in this instance the tube may be inserted orally; Patients with altered anatomy, e.g oesophageal fistula or pharyngeal pouch; Unstable cervical spinal injuries (involving vertebrae 4 or above); Basal Skull fractures; Major maxillofacial surgery; Oesophageal varices, ulceration, haemangioma, tumours or surgery; Patients who are known to have coagulopathy or receiving anticoagulant medication.

Infants on Nasal CPAP; Choanal atresia; Respiratory distress.

The device should not be used with solutions of mixed medication and feed or crushed enteric coated or sustained controlled release medications.

WARNINGS & PRECAUTIONS

- Components are supplied sterile.
- Do not re-use or re-sterilise.
- Do not use if packaging is damaged.
- Re-use can result in a loss of device performance and the risk of infection.
- Routinely check flow and tube impairment. If the tube is used for drainage, failure to clear the obstruction may cause gas & fluid buildup in stomach, aspiration of gastric contents, aspiration pneumonia and other complications.
- The device should not be used for high volume decompression in emergency situations, and/or decompression of highly viscous fluids other than gastric secretions.
- If medication is not available in a dispersible/soluble formulation where crushing tablets or opening capsules is required, the prescriber and practitioner accept liability for any adverse effects resulting from this administration.

PACK CONTENTS

Nutricare™-ISOSAF™ Long Term NG Tubes

1x Nutricare™-ISOSAF™ Nasogastric Enteral Feeding Tube

1x Fixation plaster

DIRECTIONS

- Explain procedure to patient.
- Place patient in a supported, upright position. **Warning:** patient should not lean forward and the head and neck should not be extended.
- Remove enteral feeding tube from package. Make sure the guidewire connector stays firmly connected to the administration port of the enteral feeding tube during insertion (N/a to Nutricare™-ISOSAF™ 14Fr and 16Fr tubes without guidewire).
- The tube length should be estimated before insertion.
- Check nostrils to identify any obstructions or abnormalities liable to prevent insertion, and determine the preferred nostril for insertion.
- Insert the rounded end of the tube into the nostril, aiming the tip parallel to the nasal septum, and advance the tube towards the nasopharynx. If any obstruction is felt, withdraw the tube and try again in a slightly different direction - or try the other nostril.
- If the patient has been assessed to have a safe swallow offer sips of water, advance the tube until the predetermined length of the tube has been reached.
- Warning:** If the patient shows signs of distress, coughing, gasping, or cyanosis this could indicate that the tube has passed into the trachea and the tube should be removed immediately! N.B. The absence of these symptoms does not necessarily indicate safe placement of the tube.
- When the predetermined length of tube has been inserted, correct tube position must be confirmed. This can be done with or without the guidewire. Ensure the guidewire is removed after confirmation. The tube is fully radiopaque without the guidewire if X-ray confirmation is chosen (N/a to Nutricare™-ISOSAF™ 14Fr and 16Fr tubes without guidewire).
- Secure the tube in line with local policy.
- Enteral feeding tubes are not flushed, nor any liquid/feed introduced through the tube following initial placement, until the tube tip is confirmed by pH testing or chest X-ray to be in the stomach. (NPSA/2011/PSA002)
- Warning:** For duodenal or jejunal placement X-Ray confirmation is usually indicated.
- If tube position is incorrect the tube must be removed!
- Never re-insert the guidewire when the tube is in a patient!
- Document the placement procedure in line with local policy.

RISK ASSESSMENT

Individual risk assessments should be conducted and documented by a competent person. For patients who may have a high risk of incorrect tube positioning, dislodgement or aspiration appropriate specialist advice should be taken.

CLEANING & DISINFECTION

N/A - Single use device. Dispose after use.

DISPOSAL

Dispose of waste in line with local policy.

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