ENTERAL FEEDING TUBES

GB Nutricare™ Infant Clear-ISOSAF™ Polyurethane Nasogastric Feeding Tube with ENFit® ISO 80369-3 Connector

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€1639

GBUK Group Ltd.
Woodland House
Blackwood Hall Business Park
North Duffield, Selby
North Yorkshire
YO8 5DD, UK

Enteral UK is a trading name of
GBUK Enteral Ltd.

www.gbukgroup.com

NCI-03 09

GB

ENTERAL FEEDING TUBES

Instructions For Use

Please read these instructions before use.

FOR ENTERAL USE ONLY

Nutricare™ Infant Clear-ISOSAF™ is a single-use medical device, sterilised by ethylene oxide. Nutricare™ Infant Clear-ISOSAF™ fitted with ENFit® safety connector in accordance with ISO 80369-3. Radio-opaque stripe printed along the length of the tube. Graduations printed every centimetre on the tube to aid insertion and checking for signs of migration. The device is intended for infant patients who require enteral administration of medication, fluid and/or nutrition.

INTENDED USE

The Nutricare feeding tube is intended for insertion into the alimentary canal through either the nasal or oral cavity for the purpose of supplying fluids, medication or nutrition.

INTENDED USER

Clinical personnel and others with appropriate training and knowledge.

INDICATIONS FOR USE

The Nutricare devices are indicated for use with liquids and dispersible / soluble formulations.

CONTRAINDICATIONS

Severe facial trauma (cribriform plate disruption), in this instance the tube may be inserted orally; Patients with altered anatomy, e.g oesophageal fistula or pharyngeal pouch; Unstable cervical spinal injuries (involving vertebrae 4 or above); Basal Skull fractures; Major maxillofacial surgery; Oesophageal varices, ulceration, haemangioma, tumours or surgery; Patients who are known to have coagulopathy or receiving anticoagulant medication; Infants on Nasal CPAP; Choanal atresia; Respiratory distress.

The device should not be used with solutions of mixed medication and feed or crushed enteric coated or sustained controlled release medications.

WARNINGS & PRECAUTIONS

- Components are supplied sterile.
- Do not re-use or re-sterilise.
- Do not use if packaging is damaged.
- Re-use can result in a loss of device performance and the risk of infection.
- Routinely check flow and tube impairment. If the tube is used for drainage, failure to clear the obstruction may cause gas & fluid buildup in stomach, aspiration of gastric contents, aspiration pneumonia and other complications.
- If medication is not available in a dispersible/soluble formulation and crushing tablets or opening capsules is required, the prescriber and practitioner accept liability for any adverse effects resulting from this administration.

PACK CONTENTS

1 x Nutricare™ Infant Clear-ISOSAF™ Nasogastric Enteral Feeding Tube
1 x Fixation plaster

DIRECTIONS

- Explain procedure to patient.
- Place patient in a supported, upright position. **Warning:** patient should not lean forward and the head and neck should not be extended.
- Remove enteral feeding tube from package.
- The tube length should be estimated before insertion.
- Check nostrils to identify any obstructions or abnormalities liable to prevent insertion, and determine the preferred nostril for insertion.
- Insert the rounded end of the tube into the nostril, aiming the tip parallel to the nasal septum, and advance the tube towards the nasopharynx. If any obstruction is felt, withdraw the tube and try again in a slightly different direction - or try the other nostril.
- If the patient has been assessed to have a safe swallow offer sips of water, advance the tube until the predetermined length of the tube has been reached.
- Warning:** If the patient shows signs of distress, coughing, gasping, or cyanosis this could indicate that the tube has passed into the trachea and the tube should be removed immediately! N.B. The absence of these symptoms does not necessarily indicate safe placement of the tube.
- When the predetermined length of tube has been inserted, correct tube position must be confirmed.
- Secure the tube in line with local policy.
- Enteral feeding tubes are not flushed, nor any liquid/feed introduced through the tube following initial placement, until the tube tip is confirmed by pH testing or chest X-ray to be in the stomach. (NPSA/2011/PSA002)
- If tube position is incorrect the tube must be removed!
- Document the placement procedure in line with local policy.

RISK ASSESSMENT

Individual risk assessments should be conducted and documented by a competent person. For patients who may have a high risk of incorrect tube positioning, dislodgement or aspiration appropriate specialist advice should be taken.

CLEANING & DISINFECTION

N/A - Single use device. Dispose after use.

DISPOSAL

Dispose of waste in line with local policy.

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