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Cerebral Ventricular Catheter, Atrial Catheter, Peritoneal Catheter, Connector and Reservoir


User's Manual Description

Cerebral Ventricular Catheter:

Cerebral Catheters are made of medical grade silicone with radiopaque fillet, in Adult, Infantile and Neonatal sizes. Catheters are presented in straight, right-angled and right-angled with reservoir models. The reservoir is pre-assembled and is indicated for puncture.

The straight models have a 90° angle former with fixation lugs.


Refer to table for dimensions and models.

Peritoneal Catheter:

Peritoneal Catheters are made of medical grade silicone, in Adult and Infantile sizes and different lengths. All catheters show marks to direct introduction at every 10 cm. The catheters have opened distal tip.


Refer to table for dimensions and models.

Atrial Catheter:

Atrial Catheters are made of medical grade silicone, in five models combining two tubes of different diameter and length. Models C, and and H catheter tubes are provided preconnected and model B also presents a straight poliactal connector for tube connection.

Refer to table for dimensions and models.

Reservoirs:

The reservoirs are made of medical grade silicone (dome), and poliactal (base), in two models; one for drainage and puncture and the other for puncture only:

Model RHM for drainage must be interspersed between a straight cerebral catheter and the valve, and is available in Adult and Infantile sizes.

Model RCO-V for puncture with vertical output is available in three sizes, and must be connected to a straight-model cerebral catheter.


Model RCO-L for puncture with lateral output is available in three sizes, and must be connected to a straight-or angular-model cerebral catheter. The connection is made through the straight connector that comes with the reservoir.

The reservoirs have fixation holes for suture to the subcutaneous tissue.

The reservoirs do not have metal components that may cause interference in CT scan or magnetic resonance imaging.

The reservoirs have a puncture chamber with a hard poliactal base, which protects against excessive penetration of the needle during puncture. The thickness of the dome is designed to support repeated punctures. Its high profile allows for easy localization under the cranium skin.

Refer to table for dimensions and models.


Connectors:

The connectors are made of poliactal in the Straight, Y (short and long sizes) and Angular models. The external diameter of the connectors is compatible with all cerebral and peritoneal catheters.

The model B connector makes the connection with atrial catheters.

The Y connector enables the connection between two valvular systems to a peritoneal catheter in the derivation of two ventricles.

Refer to table for dimensions and models.

If the integrity of the package is maintained the components remain sterile, during validity term of sterilization. They are packaged in a two-layer envelop of surgical grade paper and film PET/PP and sterilized with ethylene oxide. The external envelop, after opening, allows access to and removal of the internal envelop to the sterile surgical site. The latter, after opening, gives access to the components.

The silicone components are impregnated with bari-um, allowing visualization through x-ray.

The components have universal connections.

Indication

The products are designed as components for the intracranial pressure control system in post-surgical drainage and Hydrocephalus, with ventricle-peritoneal or ventricle-atrial derivation.

The components individually supplied are designed both for the new Shunt System assembly and for components replacement during checking of the implanted system.

Cerebral Reservoirs are designed for puncture in cerebral liquor aspiration and/or for administering medication to the cerebral ventricle.

Contraindication

The implant of the Shunt System and its components is contraindicated in patients with knowledge or suspicion of infection on the implant region of anyone of the components (meningitis, ventriculitis, skin infections, bacteremia, septicemia, peritonitis, etc).

The implant of the ventricle-atrial Shunt System is contraindicated in cases of pulmonary embolism or pressure increase on the right atrium, and congenital cardiopathies.

Warnings

The implant of the Derivation System in cases with hemorrhagic liquor or high protein content must be carefully assessed due to increased risk of obstruction.

The post-surgical of patients implanted with derivation system and its components must be closely observed for detection and evaluation of obstruction symptoms or hyper-drainage.

The overdrainage may be discerned by the signs and symptoms of reduced intracranial pressure and to the development of a subdural hematoma, subdural hygroma, or collapse of the lateral ventricular wall.

The obstruction may lead to the signs and symptoms of increased intracranial pressure. In children the signs may be increased tension of the fontanelle, congestion of the scalp veins, apathy, drowsiness, rigidity of the nape of the neck, irritability and vomiting. In adult patients, the signs may be: vomiting, headaches, loss of conscience, rigidity of the nape of the neck, blurred vision.

If some of these occur, immediate evaluation for surgical intervention is necessary. The choice of the procedure to be used is left at the surgeon's discretion.

The product implant is restricted to surgeon doctors with specific training and qualification in this procedure.

Do not use this device for other ends rather than the foreseen and estipulated by the manufacturer.

This product must be handled by qualified personnel only.

Reesterilization

This is a single-use product.

Do not reesterilize this device.

Disposal

The product should be disposed as hospital waste, or according regional regulation.

If the sterilization validity of the product has expired, or in case of damaged product packaging, please return to the manufacturer for processing or discard as convenient.

Usage Instructions

Implant Procedures

General

The surgical procedure for components implant, including positioning of the ventricular catheter and the reservoir over the cranium can be directed by different implant techniques, and the choice of the most suitable one may be left at the discretion of the surgeon.

The components assembly in the Shunt System is made through the connection between proximal (ventricular) catheter, valve and distal (peritoneal, atrial) catheter. A reservoir for cerebral drainage (models RHM-1, RHM-3) can be placed between the ventricular catheter and the valve.

Reservoir

The Cerebral Reservoir assembly is made through the connection between ventricular catheter and reservoir.

Atrial Catheter

The Atrial Catheter has a guide thread to assist in the introduction of the distal tip of the catheter in the patient's vein.

To remove the guide thread after introduction: Connect the catheter to a syringe with alcohol passing the tip of the thread through the syringe, and lightly inject an amount of alcohol of no more than 0.5ml. Disconnect the syringe and pull the thread.

Precautions

Patients or the persons in charge of them must be informed of possible complications and warned about the precautions involving the implant of this product.

Handle the material under strict aseptic technique.

The silicone elastomer is very electrostatic and susceptible of contamination from air or surface particles. Avoid contact with surfaces releasing threads or particles. Handle the material with sterile and clean gloves, avoiding contact with frayed fabrics, talc or other contaminants.

The safety of the connections can be reinforced with ligatures. Care must be taken to make the connections and ligatures in such a way that the catheters do not suffer cuts or occlusion. Care must be taken also, not to squeeze catheters at any point, because there is risk of obstruction of the system.

Avoid contact with piercing or cutting materials when opening the package and handling the components of the Shunt System.

When using hypodermic needles for injection or puncturing of the reservoir, calibers greater than 0.5mm should not be used. Assure the absence of air and the correct introduction depth of the ventricular catheter through observation of liquor output on the proximal tip.

Complications

The use of this device may cause complications similar to those associated with other surgical procedures, or to any foreign object implanted in the body, depending on the degree of reaction and sensitivity of the patient.

The most common complications when using shunt systems and their components are infection, obstruction, system disconnection, mechanical failure of components, overdrainage and insufficient drainage.

The use of the ventricular catheter may cause infection, subcutaneous infiltration of liquor, neurological sequela or intraventricular obstruction by particles (including hemorrhagic liquor or with high protein content).

The ventricular catheter can be obstructed by the adhesion to choroid plexus, if inserted in the brain or by the coaptation of the ventricular walls.

The ventricle-atrial derivation may cause infection and predisposition to infection dissemination to other areas or organs, skin erosion, pulmonary embolism, thrombosis, cardiac perforation and cardiac tamponade, pulmonary hypertension.

The ventricle-peritoneal derivation may cause inflammation of the abdominal cavity, skin erosion, intestinal perforation, cysts or pseudocysts, umbilical fistula and ascitis. The Peritoneal Catheter is subject to obstruction by blood clots, debris, bacterial colonization, aggregated cells of tumor, and may suffer fracture.

The squeezing of the catheters may cause system obstruction or drainage decrease.

Growth of the patient may cause displacement of the Atrial Catheter from the venous system or of the Peritoneal Catheter from the peritoneal cavity, deriving liquor to areas of difficult absorption.

The disconnection of the catheters or fracture of a component can cause migration of the components into the right atrium, peritoneal cavity, lateral ventricle of the brain or adjacent areas.

The use of a hydrocephalus shunt system has been related to the incidence of epilepsy, and incidence of apoplexy with the constant revision of the catheter.

In children, the extreme reduction of the intracranial pressure will cause marked depression of the anterior fontanelle, overriding of the cranial bones and can transform a communicating into a non-communicating hydrocephalus.

Cited References

This list of sources used as reference in this manual can be obtained by written request to be sent to Hpbio Address or through the website www.hpbio.com.br.

Communication Form

Please fill the communication form provided with the packing, if any comment on the product is necessary before or after the implant, and send it directly to Hpbio or through the supplier or representative. If possible, attach the part or component that generated the comment with serial number and lot.

Product		Features		
Ventricular Cerebral Catheter	Length (mm)	Ø Internal (mm)	Ø External (mm)	Catalog no.
Neonatal Straight	70	1,0	2,0	VN7
	100	1,0	2,0	VN10
	150	1,0	2,0	VN15
	220	1,0	2,0	VN22
Infantile Straight	150	1,2	2,5	VI15
	200	1,2	2,5	VI20
	235	1,2	2,5	VI23.5
	250	1,2	2,5	VI25
	300	1,2	2,5	VI30
	350	1,2	2,5	VI35
Adult Straight	150	1,4	3,0	VA15
	200	1,4	3,0	VA20
	235	1,4	3,0	VA23.5
	250	1,4	3,0	VA25
	300	1,4	3,0	VA30
	350	1,6	3,0	VA35
Neonatal Right-Angled	50	1,0	2,0	VNG6
	60	1,0	2,0	VNG5
	70	1,0	2,0	VNG7
Neonatal Right-Angled with Reservoir	50	1,0	2,0	VNG5R
	60	1,0	2,0	VNG6R
	70	1,0	2,0	VNG7R
Infantile Right-Angled	50	1,2	2,5	VIG5
	60	1,2	2,5	VIG6
	70	1,2	2,5	VIG7
Infantile Right-Angled with Reservoir	50	1,2	2,5	VIG5R
	60	1,2	2,5	VIG6R
	70	1,2	2,5	VIG7R
Adult Right-Angled	70 x 30	1,4	3,0	VAG7
	70 x 20	1,4	3,0	VAG7.2
	90 x 30	1,4	3,0	VAG9
	90 x 20	1,4	3,0	VAG9.2
	100 x 30	1,4	3,0	VAG10
Adult Right-Angled with Reservoir	70 x 30	1,4	3,0	VAG7R
	100 x 30	1,4	3,0	VAG10R
Peritoneal Catheter	Length (mm)	Ø Internal (mm)	Ø External (mm)	Catalog no.
Infantile Peritoneal Catheter	520	1,0	2,0	PI 52
	900	1,0	2,0	PI 90
Adult Peritoneal Catheter	900	1,2	2,5	PA 90
	1.020	1,2	2,5	PA 102
Cerebral Reservoir	Length (mm)	Height	Ø Burr hole	Catalog no.
Infantile Ventricular Drainage Reservoir	13	4	6	RHM-1
Adult Ventricular Drainage Reservoir	17	5,5	9,5	RHM-3
Small Cerebral Reservoir – Horizontal outlet	13	6	—	RCO-1L
Medium Cerebral Reservoir – Horizontal outlet	17	7	—	RCO-2L
Large Cerebral Reservoir – Horizontal outlet	21	7	—	RCO-3L
Small Cerebral Reservoir – Vertical outlet	13	6	6	RCO-1V
Medium Cerebral Reservoir – Vertical outlet	17	7	9,5	RCO-2V
Large Cerebral Reservoir – Vertical outlet	21	7	13,5	RCO-3V
Atrial Catheter (venous)	Length (mm)	Length (mm)		Catalog no.
	Tube of 2,5 x 1,2	Tube of 1,2 x 0,64		
Atrial Catheter type A	450	—		AA
Atrial Catheter type B	220	150		AB
Atrial Catheter type C	380	30		AC
Atrial Catheter type E	220	220		AE
Atrial Catheter type H	30	380		AH
Conector	Length (mm)	Ø Internal (mm)	Ø External (mm)	Catalog no.
Straight A short	11	1,0	1,8	CAC
Straight A long	17	1,0	1,8	CAL
Straight B	08	0,8	1,4 / 1,8	CAB
90°	08	1,0	1,8	C 90
Y short	09	1,0	1,8	CYC
Y long	12	1,0	1,8	CYL

Catéter Ventricular Cerebral, Catéter Peritoneal, Catéter Atrial, Conector y Reservorio

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