

Data Sheet

☐ Single use sterile device for use during surgical procedures.		
Unique simple solution for cleaning the tips of sharps instrument		
and electrosurgical pencil tips.		
Latex free.		
Easy and friendly to use.		

Product Description:

Product	Tip Cleaner
The function explaining	The tip cleaner is used during surgical procedures to clean the tips of surgical instruments from the residuals as, tissue, soft tissue, residuals, etc. that stocked on the tips.
Products material	The product is made of: Square foam pad coated with abrasive material on top of it.
Sterilization method	ETO sterilization
Shelf life Validity	5 years
Product specification	The tip cleaner is a complementary device used during the different surgical procedures, in order to keep clean the tips of the surgical instruments, diathermy pencils and other surgical devices.