




Target™ 360

DETACHABLE COILS

Conformability from Frame to Finish

Open loop configuration designed for uniform distribution, concentric filling and neck coverage.




Target XL™

DETACHABLE COILS

Larger and Longer

2x more volume per coil than Target-10 size coils and up to 50cm in length; all consistently deliverable through an Excelsior™ SL-10™ Microcatheter.



Target™ Long

DETACHABLE COILS

Length for Efficiency

Longer 10-size Target coils for more efficient treatment.



Target™ Nano™

DETACHABLE COILS

Smallest and Softest

Precisely engineered and extremely soft* 1mm and 1.5mm 360 and helical coils provide a new treatment perspective.



InZone™

DETACHMENT SYSTEM

Smarter Detachment

Electrolytic detachment affords a small and flexible coil junction, enhancing coil delivery.

Target™ Detachable Coil

See package insert for complete indications, contraindications, warnings and instructions for use.

INTENDED USE / INDICATIONS FOR USE

Target Detachable Coils are intended to endovascularly obstruct or occlude blood flow in vascular abnormalities of the neurovascular and peripheral vessels.

Target Detachable Coils are indicated for endovascular embolization of:

- Intracranial aneurysms
- Other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae
- Arterial and venous embolizations in the peripheral vasculature

InZone™ Detachment System


See package insert for complete indications, contraindications, warnings and instructions for use.

INDICATIONS FOR USE

The InZone Detachment System is intended for use with all versions of Stryker Neurovascular detachable coils in the embolization of intracranial aneurysms and other vascular malformations of the neuro and peripheral vasculature.

THIS DOCUMENT IS INTENDED SOLELY FOR THE USE OF HEALTHCARE PROFESSIONALS.


A physician must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that physicians be trained in the use of any particular product before using it in a procedure. The information presented is intended to demonstrate the breadth of Stryker product offerings. A physician must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. The Stryker products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Target™

DETACHABLE COILS

Smooth and Stable



*Bench test results. n=5. Bench test results may not necessarily be indicative of clinical performance. Data on file.





Target™ DETACHABLE COILS

Smooth and Stable

Consistently smooth deployment and exceptional microcatheter stability.

Selection Guide



360 Shape



Helical Shape

Target XL™ 360 Standard Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
6	20	M0036106200
7	20	M0036107200
8	30	M0036108300
9	30	M0036109300
10	40	M0036101040
12	45	M0036101240
14	50	M0036101450
16	50	M0036101650
18	50	M0036101850
20	50	M0036102050
22	50	M0036102250
24	50	M0036102450

Target XL 360 Soft Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
6	20	M0036126200
7	20	M0036127200
8	30	M0036128300
9	30	M0036129300
10	40	M0036121040
12	45	M0036121240
14	50	M0036121450
16	50	M0036121650
18	50	M0036121850
20	50	M0036122050

Target™ 360 Standard Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
3	6	M0035463060
	8	M0035463080
4	8	M0035464080
	10	M0035464100
	15*	M0035464150
5	15	M0035465150
	20*	M0035465200
6	15	M0035466150
	20*	M0035466200
	30*	M0035466300
7	15	M0035467150
	20	M0035467200
8	30*	M0035467300
	20	M0035468200
	30*	M0035468300
9	20	M0035469200
	30	M0035469300
10	30	M0035461030
11	30	M0035461130
12	30	M0035461230
13	30	M0035461330
14	30	M0035461430
15	40	M0035461540

Target 360 Soft Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
2	4	M0035472040
3	6	M0035473060
	8	M0035473080
	10*	M0035473100
	6	M0035474060
	8	M0035474080
4	10	M0035474100
	15*	M0035474150
5	10	M0035475100
	15	M0035475150
	20*	M0035475200
	10	M0035476100
6	15	M0035476150
	20*	M0035476200
	30*	M0035476300
7	15	M0035477150
	20	M0035477200
	30*	M0035477300
8	20	M0035478200
	30*	M0035478300
9	20	M0035479200
	30	M0035479300
10	30	M0035471030
11	30	M0035471130
12	30	M0035471230
13	30	M0035471330
14	30	M0035471430

Target™ 360 Nano™ Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
1	2	M0035421020
	3	M0035421030
1.5	2	M0035421520
	3	M0035421530
	4	M0035421540

Target 360 Ultra Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
2	3	M0035422030
	4	M0035422040
2.5	6*	M0035422060
	4	M0035422540
3	4	M0035423040
	6	M0035423060
	8	M0035423080
	10*	M0035423100
4	6	M0035424060
	8	M0035424080
	10	M0035424100
5	15*	M0035424150
	10	M0035425100
	15*	M0035425150

Target XL™ Helical Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
5	15†	M0036115150
6	20†	M0036116200
7	20†	M0036117200
8	30†	M0036118300
9	30†	M0036119300
10	40†	M0036111040
12	50	M0036111250
14	50	M0036111450
16	50	M0036111650
18	50	M0036111850
20	50	M0036112050
22	50	M0036112250
24	50	M0036112450

Target Coil Accessories	
Description	Product Number
InZone™ 2.0 Detachment System	M00345100950
IZDS™ Connecting Cable	M00345110250

Target Helical Nano Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
1	1	M0035431010
	2	M0035431020
1.5	3	M0035431030
	1	M0035431510
	2	M0035431520
	3	M0035431530
	4	M0035431540

Target Helical Ultra Coils		
Secondary Diameter (mm)	Length in Introducer (cm)	Product Number
2	1	M0035432010
	2	M0035432020
	3	M0035432030
	4	M0035432040
	6	M0035432060
2.5	8*	M0035432080
	3	M0035432530
	4	M0035432540
	6	M0035432560
3	4	M0035433040
	6	M0035433060
	8	M0035433080
	10*	M0035433100
4	6	M0035434060
	8	M0035434080
	15*	M0035434150

* Indicates Target Long Detachable Coils
† Comparable to Soft Coils
Not all products are available for sale in all countries due to regulatory requirements.