

PÁGINA: 1 de 4

## PINZA DESECHABLE PARA TOMA DE BIOPSIA



#### Descripción

- Pinza de revestimiento liso lo cual facilita el paso a través del endoscopio,
- Mango ergonómico codificado con color lo que permite identificar rápidamente la longitud del dispositivo.
- Bordes quirúrgicos y cazoletas fenestradas de 2.4mm de diámetro para la recogida de tejido
- Longitudes diferentes en caso de uso en endoscopia o colonoscopia

#### **Indicaciones**

- Obtención endoscópica de biopsias de tejido mucoso
- Extracción de cuerpos extraños

#### **Dimensiones**

Longitud: 160 ó 230cm

Diámetro: de la copa: 2.4mm

#### Características

- Sin espícula o con espícula
- Borde liso o serrado

#### Canal de trabajo

Canal de accesorios mínimo: 2.8mm

#### Contraindicaciones

• Coagulopatia y preparación insuficiente del intestino



PÁGINA: 2 de 4

### **Complicaciones posibles**

 Asociadas a la endoscopia gastrointestinal: perforación, sangrado o hemorragia, aspiración, fiebre, infección, hipotensión, paradas respiratorias y arritmias.

#### **Precauciones**

- Endoscopio debe permanecer recto como sea posible al introducir o extraer las pinzas
- Si se acciona el mango mientras las pinzas están enrolladas, están pueden resultar dañadas

#### **Fabricante**

Wilson-Cook Medical Inc

#### **Almacenamiento**

• Lugar seco y protegido de temperaturas extremas

#### Presentación

Sobre con una unidad

#### Esterilización

Desechable, único uso

#### Disposición final

Residuo contaminado peligroso

#### Clasificación de riesgo

Riesgo Moderado (IIA)

#### Registro sanitario

INVIMA 2018 DM-0018931



PÁGINA: 3 de 4

### Referencias comercializadas

| Order<br>Number | Reference<br>Part Number | Cup<br>Diameter<br>mm | Shaft<br>Length<br>cm | Forceps<br>Type | Cup<br>Edge | Box<br>Quantity | Scope<br>Type | Minimum<br>Accessory<br>Channel<br>mm | Color<br>Code |
|-----------------|--------------------------|-----------------------|-----------------------|-----------------|-------------|-----------------|---------------|---------------------------------------|---------------|
| Mini Ford       | eps                      |                       |                       |                 |             |                 |               |                                       |               |
| G53006          | DBF-1.8-160-S            | 1.8                   | 160                   | non-spiked      | surgical    | 10              | gastroscope   | 2.0                                   | green         |
| Standard        | l Forceps                |                       |                       |                 |             |                 |               |                                       |               |
| G31577          | DBF-2.4-160-S            | 2.4                   | 160                   | non-spiked      | surgical    | 10              | gastroscope   | 2.8                                   | green         |
| G56061          | DBF-2.4-160-20-S         | 2.4                   | 160                   | non-spiked      | surgical    | 20              | gastroscope   | 2.8                                   | green         |
| G31580          | DBF-2.4-160SP-S          | 2.4                   | 160                   | spiked          | surgical    | 10              | gastroscope   | 2.8                                   | green         |
| G31576          | DBF-2.4-230-S            | 2.4                   | 230                   | non-spiked      | surgical    | 10              | colonoscope   | 2.8                                   | red           |
| G56065          | DBF-2.4-230-20-S         | 2.4                   | 230                   | non-spiked      | surgical    | 20              | colonoscope   | 2.8                                   | red           |
| G31585          | DBF-2.4-230SP-S          | 2.4                   | 230                   | spiked          | surgical    | 10              | colonoscope   | 2.8                                   | red           |
| G56074          | DBF-2.4-230SP-20-S       | 2.4                   | 230                   | spiked          | surgical    | 20              | colonoscope   | 2.8                                   | red           |
| G56122          | DBF-2.4SN-160-S          | 2.4                   | 160                   | non-spiked      | serrated    | 10              | gastroscope   | 2.8                                   | green         |
| G56123          | DBF-2.4SN-160-20-S       | 2.4                   | 160                   | non-spiked      | serrated    | 20              | gastroscope   | 2.8                                   | green         |
| G56124          | DBF-2.4SN-160-40-S       | 2.4                   | 160                   | non-spiked      | serrated    | 40              | gastroscope   | 2.8                                   | green         |
| G56128          | DBF-2.4SN-160SP-S        | 2.4                   | 160                   | spiked          | serrated    | 10              | gastroscope   | 2.8                                   | green         |
| G56129          | DBF-2.4SN-160SP-20-S     | 2.4                   | 160                   | spiked          | serrated    | 20              | gastroscope   | 2.8                                   | green         |
| G56130          | DBF-2.4SN-160SP-40-S     | 2.4                   | 160                   | spiked          | serrated    | 40              | gastroscope   | 2.8                                   | green         |
| G56125          | DBF-2.4SN-230-S          | 2.4                   | 230                   | non-spiked      | serrated    | 10              | colonoscope   | 2.8                                   | red           |
| G56126          | DBF-2.4SN-230-20-S       | 2.4                   | 230                   | non-spiked      | serrated    | 20              | colonoscope   | 2.8                                   | red           |
| G56127          | DBF-2.4SN-230-40-S       | 2.4                   | 230                   | non-spiked      | serrated    | 40              | colonoscope   | 2.8                                   | red           |
| G56131          | DBF-2.4SN-230SP-S        | 2.4                   | 230                   | spiked          | serrated    | 10              | colonoscope   | 2.8                                   | red           |
| G56132          | DBF-2.4SN-230SP-20-S     | 2.4                   | 230                   | spiked          | serrated    | 20              | colonoscope   | 2.8                                   | red           |
| G56133          | DBF-2.4SN-230SP-40-S     | 2.4                   | 230                   | spiked          | serrated    | 40              | colonoscope   | 2.8                                   | red           |
| Large Fo        | rceps                    |                       |                       |                 |             |                 |               |                                       |               |
| G53630          | DBF-2.4SL-160-S          | 2.4                   | 160                   | non-spiked      | serrated    | 10              | gastroscope   | 2.8                                   | green         |
| G56057          | DBF-2.4SL-160-20-S       | 2.4                   | 160                   | non-spiked      | serrated    | 20              | gastroscope   | 2.8                                   | green         |
| G56056          | DBF-2.4SL-160-40-S       | 2.4                   | 160                   | non-spiked      | serrated    | 40              | gastroscope   | 2.8                                   | green         |
| G53638          | DBF-2.4SL-160SP-S        | 2.4                   | 160                   | spiked          | serrated    | 10              | gastroscope   | 2.8                                   | green         |
|                 |                          |                       |                       |                 |             |                 |               |                                       |               |



PÁGINA: 4 de 4

| G56070   | DBF-2.4SL-160SP-20-S  | 2.4 | 160 | spiked     | serrated | 20 | gastroscope | 2.8 | green |
|----------|-----------------------|-----|-----|------------|----------|----|-------------|-----|-------|
| G56073   | DBF-2.4SL-160SP-40-S  | 2.4 | 160 | spiked     | serrated | 40 | gastroscope | 2.8 | green |
| G53633   | DBF-2.4SL-230-S       | 2.4 | 230 | non-spiked | serrated | 10 | colonoscope | 2.8 | red   |
| G56052   | DBF-2.4SL-230-20-S    | 2.4 | 230 | non-spiked | serrated | 20 | colonoscope | 2.8 | red   |
| G56059   | DBF-2.4SL-230-40-S    | 2.4 | 230 | non-spiked | serrated | 40 | colonoscope | 2.8 | red   |
| G53641   | DBF-2.4SL-230SP-S     | 2.4 | 230 | spiked     | serrated | 10 | colonoscope | 2.8 | red   |
| G56087   | DBF-2.4SL-230SP-20-S  | 2.4 | 230 | spiked     | serrated | 20 | colonoscope | 2.8 | red   |
|          |                       |     |     |            |          |    |             |     |       |
| G56076   | DBF-2.4SL-230SP-40-S  | 2.4 | 230 | spiked     | serrated | 40 | colonoscope | 2.8 | red   |
| Jumbo Fo | orceps                |     |     |            |          |    |             |     |       |
| G53628   | DBF-2.8SJ-160-S       | 2.8 | 160 | non-spiked | serrated | 10 | gastroscope | 3.2 | green |
| G53650   | DBF-2.8SJ-160SP-S     | 2.8 | 160 | spiked     | serrated | 10 | gastroscope | 3.2 | green |
| G53642   | DBF-2.8SJ-230-S       | 2.8 | 230 | non-spiked | serrated | 10 | colonoscope | 3.2 | red   |
| G53042   | DBF-2.8SJ-230SP-S     | 2.8 | 230 | spiked     | serrated | 10 | colonoscope | 3.2 | red   |
| G56721   | DBF-2.8SJN-230-20-S   | 2.8 | 230 | non-spiked | serrated | 20 | colonoscope | 3.2 | red   |
| G56722   | DBF-2.8SJN-230-40-S   | 2.8 | 230 | non-spiked | serrated | 40 | colonoscope | 3.2 | red   |
| G56723   | DBF-2.8SJN-230SP-20-S | 2.8 | 230 | spiked     | serrated | 20 | colonoscope | 3.2 | red   |
| G56724   | DBF-2.8SJN-230SP-40-S | 2.8 | 230 | spiked     | serrated | 40 | colonoscope | 3.2 | red   |
| Max Forc | eps                   |     |     |            |          |    |             |     |       |
| G54166   | DBF-3.3SM-160-S       | 3.3 | 160 | non-spiked | serrated | 10 | gastroscope | 3.8 | green |
| G56054   | DBF-3.3SM-160-20-S    | 3.3 | 160 | non-spiked | serrated | 20 | gastroscope | 3.8 | green |
| G54168   | DBF-3.3SM-160SP-S     | 3.3 | 160 | spiked     | serrated | 10 | gastroscope | 3.8 | green |
| G56053   | DBF-3.3SM-160SP-20-S  | 3.3 | 160 | spiked     | serrated | 20 | gastroscope | 3.8 | green |
|          |                       |     |     |            |          |    |             |     |       |
| Max Forc | eps, continued        |     |     |            |          |    |             |     |       |
| G54167   | DBF-3.3SM-230-S       | 3.3 | 230 | non-spiked | serrated | 10 | colonoscope | 3.8 | red   |
| G56058   | DBF-3.3SM-230-20-S    | 3.3 | 230 | non-spiked | serrated | 20 | colonoscope | 3.8 | red   |
| G54170   | DBF-3.3SM-230SP-S     | 3.3 | 230 | spiked     | serrated | 10 | colonoscope | 3.8 | red   |
| G56055   | DBF-3.3SM-230SP-20-S  | 3.3 | 230 | spiked     | serrated | 20 | colonoscope | 3.8 | red   |
| G56060   | DBF-3.3SM-230SP-40-S  | 3.3 | 230 | spiked     | serrated | 40 | colonoscope | 3.8 | red   |