

[illegible]



Contents:

- 1: Catheter: 20 Ga. x 8 cm Radiopaque PEBA
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Both Tip on Both Ends)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

ARTERIAL CATHETERIZATION SYSTEMS - SISTEMA
DE CATETERES DE ACCESO ARTERIAL / SISTEMA
PARA CATETERIZACION

ARTERIAL
Importado por

TELEFLEX MEDICAL COLOMBIA S.A.S

Carrera 16 #97 - 46, Piso 6,

Bogotá, COLOMBIA

Registro Sanitario N° INVIMA 2016DM-0015463

ARROW

REF SAC-00820-PBX
LOT 71F22D1056
2027-03-31

(IPN820735)
2022-04-22

ARROW

REF SAC-00820-PBX
LOT 71F22D1056
2027-03-31

(IPN820735)
2022-04-22

