



FICHA TÉCNICA DISPOSITIVOS MÉDICOS

PRODUCTO	PREMICRON ® B.BRAUN					
NOMBRE GENERICO	PREMICRON ® B.BRAUN					
DESCRIPCIÓN	PremiCron® es una sutura quirúrgica sintética estéril compuesta de filamentos de polietilentereftalato (PET). Los filamentos de PET son trenzados y recubiertos de silicona según las especificaciones de calibre marcadas por USP. La sutura de poliéster PremiCron® está disponible en color blanco (sin teñir) o teñida de color verde (D&C green nº 6) y en diámetros comprendidos entre USP5 y USP8/0. Se dispone también de presentaciones de PremiCron® con pledgets. Los pledgets son 100% politetrafluoretileno (PTFE). PremiCron® cumple todos los requisitos establecidos por la Farmacopea Europea y la Farmacopea Americana (USP) para suturas no absorbibles estériles. Para mayor información ver inserto del producto.					
COMPOSICIÓN	Polietileno blanco o teñido de verde D & C green N° 6). Aguja en acero inoxidable.					
INDICACIONES DE USO	Diseñado para ser usado en cirugía general y en indicaciones cardiovasculares y vasculares. Cirugía cardiovascular o vasculares ortopedia.					
CONTRAINDICACIONES	No se han descrito.					
ADVERTENCIAS	PremiCron® no debe reesterilizarse. Desechar los envases abiertos, no utilizados o dañados. Un solo uso: peligro de infección para el paciente o el usuario y limitación de la función de los productos en caso de reutilización. La contaminación o la función limitada de los productos pueden provocar lesiones, enfermedades e incluso la muerte. Para mayor información ver inserto del producto.					
VIDA ÚTIL	5 AÑOS			CLASIFICACION DE RIESGO		III
REGISTRO SANITARIO	INVIMA 2021DM-0000118-R2		MODALIDAD	IMPORTAR Y VENDER	VIGENCIA	2-mar.-2031
TITULAR	B.BRAUN SURGICAL S.A		FABRICANTE	B.BRAUN SURGICAL S.A		
REFERENCIAS	VER ANEXO					
PRESENTACION COMERCIAL	Caja x 12, 24 y 36 unidades					
ENVASE	TIPO	MATERIAL	UNIDADES CONTENIDAS	INFORMACIÓN CONTENIDA		
	Blister	Papel grado medico Polietileno	1	Nombre.		
				Descripción del producto.		
				Lote, fecha de fabricación y expiración.		
				Titular		
EMPAQUE	TIPO	MATERIAL	UNIDADES CONTENIDAS	INFORMACIÓN CONTENIDA		
	Caja	Cartón Celofan	Caja x 12, 24 y 36 unidades	Nombre.		
				Descripción del producto.		
				Lote, fecha de fabricación y expiración.		
				Titular		
EMBALAJE	TIPO	MATERIAL	UNIDADES CONTENIDAS	INFORMACIÓN CONTENIDA		
	N/A	N/A	N/A	Nombre.		
				Descripción del producto.		
				Lote, fecha de fabricación y expiración.		
				Titular, Importador, registro sanitario.		
CONDICIONES ALMACENAMIENTO	Almacenar a condiciones normales de temperatura.					
MÉTODO DE DESECHO O DISPOSICIÓN FINAL	La disposición final del dispositivo una vez utilizado, se debe realizar de acuerdo a la normativa ambiental aplicable en materia de residuos peligrosos y residuos hospitalarios.					
BENEFICIOS AMBIENTALES	Los empaques y embalajes del producto poseen características de aprovechamiento, tras un adecuado proceso de separación en la fuente.					



ANEXO -FICHA TÉCNICA DISPOSITIVOS MEDICOS

Código	Descripción
B0026019	PREMICRON GRN 0(3,5)4X45 HR26S TO CVDDP
B0026020	PREMICRON GRN 0(3,5)4X45 HR26SS TOCVDDP
B0026044	PREMICRON GREEN 2 (5) 4X75CM HRT40(M)DDP
B0026053	PREMICRON GREEN 3 (6) 75CM HR26S CV DDP
B0026065	PREMICRON GREEN 2/0(3)4X45 HR48 TO(M)DDP
B0026075	PREMICRON GREEN 2 (5) 90CM HR40 (M) DDP
B0026082	PREMICRON GREEN 3 (6) 75CM HR48 (M) DDP
B0026083	PREMICRON GREEN 5 (7) 75CM HR40 (M) DDP
B0026097	PREMICRON GREEN 2 (5) 75CM HRT48 (M) DDP
B0026153	PREMICRON GREEN 5 (7) 45CM 2XHS40(M) DDP
B0026182	PREMICRON GREEN 2 (5) 75CM HR48 (M) DDP
B0026191	PREMICRON GREEN 5 (7) 75CM HR37SS (M)DDP
B0026205	PREMICRON GREEN 2 (5) 50CM DS76(M) DDP
B0026208	PREMICRON GREEN 5(7)50CM DS76(M)PEEL CE
B0026220	PREMICRON GRN 2(5)75CM 2XDS76(M)PEEL CE
B0026223	PREMICRON GRN 5(7)75CM 2XDS76(M)PEEL CE
B0026278	PREMICRON GREEN 2 (5) 75CM DS60(M) DDP
B0026279	PREMICRON GREEN 3 (6) 75CM DS60(M) DDP
B0026291	PREMICRON GREEN 5(7)75CM DS76(M)PEEL CE
B0026303	PREMICRON GREEN 3 (6) 75CM HR65(M) DDP
B0026325	PREMICRON GREEN 3 (6)75CM 2XHR26SSCVDDP
B0026333	PREMICRON GREEN 3 (6) 75CM DS39(M) DDP
B0026355	PREMICRON GREEN 2 (5) 75CM DS48(M) DDP
B0026376	PREMICRON GREEN 2 (5) 75CM HS48(M) DDP
B0026377	PREMICRON GREEN 3 (6) 75CM HS48(M) DDP
B0026399	PREMICRON GREEN 5 (7) 75CM HS48(M) DDP
B0026409	PREMICRON GREEN 5 (7) 75CM HS60(M) DDP
B0026425	PREMICRON GREEN 3 (6) 75CM HR40(M) DDP
B0026481	PREMICRON GREEN 3 (6) 75CM HS40(M) DDP
B0026482	PREMICRON GREEN 5 (7) 75CM HS40(M) DDP
B0026509	PREMICRON GREEN 2 (5) 75CM GS51(M) DDP
B0026513	PREMICRON GREEN 2 (5) 75CM HR40S(M) DDP
B0026545	PREMICRON GREEN 3 (6) 75CM HS30(M) DDP
B0026696	PREMICRON GREEN 0(3,5)75CM 2XHR70 (M)DDP
B0026697	PREMICRON GREEN 1 (4) 75CM 2XHR70 (M)DDP
B0026717	PREMICRON GREEN 2 (5) 4X75CM HRT37(M)DDP
B0026718	PREMICRON GREEN 2 (5) 4X75CM HRT43(M)DDP
B0026720	PREMICRON GREEN 2 (5) 75CM 2XHR48 (M)DDP
B0026721	PREMICRON GREEN 0(3,5)4X45 HR48 TO(M)DDP
B0026727	PREMICRON GREEN 0(3,5)4X45 HR37S TO(MDDP
B0026729	PREMICRON GRN 0(3,5)4X45 HR22SS TOCVDDP
B0026759	PREMICRON GREEN 2 (5) 90CM 2XHR37S(M)DDP
B0026766	PREMICRON GREEN 2 (5) 4X75CM HRS40 DDP
B0026767	PREMICRON GREEN 2 (5) 4X75CM HRS48 DDP
B0026778	PREMICRON GREEN 2 (5) 75CM HS26S(M) DDP
B0026880	PREMICRON GREEN 2 (5) 4X75CM HRT48(M)DDP
B0026898	PREMICRON GRN 0(3,5)90CM2XHRC22BCVPF6DDP
B0026899	PREMICRON GRN 0(3,5)90CM2XHRC22(M)PF6DDP
B0027082	PREMICRON WHITE 3 (6) 75CM HR48 (M) DDP
B0027148	PREMICRON G/W 2 (5) 4X75CM HR26SS CVDDP
B0027642	PREMICRON G/W 2 (5) 2X60CM HR22SS CVDDP
B0027718	PREMICRON G/W 2/0 2X75 2XHR22 CVPF6 DDP
B0027875	PREMICRON G/W 2/0(3)2X90CM2XHRC17MPF6DDP
B0027885	PREMICRON G/W 2/0(3)2X90CM2XHRC26MPF6DDP
B0027978	PREMICRON G/W 2 (5) 4X90CM HRC22 (M) DDP
B0102030	PREMICRON WHITE 2 (5) 5X45CM DDP
B0102035	PREMICRON WHITE 3 (6) 150CM DDP
B0120009	PREMICRON GREEN 0 (3,5) 10X45CM DDP
B0120010	PREMICRON GREEN 1 (4) 10X45CM DDP
B0120011	PREMICRON GREEN 2 (5) 10X45CM DDP
B0120022	PREMICRON GREEN 0 (3,5) 10X75CM DDP
B0120024	PREMICRON GREEN 2 (5) 10X75CM DDP
B0120031	PREMICRON GREEN 3 (6) 150CM DDP
B0120032	PREMICRON GREEN 4 (6) 150CM DDP
B0120033	PREMICRON GREEN 5 (7) 150CM DDP
B0120070	PREMICRON GREEN 1 (4) 200CM HS48(M) DDP
B0120085	PREMICRON GREEN 3 (6) 2X75CM DDP
B0120087	PREMICRON GREEN 5 (7) 2X75CM DDP
B0120105	PREMICRON GREEN 1 (4) 5X45CM DDP
B0120106	PREMICRON GREEN 2 (5) 5X45CM DDP
B0120111	PREMICRON GREEN 1 (4) 5X60CM DDP
B0120122	PREMICRON GREEN 1 (4) 5X75CM DDP

B0120123	PREMICRON GREEN 2 (5) 5X75CM DDP
C0026001	PREMICRON GREEN 0 (3,5) 120CM HR26CVDDP
C0026002	PREMICRON GREEN 3/0 (2) 120CM HR26CVDDP
C0026003	PREMICRON GREEN 5/0 (1) 75CM HR17 CVDDP
C0026004	PREMICRON GREEN 4/0(1,5)75CM HR17 CVDDP
C0026005	PREMICRON GREEN 3/0 (2) 75CM HR17 CVDDP
C0026006	PREMICRON GREEN 2/0 (3) 75CM HR17 CVDDP
C0026007	PREMICRON GREEN 0 (3,5) 45CM HR22 CVDDP
C0026008	PREMICRON GREEN 6/0(0,7)45CM HR10 CVDDP
C0026009	PREMICRON GREEN 5/0 (1) 45CM HR10 CVDDP
C0026010	PREMICRON GREEN 4/0(1,5)45CM HR10 CVDDP
C0026011	PREMICRON GREEN 2 (5) 90 2XHS23SS (M)DDP
C0026014	PREMICRON GREEN 4/0(1,5)75CM HR26 CVDDP
C0026015	PREMICRON GREEN 3/0 (2) 75CM HR26 CVDDP
C0026016	PREMICRON GREEN 2/0 (3) 75CM HR26 CVDDP
C0026017	PREMICRON GREEN 0 (3,5) 75CM HR26 CVDDP
C0026018	PREMICRON GREEN 1 (4) 75CM HR26 CVDDP
C0026024	PREMICRON GREEN 4/0(1,5)75CM HR22 CVDDP
C0026025	PREMICRON GREEN 3/0 (2) 75CM HR22 CVDDP
C0026026	PREMICRON GREEN 2/0 (3) 75CM HR22 CVDDP
C0026029	PREMICRON GREEN 5/0 (1) 45CM HR13 CVDDP
C0026030	PREMICRON GREEN 4/0(1,5)45CM HR13 CVDDP
C0026031	PREMICRON GREEN 3/0 (2) 45CM HR13 CVDDP
C0026032	PREMICRON GREEN 3/0 (2) 75CM HR13 CVDDP
C0026035	PREMICRON GREEN 3/0 (2) 75CM HR37 (M)DDP
C0026036	PREMICRON GREEN 2/0 (3) 75CM HR37 (M)DDP
C0026037	PREMICRON GREEN 0 (3,5) 75CM HR37 (M)DDP
C0026038	PREMICRON GREEN 1 (4) 75CM HR37 (M)DDP
C0026039	PREMICRON GREEN 2 (5) 75CM HR37 (M) DDP
C0026041	PREMICRON GREEN 3/0 (2) 45CM GR19 (M)DDP
C0026045	PREMICRON GREEN 2/0 (3) 120CM GR19(M)DDP
C0026046	PREMICRON GREEN 2/0 (3)120CM SKR26(M)DDP
C0026049	PREMICRON GREEN 2/0 (3) 90CM HR26 CVDDP
C0026050	PREMICRON GREEN 0 (3,5) 75CM HR26SCVDDP
C0026051	PREMICRON GREEN 1 (4) 75CM HR26S CVDDP
C0026052	PREMICRON GREEN 2/0 (3) 75CM HR26SCVDDP
C0026054	PREMICRON GREEN (2,5) 75CM HR30 (M) DDP
C0026056	PREMICRON GREEN 3/0 (2) 75CM HR30 (M)DDP
C0026057	PREMICRON GREEN 2/0 (3) 75CM HR30 (M)DDP
C0026058	PREMICRON GREEN 0 (3,5) 75CM HR30 (M)DDP
C0026059	PREMICRON GREEN 1 (4) 75CM HR30 (M)DDP
C0026061	PREMICRON GREEN 3/0 (2) 75CM HR26SCVDDP
C0026062	PREMICRON GREEN 2/0 (3)75CM HR13SSCVDDP
C0026064	PREMICRON GRN 2/0(3)4X45 HR26SS TOCVDDP
C0026066	PREMICRON GREEN 2/0 (3) 75CM HR37S(M)DDP
C0026067	PREMICRON GREEN 0 (3,5) 75CM HR37S(M)DDP
C0026068	PREMICRON GREEN 1 (4) 75CM HR37S (M)DDP
C0026069	PREMICRON GREEN 2 (5) 75CM HR37S (M) DDP
C0026070	PREMICRON GREEN 2 (5)75CM 2XHR22SSCVDDP
C0026071	PREMICRON GREEN 0 (3,5)75CM 2XDR48(M)DDP
C0026072	PREMICRON GREEN 2/0 (3) 75CM HR40 (M)DDP
C0026074	PREMICRON GREEN 1 (4) 75CM HR40 (M)DDP
C0026079	PREMICRON GREEN 2 (5) 75CM HR37SS (M)DDP
C0026080	PREMICRON GREEN 1 (4) 75CM HR48 (M)DDP
C0026086	PREMICRON GREEN 2/0(3)100CM HR22 CV DDP
C0026088	PREMICRON GREEN 1 (4) 75CM FR26 (M) DDP
C0026089	PREMICRON GREEN 1 (4) 90CM FR36 (M) DDP
C0026096	PREMICRON GREEN 0 (3,5) 100CM HR26CVDDP
C0026098	PREMICRON GREEN 0 (3,5) 75CM HR48 (M)DDP
C0026099	PREMICRON GREEN 0(3,5)120CM PRC24(M)DDP
C0026102	PREMICRON GREEN 2/0(3)120CM PRC24(M)DDP
C0026104	PREMICRON GREEN 4/0(1,5)75CM DR18 CVDDP
C0026105	PREMICRON GREEN 3/0 (2) 75CM DR18 CVDDP
C0026111	PREMICRON GREEN 2/0 (3) 120CM HR26CVDDP
C0026113	PREMICRON GREEN 4/0(1,5)75CM GR19 (M)DDP
C0026119	PREMICRON GREEN 2/0 (3) 20CM GR35 (M)DDP
C0026120	PREMICRON GREEN 6/0(0,7)45CM DR12 CVDDP
C0026121	PREMICRON GREEN 5/0 (1) 45CM DR12 CVDDP
C0026130	PREMICRON GREEN 4/0(1,5)75CM DR15 CVDDP
C0026131	PREMICRON GREEN 3/0 (2) 75CM DR15 CVDDP
C0026136	PREMICRON GREEN 3/0 (2) 75CM DR26 CVDDP
C0026137	PREMICRON GREEN 2/0 (3) 75CM DR26 CVDDP
C0026139	PREMICRON GREEN 2/0(3)100CM DR26 CV DDP
C0026141	PREMICRON GREEN 2/0 (3) 75CM DR32 (M)DDP
C0026142	PREMICRON GREEN 0 (3,5) 75CM DR32 (M)DDP
C0026143	PREMICRON GREEN 2/0 (3)75CM 2XDR22CVDDP
C0026145	PREMICRON GREEN 2 (5) 75CM HR30 (M)DDP
C0026152	PREMICRON GREEN 4/0(1,5)75CM HR13 CVDDP

C0026156	PREMICRON GREEN 4/0(1,5)75CM HRC22CV DDP
C0026158	PREMICRON GREEN 0(3,5)75CM 2XHRC17CVDDP
C0026159	PREMICRON GREEN 0(3,5)75CM 2XHRC17BCVDDP
C0026161	PREMICRON GRN 2/0(3)75CM2XDRC18BCVPF6DDP
C0026162	PREMICRON GRN 0(3,5)75CM 2XHRC26CVPF6DDP
C0026163	PREMICRON GREEN 0(3,5)90CM 2XHRC22 CVDDP
C0026165	PREMICRON GRN 0(3,5)90CM 2XHRC26CVPF6DDP
C0026168	PREMICRON GREEN 0(3,5)90CM 2XHRC26BCVDDP
C0026169	PREMICRON GRN0(3,5)90CM 2XHRC26BCVPF6DDP
C0026171	PREMICRON GREEN 6/0 (0,7)45CM DS12(M)DDP
C0026172	PREMICRON GREEN 5/0 (1) 45CM DS12(M) DDP
C0026173	PREMICRON GREEN 4/0 (1,5)45CM DS12(M)DDP
C0026174	PREMICRON GRN 2/0(3)90CM2XHRC22BCVPF6DDP
C0026176	PREMICRON GREEN 5/0 (1) 45CM DS16(M) DDP
C0026177	PREMICRON GREEN 4/0 (1,5)45CM DS16(M)DDP
C0026179	PREMICRON GREEN 2/0(3)90CM 2XHRC22BCVDDP
C0026180	PREMICRON GREEN 2/0(3)90CM 2XHRC26BCVDDP
C0026181	PREMICRON GRN 2/0(3)90CM2XHRC26BCVPF6DDP
C0026184	PREMICRON GRN 2/0(3)75CM2XDRC18BCVPS6DDP
C0026185	PREMICRON GREEN 4/0(1,5)4X45CM DS19(MDDP
C0026189	PREMICRON GREEN 3/0 (2) 45CM DS24(M) DDP
C0026190	PREMICRON GREEN 2/0 (3) 45CM DS24(M) DDP
C0026192	PREMICRON GREEN 6/0 (0,7) 45CM DS9(M)DDP
C0026193	PREMICRON GREEN 4/0(1,5)45CM DSMP13(MDDP
C0026194	PREMICRON GREEN 1 (4) 120CM HR26 (M) DDP
C0026196	PREMICRON GREEN 5/0 (1) 45CM HS15(M) DDP
C0026197	PREMICRON GREEN 4/0 (1,5)45CM HS15(M)DDP
C0026198	PREMICRON GREEN 3/0 (2) 45CM HS15(M) DDP
C0026199	PREMICRON GREEN 4/0 (1,5)45CM HS21(M)DDP
C0026200	PREMICRON GREEN 3/0 (2) 45CM HS21(M) DDP
C0026202	PREMICRON GREEN 3/0 (2) 45CM HS23(M) DDP
C0026204	PREMICRON GREEN 3/0 (2) 75CM DR32(M) DDP
C0026211	PREMICRON GRN 2/0(3)75CM 2XHR17B CVDDP
C0026213	PREMICRON GREEN 5/0(1)75CM 2XHR10FCVDDP
C0026214	PREMICRON GREEN 2 (5) 75CM 2XHS30 (M)DDP
C0026215	PREMICRON GREEN 3/0(2)75CM 2XHRC26CVDDP
C0026226	PREMICRON GREEN 3/0 (2)75CM DR22 CV DDP
C0026227	PREMICRON GREEN 2/0 (3)75CM DR22 CV DDP
C0026228	PREMICRON GREEN 0 (3,5)75CM DR22 CV DDP
C0026231	PREMICRON GREEN 4/0 (1,5)75CM DS19(M)DDP
C0026232	PREMICRON GREEN 3/0 (2) 75CM DS19(M) DDP
C0026233	PREMICRON GREEN 2/0 (3) 75CM DS19(M) DDP
C0026234	PREMICRON GREEN 0 (3,5) 75CM DS19(M) DDP
C0026235	PREMICRON GREEN 1 (4) 75CM DS19(M) DDP
C0026238	PREMICRON GREEN 1 (4) 75CM DR32(M) DDP
C0026239	PREMICRON GREEN 2 (5) 75CM DR32(M) DDP
C0026240	PREMICRON GRN 3/0(2)75CM 2XHRC17CVPF6DDP
C0026242	PREMICRON GREEN 5/0 (1) 75CM DS12(M) DDP
C0026245	PREMICRON GREEN 4/0 (1,5)75CM DS16(M)DDP
C0026246	PREMICRON GREEN 3/0 (2) 75CM DS16(M) DDP
C0026248	PREMICRON G 2/0(3)75CM 2XHRC17CVPF6 DDP
C0026249	PREMICRON GRN 2/0(3)75CM 2XHRC17CVPF3DDP
C0026250	PREMICRON GREEN 2/0(3)75CM 2XHRC17CVDDP
C0026251	PREMICRON GRN 0(3,5)75CM 2XHRC17(MPF6DDP
C0026252	PREMICRON GRN 2/0(3)75CM 2XHR17BCVPF3DDP
C0026253	PREMICRON GR 0(3,5)75CM 2XHRC17BCVPF6DDP
C0026255	PREMICRON GREEN 2 (5) 90CM HR26SS CVDDP
C0026257	PREMICRON GRN 2/0(3)75CM 2XHR17BCVPS3DDP
C0026260	PREMICRON GREEN 2 (5) 75CM HRC40 (M) DDP
C0026263	PREMICRON G 2/0(3)75CM 2XHRC22CVPF6 DDP
C0026265	PREMICRON GRN 2/0(3)75CM 2XHR22B CVDDP
C0026266	PREMICRON GREEN 0 (3,5) 90CM GS60 (M)
C0026267	PREMICRON GREEN 0 (3,5) 90CM HR37S(M)DDP
C0026268	PREMICRON GREEN 0 (3,5) 90CM FR36 (M)DDP
C0026275	PREMICRON GREEN 2 (5) 75CM HR40SS (M)DDP
C0026277	PREMICRON GREEN 1 (4) 90CM HR37S(M) DDP
C0026281	PREMICRON GREEN 2/0 (3) 75CM DS35(M) DDP
C0026282	PREMICRON GREEN 0 (3,5) 75CM DS35(M) DDP
C0026283	PREMICRON GREEN 1 (4) 75CM DS35(M) DDP
C0026284	PREMICRON GREEN 2 (5) 75CM DS35(M) DDP
C0026286	PREMICRON GREEN 3/0 (2) 75CM GS65(M) DDP
C0026287	PREMICRON GREEN 2/0 (3) 75CM GS65(M) DDP
C0026288	PREMICRON GREEN 0 (3,5) 75CM GS65(M) DDP
C0026289	PREMICRON GREEN 1 (4) 75CM GS65(M) DDP
C0026296	PREMICRON GREEN 2/0(3)100CM HRC22 CVDDP
C0026299	PREMICRON GREEN 3/0(2)100CM HRC26 CVDDP
C0026300	PREMICRON GREEN 2/0(3)100CM HRC26 CVDDP
C0026301	PREMICRON GREEN 3/0 (2) 75CM HRC26CVDDP

C0026302	PREMICRON GREEN 2/0 (3) 75CM HRC26CVDDP
C0026304	PREMICRON GR 4/0(1,5)90CM 2XHRC17 CVDDP
C0026305	PREMICRON GRN 3/0(2)90CM 2XHRC17 CV DDP
C0026306	PREMICRON GRN 2/0(3)90CM 2XHRC17 CV DDP
C0026307	PREMICRON GREEN 3/0(2)75 2XHR22CVPF6DDP
C0026308	PREMICRON GRN 2/0(3)90CM 2XHRC17CVPF6DDP
C0026312	PREMICRON GRN 3/0(2)90CM 2XHRC26 CV DDP
C0026313	PREMICRON GRN 2/0(3)90CM 2XHRC26 CV DDP
C0026314	PREMICRON GREEN 3/0 (2) 75CM GS60 (M)DDP
C0026315	PREMICRON GREEN 2/0 (3) 75CM GS60(M) DDP
C0026316	PREMICRON GRN 2/0(3)75CM 2XHRC26 CV DDP
C0026317	PREMICRON GREEN 0(3,5)75CM 2XHRC26CVDDP
C0026318	PREMICRON G 2/0(3)90CM 2XHRC26CVPF6 DDP
C0026319	PREMICRON G 2/0(3)75CM 2XHRC26CVPF6 DDP
C0026320	PREMICRON GREEN 0(3,5)75CM 2XHRC26BCVDDP
C0026321	PREMICRON GREEN 2/0(3)90CM2XHRC22MPF6DDP
C0026323	PREMICRON GRN 2/0(3)90CM 2XHRC22 CV DDP
C0026328	PREMICRON GREEN 2/0 (3)75CM 2XDR48(M)DDP
C0026329	PREMICRON GREEN 2/0 (3) 75CM DS39(M) DDP
C0026330	PREMICRON GREEN 0 (3,5) 75CM DS39(M) DDP
C0026331	PREMICRON GREEN 1 (4) 75CM DS39(M) DDP
C0026332	PREMICRON GREEN 2 (5) 75CM DS39(M) DDP
C0026337	PREMICRON GREEN 0 (3,5) 75CM HRT26(M)DDP
C0026339	PREMICRON GREEN 1 (4) 75CM HS26S(M) DDP
C0026340	PREMICRON GREEN 0 (3,5) 90CM DS24 (M)DDP
C0026341	PREMICRON GREEN 4/0 (1,5)75CM HS15(M)DDP
C0026342	PREMICRON GREEN 3/0 (2) 75CM HS15(M) DDP
C0026352	PREMICRON GREEN 6/0(0,7)75CM DRC12CVDDP
C0026353	PREMICRON GREEN 0 (3,5) 75CM DS48(M) DDP
C0026354	PREMICRON GREEN 1 (4) 75CM DS48(M) DDP
C0026365	PREMICRON GREEN 3/0 (2) 75CM HS21(M) DDP
C0026366	PREMICRON GREEN 2/0 (3) 75CM HS21(M) DDP
C0026367	PREMICRON GREEN 0 (3,5) 75CM HS21(M) DDP
C0026368	PREMICRON GREEN 6/0 (0,7) 45CM DRTM8 DDP
C0026371	PREMICRON GREEN 2 (5) 75CM HS26SS(M) DDP
C0026373	PREMICRON GREEN 2 (5) 75CM HS37SS(M) DDP
C0026375	PREMICRON GREEN 1 (4) 75CM HS48(M) DDP
C0026390	PREMICRON GREEN 2 (5) 75CM HS40S(M) DDP
C0026392	PREMICRON GREEN 0 (3,5) 75CM HRT37(M)DDP
C0026393	PREMICRON GREEN 1 (4) 75CM HRT37 (M)DDP
C0026394	PREMICRON GREEN 2 (5) 75CM HRT37 (M) DDP
C0026395	PREMICRON GREEN 2 (5) 75CM HRT40S (M)DDP
C0026397	PREMICRON GREEN 1 (4) 75CM HS40SS(M) DDP
C0026401	PREMICRON GREEN 0 (3,5) 75CM HRT40(M)DDP
C0026403	PREMICRON GREEN 0 (3,5)75CM HS23SS(M)DDP
C0026404	PREMICRON GREEN 1 (4) 75CM HS23SS(M) DDP
C0026405	PREMICRON GREEN 2 (5) 75CM HS23SS(M) DDP
C0026419	PREMICRON GREEN 0(3,5)90CM 2XHRC26 CVDDP
C0026420	PREMICRON GREEN 3/0 (2) 75CM HS26(M) DDP
C0026421	PREMICRON GREEN 2/0 (3) 75CM HS26(M) DDP
C0026422	PREMICRON GRN 2/0(3)75CM 2XHRC26CVPF3DDP
C0026423	PREMICRON GREEN 0 (3,5) 75CM HS26 (M)DDP
C0026441	PREMICRON GREEN 2/0 (3) 90CM GS60(M) DDP
C0026449	PREMICRON GRN 4/0(1,5)75CM 2XGR19 (M)DDP
C0026458	PREMICRON GREEN 2/0(3)90CM 2XHRT37(M)DDP
C0026460	PREMICRON GRN 2/0(3)75CM 2XHR22BCVPF6DDP
C0026465	PREMICRON GREEN 3/0(2)75CM 2XDRC18CVDDP
C0026466	PREMICRON GREEN 2/0(3)75CM 2XDRC18CVDDP
C0026467	PREMICRON GREEN 0 (3,5) 75CM HS37S(M)DDP
C0026468	PREMICRON GREEN 1 (4) 75CM HS37S(M) DDP
C0026469	PREMICRON GREEN 2 (5) 75CM HS37S(M) DDP
C0026471	PREMICRON GREEN 0 (3,5) 75CM HS37(M) DDP
C0026472	PREMICRON GREEN 1 (4) 75CM HS37(M) DDP
C0026473	PREMICRON GREEN 2 (5) 75CM HS37(M) DDP
C0026476	PREMICRON GRN 2/0(3)75CM 2XDRC18CVPF3DDP
C0026477	PREMICRON GREEN 4/0(1,5)4X45 HR13 TOMDDP
C0026479	PREMICRON GREEN 1 (4) 75CM HS40(M) DDP
C0026480	PREMICRON GREEN 2 (5) 75CM HS40(M) DDP
C0026493	PREMICRON GRN 2/0(3)75CM 2XHR22BCVPS6DDP
C0026500	PREMICRON GRN 2/0(3)75CM 2XHR26B CVDDP
C0026505	PREMICRON GREEN 3/0 (2) 75CM GS51(M) DDP
C0026506	PREMICRON GREEN 2/0 (3) 75CM GS51(M) DDP
C0026507	PREMICRON GREEN 0 (3,5) 75CM GS51(M) DDP
C0026508	PREMICRON GREEN 1 (4) 75CM GS51(M) DDP
C0026512	PREMICRON GREEN 1 (4) 75CM HR40S(M) DDP
C0026514	PREMICRON GREEN 4/0 (1,5)75CM DS24(M)DDP
C0026515	PREMICRON GREEN 3/0 (2) 75CM DS24(M) DDP
C0026516	PREMICRON GREEN 2/0 (3) 75CM DS24(M) DDP

C0026517	PREMICRON GREEN 0 (3,5) 75CM DS24(M) DDP
C0026518	PREMICRON GREEN 1 (4) 75CM DS24(M) DDP
C0026522	PREMICRON GRN 2/0(3)75CM 2XHR26BCVPF6DDP
C0026527	PREMICRON GRN 2/0(3)75CM 2XHR26BCVPS6DDP
C0026531	PREMICRON GREEN 2/0 (3) 75CM DS30(M) DDP
C0026532	PREMICRON GREEN 0 (3,5) 75CM DS30(M) DDP
C0026533	PREMICRON GREEN 1 (4) 75CM DS30(M) DDP
C0026534	PREMICRON GREEN 2 (5) 75CM DS30(M) DDP
C0026535	PREMICRON GREEN 3/0 (2) 30CM SKR19(M)DDP
C0026536	PREMICRON GREEN 2/0 (3) 30CM SKR19(M)DDP
C0026537	PREMICRON GREEN 0 (3,5) 90CM DS30(M) DDP
C0026538	PREMICRON GRN 2/0(3)75CM 2XHRC17BCVDDP
C0026540	PREMICRON GREEN 3/0(2)30CM SKR26(M)DDP
C0026541	PREMICRON GREEN 2/0 (3) 75CM HS30(M) DDP
C0026542	PREMICRON GREEN 0 (3,5) 75CM HS30(M) DDP
C0026543	PREMICRON GREEN 1 (4) 75CM HS30(M) DDP
C0026544	PREMICRON GREEN 2 (5) 75CM HS30(M) DDP
C0026546	PREMICRON GREEN 3/0 (2)120CM SKR26(M)DDP
C0026548	PREMICRON GREEN 2/0(3)45CM SKR26 (M) DDP
C0026549	PREMICRON GREEN 0(3,5)120CM SKR26 (M)DDP
C0026550	PREMICRON GREEN 2/0 (3) 30CM SKR26(M)DDP
C0026551	PREMICRON GREEN 5/0 (1) 45CM DS19(M) DDP
C0026552	PREMICRON GREEN 4/0 (1,5)45CM DS19(M)DDP
C0026553	PREMICRON GREEN 3/0 (2) 45CM DS19(M) DDP
C0026554	PREMICRON GREEN 2/0 (3) 45CM DS19(M) DDP
C0026556	PREMICRON GRN 2/0(3)75CM2XHRC17BCVPF3DDP
C0026560	PREMICRON GRN 2/0(3)75CM2XHRC17CVPS3DDP
C0026563	PREMICRON GRN 2/0(3)75CM2XHRC17CVPS6DDP
C0026568	PREMICRON GRN 2/0(3)75CM2XHRC17BCVPS3DDP
C0026579	PREMICRON GRN 2/0(3)75CM2XHRC17BCVPF6DDP
C0026585	PREMICRON GRN 2/0(3)75CM2XHRC22CVPS6DDP
C0026587	PREMICRON GRN 2/0(3)75CM2XHRC17BCVPS6DDP
C0026592	PREMICRON GRN 2/0(3)75CM 2XHRC18BCVDDP
C0026594	PREMICRON GRN 2/0(3)75CM2XHRC26CVPS6DDP
C0026596	PREMICRON GRN 2/0(3)75CM2XHRC18BCVPF3DDP
C0026600	PREMICRON GRN 2/0(3)75CM2XHRC18BCVPS3DDP
C0026606	PREMICRON GREEN 2/0 (3)45CM 2XHR17CVDDP
C0026612	PREMICRON GRN 2/0(3)75CM2XHRC18BCVPS6DDP
C0026616	PREMICRON GRN 2/0(3)75CM 2XHRC22BCVDDP
C0026620	PREMICRON GRN 2/0(3)75CM2XHRC22BCVPF6DDP
C0026625	PREMICRON GRN 2/0(3)75CM2XHRC22BCVPS6DDP
C0026634	PREMICRON GRN 2/0(3)75CM 2XHRC26BCVDDP
C0026636	PREMICRON GRN 2/0(3)75CM 2XDRC18CVPS3DDP
C0026643	PREMICRON GRN 2/0(3)75CM2XHRC26BCVPF6DDP
C0026644	PREMICRON GRN 0(3,5)75CM2XHRC26BCVPF6DDP
C0026650	PREMICRON GRN 2/0(3)75CM2XHRC26BCVPS6DDP
C0026658	PREMICRON GRN 2/0(3)75CM 2XHRC27BCVDDP
C0026665	PREMICRON GRN 2/0(3)75CM2XHRC27BCVPF6DDP
C0026674	PREMICRON GRN 2/0(3)75CM2XHRC27BCVPS6DDP
C0026676	PREMICRON GRN 2/0(3)90CM 2XHRC26CVPS6DDP
C0026689	PREMICRON GREEN 0(3,5)75CM HRT30S (M)DDP
C0026695	PREMICRON GREEN 5/0 (1) 75CM 2XHR10CVDDP
C0026700	PREMICRON GRN 4/0(1,5)4X75 HR17 TOCVDDP
C0026704	PREMICRON GRN 4/0(1,5)4X45CM HR17 TOMDDP
C0026724	PREMICRON GRN 2/0 (3)4X45 HR26S TOCVDDP
C0026730	PREMICRON GRN 4/0(1,5)90CM 2XHR22 CVDDP
C0026731	PREMICRON GREEN 3/0(2)90CM 2XHR22 CVDDP
C0026732	PREMICRON GREEN 2/0(3)90CM 2XHR22 CVDDP
C0026733	PREMICRON GRN 2/0(3)75CM 2XHR22 CVPF6DDP
C0026735	PREMICRON G 4/0(1,5)75 2XHR22 CVPF6 DDP
C0026736	PREMICRON GRN 2/0(3)4X45 HR22SS TOCVDDP
C0026737	PREMICRON GREEN 2/0(3)75CM 2XHR22 CVDDP
C0026738	PREMICRON GREEN 3/0(2)75CM 2XHR26 CVDDP
C0026739	PREMICRON GREEN 2/0(3)75CM 2XHR26 CVDDP
C0026740	PREMICRON GRN 3/0(2)4X45CM HR17 TO(M)DDP
C0026750	PREMICRON GREEN 3/0(2)90CM 2XHR37 (M)DDP
C0026751	PREMICRON GREEN 2/0(3)90CM 2XHR37 (M)DDP
C0026752	PREMICRON GREEN 0(3,5)90CM 2XHR37 (M)DDP
C0026763	PREMICRON GRN 4/0(1,5)75CM 2XHR22 CVDDP
C0026768	PREMICRON GREEN 0(3,5)2X75CM HR26SS(MDDP
C0026781	PREMICRON GREEN 2 (5) 75CM HRT43S (M)DDP
C0026792	PREMICRON GRN 5/0(1)90 2XHR10 CVPF3 DDP
C0026802	PREMICRON GREEN 6/0(0,7)75CM 2XHR17CVDDP
C0026803	PREMICRON GREEN 5/0(1)75CM 2XHR17 CVDDP
C0026804	PREMICRON GRN 4/0(1,5)75CM 2XHR17 CVDDP
C0026805	PREMICRON GREEN 3/0(2)75CM 2XHR17 CVDDP
C0026806	PREMICRON GREEN 2/0(3)75CM 2XHR17 CVDDP
C0026814	PREMICRON GRN 4/0(1,5)90CM 2XHR26 CVDDP

C0026815	PREMICRON GREEN 3/0(2)90CM 2XHR26 CVDDP
C0026816	PREMICRON GREEN 2/0(3)90CM 2XHR26 CVDDP
C0026817	PREMICRON GREEN 0(3,5)90CM 2XHR26 CVDDP
C0026822	PREMICRON GREEN 5/0(1)45CM 2XDR8 CVDDP
C0026836	PREMICRON GREEN 2/0(3)75CM 2XHR37 (M)DDP
C0026837	PREMICRON GREEN 0(3,5)75CM 2XHR37 (M)DDP
C0026840	PREMICRON GREEN 2(5)60CM 2XHR37S (M)DDP
C0026843	PREMICRON GREEN 5/0(1)75CM 2XHR13 CVDDP
C0026844	PREMICRON GRN 4/0(1,5)75CM 2XHR13 CVDDP
C0026845	PREMICRON GREEN 3/0(2)75CM 2XHR13 CVDDP
C0026852	PREMICRON GRN 6/0(0,7)75CM 2XDR12 CVDDP
C0026853	PREMICRON GREEN 5/0(1)75CM 2XDR12 CVDDP
C0026854	PREMICRON GRN 4/0(1,5)75CM 2XDR12 CVDDP
C0026858	PREMICRON GREEN 3/0(2)90CM 2XDR18 CVDDP
C0026859	PREMICRON GREEN 3/0 (2)75CM 2XDR12CVDDP
C0026863	PREMICRON GREEN 5/0(1)75CM 2XDR18 CVDDP
C0026864	PREMICRON GRN 4/0(1,5)75CM 2XDR18 CVDDP
C0026865	PREMICRON GREEN 3/0(2)75CM 2XDR18 CVDDP
C0026871	PREMICRON GRN 6/0(0,7)75CM 2XDR10 CVDDP
C0026873	PREMICRON GREEN 2/0 (3)90CM 2XDR26CVDDP
C0026883	PREMICRON GREEN 2/0 (3)90CM 2XDR18CVDDP
C0026887	PREMICR GRN 5/0(1)75CM 2XHR10 CVPF3 DDP
C0026888	PREMICR GRN 4/0(1,5)75CM2XHR10CVPF3 DDP
C0026903	PREMICRON GREEN 5/0(1)90CM 2XHR17 CVDDP
C0026904	PREMICRON GRN 4/0(1,5)90CM 2XHR17 CVDDP
C0026905	PREMICRON GREEN 3/0(2)90CM 2XHR17 CVDDP
C0026906	PREMICRON GREEN 2/0(3)90CM 2XHR17 CVDDP
C0026907	PREMICRON GR 2/0(3)75CM 2XHR20CVPF6 DDP
C0026915	PREMICRON GRN 3/0(2)90CM 2XHR26 CVPF6DDP
C0026916	PREMICRON GR 2/0(3)90CM 2XHR26CVPF6 DDP
C0026920	PREMICRON GREEN 1 (4) 90CM HRT30S (M)DDP
C0026928	PREMICRON GREEN 2/0 (3) 90CM HR30 (M)DDP
C0026936	PREMICRON GREEN 5/0 (1)75CM HRC17 CV DDP
C0026937	PREMICRON GREEN 4/0(1,5)75CM HRC17 CVDDP
C0026938	PREMICRON GREEN 3/0(2)75CM HRC17CVDDP
C0026949	PREMICRON GRN 5/0(1)75CM 2XHR13FF CV DDP
C0026954	PREMICRON G 4/0(1,5)75 2XDR12 CVPF3 DDP
C0026972	PREMICRON GREEN 1 (4) 90CM HR40 (M)DDP
C0026977	PREMICRON GREEN 1 (4) 75CM JR30S (M) DDP
C0026987	PREMICRON GREEN 1 (4) 75CM JRC30S (M)DDP
C0026992	PREMICRON GREEN 0 (3,5) 120CM GR19(M)DDP
C0027010	PREMICRON WHITE 2/0(3)75 2XHR20CVPF6DDP
C0027031	PREMICRON WHITE 5/0(1)45CM 2XDSMP11(MDDP
C0027034	PREMICRON WHITE 2/0(3)75 2XHR22CVPF6DDP
C0027059	PREMICRON WHITE 1 (4) 75CM HR30 (M)DDP
C0027074	PREMICRON WHITE 6/0(0,7)45CM2XDSMP7(MDDP
C0027086	PREMICRON WHITE 5/0 (1) 45CM DSMP13(MDDP
C0027087	PREMICRON WHITE 4/0(1,5)45CM DSMP13(MDDP
C0027096	PREMICRON WHITE 4/0(1,5)45CM DSMP19(MDDP
C0027097	PREMICRON WHITE 3/0 (2)45CM DSMP19(M)DDP
C0027306	PREMICRON WH 2/0(3)90CM 2XHR26CVPF6 DDP
C0027315	PREMICRON WHT 2/0(3)90CM 2XHRC17CVPF6DDP
C0027317	PREMICRON WHITE 2/0(3)90CM 2XHRC17CVDDP
C0027320	PREMICRON WHT 2/0(3)75CM2XHRC17 CVPS3DDP
C0027326	PREMICRON W 2/0(3)90CM 2XHRC26CV DDP
C0027363	PREMICRON WHITE 2/0(3)90CM 2XHRC22CVDDP
C0027433	PREMICRON WHITE 3/0(2)4X45CM DS19 (M)DDP
C0027455	PREMICRON WHITE 4/0(1,5)75CM 2XHR17CVDDP
C0027457	PREMICRON WHITE 2/0 (3)75CM 2XHR17CVDDP
C0027536	PREMICRON WH 0(3,5)90CM 2XHRC22BCVPF6DDP
C0027537	PREMICRON WH 0(3,5)90CM 2XHRC22 CVPF6DDP
C0027591	PREMICRON WHITE 2/0(3)90CM 2XHR20 CVDDP
C0027656	PREMICRON WHITE 2/0(3)75CM 2XHRC17CVDDP
C0027816	PREMICRON WHITE 2/0(3)90CM 2XHR26 CVDDP
C0027817	PREMICRON WHITE 0(3,5)90CM 2XHR26 CVDDP
C0027822	PREMICRON WHITE 5/0(1)45CM 2XDR8 CVDDP
C0027839	PREMICRON WHITE 2/0(3)75CM 2XHR26 CVDDP
C0027853	PREMICRON WHITE 5/0 (1)75CM 2XDR12CVDDP
C0027854	PREMICRON WHITE 4/0(1,5)75CM 2XDR12CVDDP
C0027906	PREMICRON WHITE 2/0(3)90CM 2XHR17CVDDP
C0027957	PREMICRON WHITE 2/0 (3)75CM DSMP24(M)DDP
C0027965	PREMICRON WHITE 2/0 (3) 75CM DS24 (M)DDP
C0027974	PREMICRON WHITE 3/0 (2) 45CM DS19 (M)DDP
C0100008	PREMICRON WHITE 1 (4) 2X75CM DDP
C0100017	PREMICRON WHITE 2/0 (3) 2X75CM DDP
C0102020	PREMICRON WHITE 1 (4) 150CM DDP
C0102025	PREMICRON WHITE 2 (5) 150CM DDP
C0120005	PREMICRON GREEN 5/0 (1) 10X45CM DDP

C0120006	PREMICRON GREEN 4/0 (1,5) 10X45CM DDP
C0120007	PREMICRON GREEN 3/0 (2) 10X45CM DDP
C0120008	PREMICRON GREEN 2/0 (3) 10X45CM DDP
C0120015	PREMICRON GREEN 4/0 (1,5) 10X60CM DDP
C0120016	PREMICRON GREEN 3/0 (2) 10X60CM DDP
C0120017	PREMICRON GREEN 2/0 (3) 10X60CM DDP
C0120020	PREMICRON GREEN 3/0 (2) 10X75CM DDP
C0120021	PREMICRON GREEN 2/0 (3) 10X75CM DDP
C0120025	PREMICRON GREEN 4/0 (1,5) 150CM DDP
C0120026	PREMICRON GREEN 3/0 (2) 150CM DDP
C0120027	PREMICRON GREEN 2/0 (3) 150CM DDP
C0120028	PREMICRON GREEN 0 (3,5) 150CM DDP
C0120029	PREMICRON GREEN 1 (4) 150CM DDP
C0120030	PREMICRON GREEN 2 (5) 150CM DDP
C0120036	PREMICRON GREEN 3/0 (2) 13X60CM DDP
C0120037	PREMICRON GREEN 2/0 (3) 13X60CM DDP
C0120050	PREMICRON GREEN 3/0 (2) 17X45CM DDP
C0120051	PREMICRON GREEN 2/0 (3) 17X45CM DDP
C0120080	PREMICRON GREEN 3/0 (2) 2X75CM DDP
C0120081	PREMICRON GREEN 2/0 (3) 2X75CM DDP
C0120082	PREMICRON GREEN 0 (3,5) 2X75CM DDP
C0120083	PREMICRON GREEN 1 (4) 2X75CM DDP
C0120084	PREMICRON GREEN 2 (5) 2X75CM DDP
C0120090	PREMICRON GREEN 5/0 (1) 3X45CM DDP
C0120091	PREMICRON GREEN 4/0 (1,5) 3X45CM DDP
C0120092	PREMICRON GREEN 3/0 (2) 3X45CM DDP
C0120093	PREMICRON GREEN 2/0 (3) 3X45CM DDP
C0120094	PREMICRON GREEN 0 (3,5) 3X45CM DDP
C0120095	PREMICRON GREEN 1 (4) 3X45CM DDP
C0120096	PREMICRON GREEN 2 (5) 3X45CM DDP
C0120097	PREMICRON GREEN 3 (6) 3X45CM DDP
C0120100	PREMICRON GREEN 5/0 (1) 5X45CM DDP
C0120101	PREMICRON GREEN 4/0 (1,5) 5X45CM DDP
C0120102	PREMICRON GREEN 3/0 (2) 5X45CM DDP
C0120103	PREMICRON GREEN 2/0 (3) 5X45CM DDP
C0120104	PREMICRON GREEN 0 (3,5) 5X45CM DDP
C0120110	PREMICRON GREEN 0 (3,5) 5X60CM DDP
C0120119	PREMICRON GREEN 3/0 (2) 5X75CM DDP
C0120120	PREMICRON GREEN 2/0 (3) 5X75CM DDP
C0120121	PREMICRON GREEN 0 (3,5) 5X75CM DDP
C0120139	PREMICRON GREEN 2/0 (3) 12X45CM DDP
C0120149	PREMICRON GREEN 3/0 (2) 6X45CM DDP
C2026091	PREMICRON GREEN 2/0(3)90CM 2XHR20 CV.RCP
C2026094	PREMICRON GREEN 4/0(1,5)75CM2XHR20CV.RCP
C2026116	PREMICRON GREEN 3/0(2)75CM 2XDR22 CV.RCP
C2026143	PREMICRON GREEN 2/0 (3)75CM 2XDR22CV.RCP
C2026157	PREMICRON GREEN 3/0 (2)60CM 2XHR17CV.RCP
C2026161	PREMICRON GRN2/0(3)75CM2XDRC18BCVPF6.RCP
C2026183	PREMICRON GREEN 3/0 (2) 60CM HR13 CV.RCP
C2026184	PREMICRON GRN2/0(3)75CM2XDRC18BCVPS6.RCP
C2026210	PREMICRON GRN2/0(3)75CM 2XHR17BCVPF6.RCP
C2026211	PREMICRON GRN 2/0(3)75CM 2XHR17B CV.RCP
C2026212	PREMICRON GREEN 0(3,5)75CM 2XHRC17CV.RCP
C2026213	PREMICRON GREEN 5/0(1)75CM 2XHR10FCV.RCP
C2026215	PREMICRON GREEN 3/0(2)75CM 2XHRC26CV.RCP
C2026240	PREMICRON GRN3/0(2)75CM 2XHRC17CVPF6.RCP
C2026241	PREMICRON GRN3/0(2)75CM 2XHRC17CVPF3.RCP
C2026248	PREMICRON GRN 2/0(3)75CM2XHRC17CVPF6.RCP
C2026249	PREMICRON GRN2/0(3)75CM 2XHRC17CVPF3.RCP
C2026250	PREMICRON GREEN 2/0(3)75CM 2XHRC17CV.RCP
C2026252	PREMICRON GRN 2/0(3)75CM2XHR17BCVPF3.RCP
C2026257	PREMICRON GRN 2/0(3)75CM2XHR17BCVPS3.RCP
C2026263	PREMICRON GRN 2/0(3)75CM2XHRC22CVPF6.RCP
C2026265	PREMICRON GRN 2/0(3)75CM 2XHR22B CV.RCP
C2026304	PREMICRON GR 4/0(1,5)90CM 2XHRC17 CV.RCP
C2026305	PREMICRON GRN 3/0(2)90CM 2XHRC17 CV.RCP
C2026306	PREMICRON GRN 2/0(3)90CM 2XHRC17 CV.RCP
C2026307	PREMICRON GRN 3/0(2)75CM2XHR22 CVPF6.RCP
C2026308	PREMICRON GRN2/0(3)90CM 2XHRC17CVPF6.RCP
C2026311	PREMICRON GRN 4/0(1,5)90CM2XHRC26CV.RCP
C2026313	PREMICRON GRN 2/0(3)90CM 2XHRC26 CV.RCP
C2026316	PREMICRON GRN 2/0(3)75CM 2XHRC26 CV.RCP
C2026317	PREMICRON GREEN 0(3,5)75CM 2XHRC26CV.RCP
C2026318	PREMICRON GRN 2/0(3)90CM2XHRC26CVPF6.RCP
C2026319	PREMICRON GRN 2/0(3)75CM2XHRC26CVPF6.RCP
C2026323	PREMICRON GRN 2/0(3)90CM 2XHRC22 CV.RCP
C2026419	PREMICRON GREEN 0(3,5)90CM 2XHRC26CV.RCP
C2026422	PREMICRON GRN2/0(3)75CM 2XHRC26CVPF3.RCP

C2026460	PREMICRON GRN 2/0(3)75CM2XHR22BCVPF6.RCP
C2026465	PREMICRON GREEN 3/0(2)75CM 2XDRC18CV.RCP
C2026466	PREMICRON GREEN 2/0(3)75CM 2XDRC18CV.RCP
C2026476	PREMICRON GRN2/0(3)75CM 2XDRC18CVPF3.RCP
C2026493	PREMICRON GRN 2/0(3)75CM2XHR22BCVPS6.RCP
C2026500	PREMICRON GRN 2/0(3)75CM 2XHR26B CV .RCP
C2026522	PREMICRON GRN 2/0(3)75CM2XHR26BCVPF6.RCP
C2026527	PREMICRON GRN 2/0(3)75CM2XHR26BCVPS6.RCP
C2026538	PREMICRON GRN 2/0(3)75CM 2XHRC17BCV.RCP
C2026556	PREMICRON GRN2/0(3)75CM2XHRC17BCVPF3.RCP
C2026560	PREMICRON GRN 2/0(3)75CM2XHRC17CVPS3.RCP
C2026563	PREMICRON GRN 2/0(3)75CM2XHRC17CVPS6.RCP
C2026568	PREMICRON GRN2/0(3)75CM2XHRC17BCVPS3.RCP
C2026572	PREMICRON GREEN 3/0(2)75CM 2XHRC17CV.RCP
C2026579	PREMICRON GRN2/0(3)75CM2XHRC17BCVPF6.RCP
C2026585	PREMICRON GRN 2/0(3)75CM2XHRC22CVPS6.RCP
C2026587	PREMICRON GRN2/0(3)75CM2XHRC17BCVPS6.RCP
C2026592	PREMICRON GRN 2/0(3)75CM 2XHRC18BCV.RCP
C2026594	PREMICRON GRN 2/0(3)75CM2XHRC26CVPS6.RCP
C2026596	PREMICRON GRN2/0(3)75CM2XHRC18BCVPF3.RCP
C2026600	PREMICRON GRN2/0(3)75CM2XHRC18BCVPS3.RCP
C2026606	PREMICRON GREEN 2/0 (3)45CM 2XHR17CV.RCP
C2026608	PREMICRON GRN2/0(3)75CM2XHRC18BCVPF6.RCP
C2026612	PREMICRON GRN2/0(3)75CM2XHRC18BCVPS6.RCP
C2026616	PREMICRON GRN 2/0(3)75CM 2XHRC22BCV.RCP
C2026620	PREMICRON GRN2/0(3)75CM2XHRC22BCVPF6.RCP
C2026625	PREMICRON GRN2/0(3)75CM2XHRC22BCVPS6.RCP
C2026634	PREMICRON GRN 2/0(3)75CM 2XHRC26BCV.RCP
C2026636	PREMICRON GRN 2/0(3)75CM2XDRC18CVPS3.RCP
C2026643	PREMICRON GRN2/0(3)75CM2XHRC26BCVPF6.RCP
C2026650	PREMICRON GRN2/0(3)75CM2XHRC26BCVPS6.RCP
C2026658	PREMICRON GRN 2/0(3)75CM 2XHRC27BCV.RCP
C2026665	PREMICRON GRN2/0(3)75CM2XHRC27BCVPF6.RCP
C2026674	PREMICRON GRN2/0(3)75CM2XHRC27BCVPS6.RCP
C2026676	PREMICRON GRN 2/0(3)90CM2XHRC26CVPS6.RCP
C2026695	PREMICRON GREEN 5/0 (1)75CM 2XHR10CV.RCP
C2026730	PREMICRON GRN 4/0(1,5)90CM 2XHR22CV.RCP
C2026731	PREMICRON GREEN 3/0(2)90CM 2XHR22CV.RCP
C2026732	PREMICRON GREEN 2/0(3)90CM 2XHR22CV.RCP
C2026733	PREMICRON GRN 2/0(3) 75CM2XHR22CVPF6.RCP
C2026734	PREMICRON GREEN 0(3,5)90CM 2XHR22CV.RCP
C2026735	PREMICRON GRN4/0(1,5)75CM2XHR22CVPF6.RCP
C2026737	PREMICRON GREEN 2/0(3)75CM 2XHR22CV.RCP
C2026738	PREMICRON GREEN 3/0(2)75CM 2XHR26CV.RCP
C2026739	PREMICRON GREEN 2/0(3)75CM 2XHR26CV.RCP
C2026763	PREMICRON GRN 4/0(1,5)75CM 2XHR22CV.RCP
C2026792	PREMICRON GRN 5/0 (1)90CM2XHR10CVPF3.RCP
C2026802	PREMICRON GREEN 6/0(0,7)75CM2XHR17CV.RCP
C2026803	PREMICRON GREEN 5/0(1)75CM 2XHR17 CV.RCP
C2026804	PREMICRON GREEN 4/0(1,5)75CM2XHR17CV.RCP
C2026805	PREMICRON GREEN 3/0(2)75CM 2XHR17 CV.RCP
C2026806	PREMICRON GREEN 2/0(3)75CM 2XHR17 CV.RCP
C2026814	PREMICRON GREEN 4/0(1,5)90CM2XHR26CV.RCP
C2026815	PREMICRON GREEN 3/0(2)90CM 2XHR26 CV.RCP
C2026816	PREMICRON GREEN 2/0(3)90CM 2XHR26 CV.RCP
C2026817	PREMICRON GREEN 0(3,5)90CM 2XHR26CV.RCP
C2026822	PREMICRON GREEN 5/0(1)45CM 2XDR8CV.RCP
C2026843	PREMICRON GREEN 5/0(1)75CM 2XHR13 CV.RCP
C2026844	PREMICRON GRN 4/0(1,5)75CM 2XHR13 CV.RCP
C2026845	PREMICRON GREEN 3/0(2)75CM 2XHR13 CV.RCP
C2026846	PREMICRON GREEN 3/0(2)60CM 2XHR13 CV.RCP
C2026852	PREMICRON GREEN 6/0(0,7)75CM2XDR12CV.RCP
C2026853	PREMICRON GREEN 5/0(1)75CM 2XDR12 CV.RCP
C2026854	PREMICRON GREEN 4/0(1,5)75CM2XDR12CV.RCP
C2026858	PREMICRON GREEN 3/0(2)90CM 2XDR18 CV.RCP
C2026863	PREMICRON GREEN 5/0(1)75CM 2XDR18 CV.RCP
C2026864	PREMICRON GREEN 4/0(1,5)75CM2XDR18CV.RCP
C2026865	PREMICRON GREEN 3/0(2)75CM 2XDR18 CV.RCP
C2026871	PREMICRON GRN 6/0(0,7)75CM 2XDR10 CV.RCP
C2026873	PREMICRON GREEN 2/0 (3)90CM 2XDR26CV.RCP
C2026874	PREMICRON GRN 2/0(3)90CM2XDR26 CVPF6.RCP
C2026883	PREMICRON GREEN 2/0(3)90CM 2XDR18 CV.RCP
C2026887	PREMICRON GRN 5/0(1)75CM2XHR10 CVPF3.RCP
C2026888	PREMICRON GRN4/0(1,5)75CM2XHR10CVPF3.RCP
C2026896	PREMICRON GRN4/0(1,5)75CM2XDR13CVPF6.RCP
C2026903	PREMICRON GREEN 5/0(1)90CM 2XHR17 CV.RCP
C2026904	PREMICRON GREEN 4/0(1,5)90CM2XHR17CV.RCP
C2026905	PREMICRON GREEN 3/0(2)90CM 2XHR17 CV.RCP

C2026906	PREMICRON GREEN 2/0(3)90CM 2XHR17 CV.RCP
C2026915	PREMICRON GRN 3/0(2) 90CM2XHR26CVPF6.RCP
C2026916	PREMICRON GRN 2/0 (3)90CM2XHR26CVPF6.RCP
C2026932	PREMICRON GRN 3/0 (2)120CM 2XHR26 CV.RCP
C2026940	PREMICRON GREEN 0(3,5)75CM HRC17 CV.RCP
C2026942	PREMICRON GREEN 6/0(0,7)75CM HR6(M).RCP
C2026949	PREMICRON GREEN 5/0(1)75CM2XHR13FFCV.RCP
C2026954	PREMICRON GRN4/0(1,5)75CM2XDR12CVPF3.RCP
C2026955	PREMICRON GREEN 2/0(3)120CM HRC26BCV.RCP
C2026956	PREMICRON GREEN 0(3,5)120CM HRC26BCV.RCP
C2026957	PREMICRON GREEN2/0(3)120CM2XHR26BCV.RCP
C2026960	PREMICRON GREEN0(3,5)120CM2XHR26BCV.RCP
C2026995	PREMICRON GREEN 0(3,5)90CM2XHR22BCV.RCP
C2027014	PREMICRON WHITE 6/0(0,7)75CM2XDR12CV.RCP
C2027034	PREMICRON WHITE2/0(3)75CM2XHR22CVPF6.RCP
C2027042	PREMICRON GRN 2/0(3)75CM 2XHR17CVPF6.RCP
C2027061	PREMICRON WHITE 6/0(0,7)45CM 2XHLM8M.RCP
C2027069	PREMICRON WHITE 6/0(0,7)75CM 2XDR8CV.RCP
C2027070	PREMICRON WHITE 5/0(1)75CM 2XDR8 CV.RCP
C2027071	PREMICRON WHITE 4/0(1,5)75CM 2XDR8CV.RCP
C2027306	PREMICRON WHITE2/0(3)90CM2XHR26CVPF6.RCP
C2027315	PREMICRON WHIT2/0(3)90CM2XHR26CVPF6.RCP
C2027317	PREMICRON WHITE 2/0(3)90CM 2XHR26CV.RCP
C2027320	PREMICRON WHIT2/0(3)75CM2XHR26CVPS3.RCP
C2027323	PREMICRON WHIT2/0(3)75CM2XHR26CVPS6.RCP
C2027326	PREMICRON W 2/0(3)90CM 2XHR26 CV.RCP
C2027340	PREMICRON WHIT2/0(3)75CM2XHR26CVPS6.RCP
C2027343	PREMICRON WHIT2/0(3)75CM2XHR26CVPS6.RCP
C2027363	PREMICRON WHITE 2/0(3)90CM 2XHR26CV.RCP
C2027422	PREMICRON WHIT2/0(3)75CM2XHR26CVPF3.RCP
C2027455	PREMICRON WHITE 4/0(1,5)75CM2XHR17CV.RCP
C2027457	PREMICRON WHITE 2/0 (3)75CM 2XHR17CV.RCP
C2027537	PREMICRON WHIT0(3,5)90CM2XHR22CVPF6.RCP
C2027591	PREMICRON WHITE 2/0(3) 90CM 2XHR20CV.RCP
C2027656	PREMICRON WHITE 2/0(3)75CM 2XHR26CV.RCP
C2027815	PREMICRON WHITE 3/0(2)90CM 2XHR26 CV.RCP
C2027816	PREMICRON WHITE 2/0(3)90CM 2XHR26 CV.RCP
C2027817	PREMICRON WHITE 0(3,5)90CM 2XHR26CV.RCP
C2027822	PREMICRON WHITE 5/0(1)45CM 2XDR8 CV.RCP
C2027839	PREMICRON WHITE 2/0(3)75CM 2XHR26 CV.RCP
C2027853	PREMICRON WHITE 5/0 (1)75CM 2XDR12CV.RCP
C2027854	PREMICRON WHITE 4/0(1,5)75CM2XDR12CV.RCP
C2027906	PREMICRON WHITE 2/0(3)90CM 2XHR17 CV.RCP
G0026081	PREMICRON GREEN 5 (7) 4X75CM HRT48(M)DDP
G0026114	PREMICRON GREEN 5 (7) 4X75CM HRC55(M)DDP
G0026146	PREMICRON GREEN 5 (7) 75CM HRT55 (M)DDP
G0026148	PREMICRON GREEN 5 (7) 75CM HRT48 (M)DDP
G0026155	PREMICRON GREEN 5/0 (1)45CM 2XVLM8(M)DDP
G0026203	PREMICRON GREEN 5 (7) 4X75CM HS48(M) DDP
G0026702	PREMICRON GREEN 5 (7) 4X75CM HRS48 DDP
G0026705	PREMICRON GREEN 5 (7) 4X75CM HRS60 DDP
G0026723	PREMICRON GREEN 5 (7) 4X75CM HRT55(M)DDP
G0027062	PREMICRON WHITE 5/0 (1)45CM 2XHLM8 (MDDP
G0027114	PREMICRON WHITE 6/0(0,7)45CM 2XVLM8(MDDP
G0027115	PREMICRON WHITE 5/0 (1)45CM 2XVLM8 (MDDP
G0027116	PREMICRON WHITE 4/0(1,5)45CM2XVLM8(M)DDP
G0027895	PREMICRON WHITE 5/0(1)45CM 2XVLM6(M)DDP
G0120055	PREMICRON GREEN 5/0 (1) 250CM ARO DDP
G0120056	PREMICRON GREEN 4/0 (1,5) 250CM ARO DDP
G0120060	PREMICRON GREEN 3/0 (2) 250CM ARO DDP
G0120061	PREMICRON GREEN 2/0 (3) 250CM ARO DDP
G0120062	PREMICRON GREEN 0 (3,5) 250CM ARO DDP
G0120063	PREMICRON GREEN 1 (4) 250CM ARO DDP
G0120064	PREMICRON GREEN 2 (5) 250CM ARO DDP
G0120065	PREMICRON GREEN 3 (6) 250CM ARO DDP
G0120067	PREMICRON GREEN 5 (7) 250CM ARO DDP
M0026175	PREMICRON GRN 2/0 (3)4X75 2XHR17 CV MLP
M0026225	PREMICRON GRN 2/0 (3) 4X75 2XHR26 CV MLP
M0026699	PREMICRON GRN 4/0(1,5)4X75 HR13 CV MLP
M0026708	PREMICRON GRN 4/0(1,5)4X75 2XHR17 CV MLP
M0026712	PREMICRON GRN 2/0 8X90 2XHR26 CV PF6 MLP
M0026725	PREMICRON GRN 2/0 (3) 4X75 HR26 CV MLP
M0027001	PREMICRON G/W 2/0(3)8X90CM 2XHR22(M)MLP
M0027002	PREMICRON G/W2/0(3)4X75CM2XDR18BCVPF6MLP
M0027003	PREMICRON GW2/0(3)4X90CM2XHR26BCVPF6MLP
M0027004	PREMICRON G/W 2/0(3)4X90CM2XHR22MPF6MLP
M0027008	PREMICRON G/W 0(3,5)4X90CM 2XHR17 CV MLP
M0027009	PREMICRON G/W 0(3,5)4X75CM 2XHR26 CV MLP

M0027013	PREMICRONG/W2/0(3)4X75CM2XDRC18BCVPS6MLP
M0027019	PREMICRON G/W 0(3,5)4X75CM 2XHR22 CV MLP
M0027046	PREMICRONG/W2/0(3)8X75CM2XDRC18BCVPF6MLP
M0027079	PREMICRONG/W2/0(3)8X75CM2XDRC18BCVPS6MLP
M0027088	PREMICRON G/W 2/0(3)4X75CM2XHR17CVPF6MLP
M0027089	PREMICRON G/W2/0(3)4X75CM2XHR17BCVPF6MLP
M0027090	PREMICRON G/W2/0(3)4X75CM2XHR17BCVPF3MLP
M0027091	PREMICRON G/W 0(3,5)4X75CM2XHR17CVPF6MLP
M0027100	PREMICRN G/W 2/0 8X75 2XHRC22 CV PF6 MLP
M0027102	PREMICRON G/W 2/0 (3)8X75 2XHRC22 CV MLP
M0027105	PREMICRON G/W2/0(3)4X75CM2XHR17BCVPS3MLP
M0027107	PREMICRON G/W 0(3,5)4X75 2XHRC22 CV MLP
M0027108	PREMICRON GW 0(3,5)4X75CM2XHRC22CVPF6MLP
M0027109	PREMICRONG/W0(3,5)4X75CM2XHRC22BCVPF6MLP
M0027112	PREMICRON WHT 2/0(3) 4X75CM HR26 CV MLP
M0027119	PREMICRON G/W2/0(3)8X75CM2XHR17BCVPF3MLP
M0027133	PREMICRON G/W2/0(3)8X75CM2XHR17BCVPS3MLP
M0027150	PREMICRON G/W2/0(3)8X75CM2XHR17BCVPF6MLP
M0027151	PREMICRON G/W 2/0(3)8X75CM 2XHR17B CVMLP
M0027157	PREMICRON G/W2/0(3)4X75CM2XHR22BCVPF6MLP
M0027158	PREMICRON G/W 0(3,5)4X75CM2XHR22CVPF6MLP
M0027163	PREMICRON G/W2/0(3)4X75CM2XHR22BCVPS6MLP
M0027169	PREMICRON G/W2/0(3)8X75CM2XHR22BCVPF6MLP
M0027175	PREMICRON G/W2/0(3)8X75CM2XHR22BCVPS6MLP
M0027181	PREMICRON G/W 2/0(3)8X75CM 2XHR22B CVMLP
M0027186	PREMICRON G/W 2/0(3)4X75CM 2XHR26B CVMLP
M0027187	PREMICRON G/W2/0(3)4X75CM2XHR26BCVPF6MLP
M0027193	PREMICRON G/W2/0(3)4X75CM2XHR26BCVPS6MLP
M0027199	PREMICRON G/W2/0(3)8X75CM2XHR26BCVPF6MLP
M0027205	PREMICRON G/W 2/0(3)8X75CM 2XHR26B CVMLP
M0027211	PREMICRON G/W2/0(3)8X75CM2XHR26BCVPS6MLP
M0027217	PREMICRONG/W2/0(3)4X75CM2XHRC17BCVPF3MLP
M0027223	PREMICRONG/W2/0(3)4X75CM2XHRC17BCVPS3MLP
M0027229	PREMICRONG/W2/0(3)4X75CM2XHRC17BCVPF6MLP
M0027235	PREMICRONG/W2/0(3)4X75CM2XHRC17BCVPS6MLP
M0027241	PREMICRONG/W2/0(3)8X75CM2XHRC17BCVPF3MLP
M0027247	PREMICRONG/W2/0(3)8X75CM2XHRC17BCVPS3MLP
M0027251	PREMICRON G/W 2/0(3)4X75CM2xHRC22BCVMLP
M0027253	PREMICRON G/W 2/0(3)8X75CM 2XHRC17BCVMLP
M0027259	PREMICRONG/W2/0(3)8X75CM2XHRC17BCVPF6MLP
M0027265	PREMICRONG/W2/0(3)8X75CM2XHRC17BCVPS6MLP
M0027271	PREMICRONG/W2/0(3)4X75CM2XHRC18BCVPF3MLP
M0027277	PREMICRONG/W2/0(3)4X75CM2XHRC18BCVPS3MLP
M0027283	PREMICRONG/W2/0(3)4X75CM2XHRC18BCVPF6MLP
M0027289	PREMICRONG/W2/0(3)4X75CM2XHRC18BCVPS6MLP
M0027295	PREMICRONG/W2/0(3)8X75CM2XHRC18BCVPF3MLP
M0027300	PREMICRON GW2/0(3)4X90CM2XHRC22BCVPF6MLP
M0027302	PREMICRONG/W2/0(3)8X75CM2XHRC18BCVPS3MLP
M0027310	PREMICRN G/W 2/0 8X75 2XDRC18 CV PF3 MLP
M0027311	PREMICRON G/W 2/0 (3)8X75 2XDRC18 CV MLP
M0027312	PREMICRN G/W 2/0 8X75 2XDRC18 CV PF6 MLP
M0027353	PREMICRON G/W 2/0(3)8X75CM 2XHRC18BCVMLP
M0027360	PREMICRONG/W2/0(3)8X75CM2XHRC18BCVPF6MLP
M0027368	PREMICRONG/W2/0(3)8X75CM2XHRC18BCVPS6MLP
M0027371	PREMICRON G/W 2/0(3) 4X75 2XHRC26 CV MLP
M0027375	PREMICRONG/W2/0(3)4X75CM2XHRC22BCVPF6MLP
M0027383	PREMICRONG/W2/0(3)4X75CM2XHRC22BCVPS6MLP
M0027385	PREMICRON G/W 0(3,5)4X75CM 2XHRC22BCVMLP
M0027391	PREMICRONG/W2/0(3)8X75CM2XHRC22BCVPF6MLP
M0027394	PREMICRON G/W2/0(3)4X75CM2XDRC18CVPS3MLP
M0027402	PREMICRONG/W2/0(3)8X75CM2XHRC22BCVPS6MLP
M0027414	PREMICRON G/W 2/0(3)8X75CM 2XHRC22BCVMLP
M0027424	PREMICRONG/W2/0(3)4X75CM2XHRC26BCVPF6MLP
M0027425	PREMICRON G/W 2/0(3)4X75CM 2XHRC26BCVMLP
M0027443	PREMICRONG/W2/0(3)4X75CM2XHRC26BCVPS6MLP
M0027451	PREMICRONG/W2/0(3)8X75CM2XHRC26BCVPF6MLP
M0027465	PREMICRON G/W 2/0(3)8X75CM 2XHRC26BCVMLP
M0027471	PREMICRONG/W2/0(3)8X75CM2XHRC26BCVPS6MLP
M0027483	PREMICRONG/W2/0(3)4X75CM2XHRC27BCVPS6MLP
M0027495	PREMICRONG/W2/0(3)8X75CM2XHRC27BCVPS6MLP
M0027501	PREMICRON G/W 2/0(3)8X75CM 2XHRC27BCVMLP
M0027507	PREMICRON G/W2/0(3)8X75CM2XDRC18CVPS3MLP
M0027508	PREMICRON G/W 0(3,5) 8X75CM 2XDRC18CVMLP
M0027513	PREMICRON G/W2/0(3)4X90CM2XDRC18CVPS6MLP
M0027519	PREMICRON G/W2/0(3)8X90CM2XDRC18CVPS6MLP
M0027525	PREMICRON G/W2/0(3)8X75CM 2XHR20CVPS6MLP
M0027531	PREMICRON G/W2/0(3)8X90CM2XHRC26CVPS6MLP
M0027543	PREMICRON G/W2/0(3)4X75CM2XHRC26CVPF6MLP

M0027549	PREMICRON G/W 2/0(3)4X75CM 2XHRC22 CVMLP
M0027555	PREMICRON G/W 2/0(3)4X90CM 2XHRC26CV MLP
M0027565	PREMICRON G/W2/0(3)4X90 2XHRC26BCVPS6MLP
M0027566	PREMICRON G/W 2/0(3)4X90CM 2XHRC22BCVMLP
M0027567	PREMICRON G/W2/0(3)4X90 2XHRC22BCVPS6MLP
M0027568	PREMICRON G/W 2/0(3)4X75CM 2XHRC17BCVMLP
M0027569	PREMICRON G/W 2/0(3)4X90CM 2XHRC26BCVMLP
M0027579	PREMICRON G/W0(3,5)4X75CM2XHRC17CVPF6MLP
M0027580	PREMICRON GW0(3,5)4X75CM2XHRC17BCVPF6MLP
M0027584	PREMICRON G/W 3/0(2)4X75CM 2XHRC17 CVMLP
M0027585	PREMICRON G/W 0 (3,5)4X75CM2XHRC17CV MLP
M0027586	PREMICRON G/W 0(3,5)4X75CM 2XHRC17BCVMLP
M0027600	PREMICRN G/W 2/0 4X75 2XDRC18 CV PF3 MLP
M0027601	PREMICRON G/W 2/0(3) 4X75 2XDRC18 CV MLP
M0027603	PREMICRON G/W 2/0(3)4X90CM 2XHRC22CV MLP
M0027605	PREMICRON G/W 2/0 8X75 2XHR22 CV PF6 MLP
M0027610	PREMICRON G/W 0(3,5)4X90 2XHRC22 CV MLP
M0027611	PREMICRON G/W 0(3,5)4X90CM 2XHRC22BCVMLP
M0027615	PREMICRON G/W 2/0(3)4X90CM2XHRC22 (M)MLP
M0027617	PREMICRON G/W2/0(3)4X90CM2XHRC22CVPF6MLP
M0027618	PREMICRON G/W 0 4X90 2XHRC22 CV PF6 MLP
M0027619	PREMICRONG/W0(3,5)4X90CM2XHRC22BCVPF6MLP
M0027625	PREMICRON GW2/0(3)8X75CM2XHRC20CVPF6MLP
M0027626	PREMICRON GW2/0(3)4X75CM2XHRC20CVPF6MLP
M0027628	PREMICRN G/W 2/0 4X75 2XHRC17 CV PF6 MLP
M0027660	PREMICRN G/W 2/0 8X90 2XHRC26 CV PF6 MLP
M0027661	PREMICRON G/W 0 8X90 2XHRC26 CV PF6 MLP
M0027671	PREMICRON G/W 2/0 8X75 2XDR18 CV PF6 MLP
M0027675	PREMICRON G/W0(3,5)8X75CM2XHRC17CVPF6MLP
M0027685	PREMICRN G/W 2/0 4X90 2XHRC26 CV PF6 MLP
M0027686	PREMICRON G/W 0(3,5)4X90CM 2XHRC26 CVMLP
M0027687	PREMICRON G/W0(3,5)4X90CM2XHRC26CVPF6MLP
M0027688	PREMICRON G/W 0(3,5)4X90CM 2XHRC26BCVMLP
M0027689	PREMICRON GW0(3,5)4X90CM2XHRC26BCVPF6MLP
M0027700	PREMICRON G/W 2/0(3) 4X75 2XHR26 CV MLP
M0027704	PREMICRON G/W 2/0 8X90 2XHR22 CV PF6 MLP
M0027705	PREMICRON G/W 2/0 (3) 8X75 2XHR22 CV MLP
M0027706	PREMICRON G/W 2/0 8X75 2XHR20 CV PF6 MLP
M0027707	PREMICRON G/W 2/0 (3) 8X75 2XHR20 CV MLP
M0027712	PREMICRON G/W 2/0 8X90 2XHR26 CV PF6 MLP
M0027713	PREMICRON G/W 2/0 (3) 8X90 2XHR26 CV MLP
M0027715	PREMICRON G/W 2/0 (3) 8X75 2XHR26 CV MLP
M0027716	PREMICRON G/W 2/0 8X75 2XHR26 CV PF3 MLP
M0027717	PREMICRON G/W 2/0 8X90 2XHR26 CV PF3 MLP
M0027722	PREMICRON G/W 2/0(3) 8X75 2XHRC17 CV MLP
M0027723	PREMICRN G/W 2/0 8X75 2XHRC17 CV PF6 MLP
M0027738	PREMICRN G/W 2/0 8X75 2XHRC17 CV PF3 MLP
M0027739	PREMICRON GW0(3,5)4X75CM 2XHRC17CVPF3MLP
M0027752	PREMICRON G/W 3/0 4X90 2XHR26 CV PF6 MLP
M0027759	PREMICRON G/W 0 4X75 2XHR26 CV PF6 MLP
M0027764	PREMICRON G/W 2/0 4X75 2XHR17 CV PF3 MLP
M0027766	PREMICRON G/W 2/0 8X75 2XHR26 CV PF6 MLP
M0027769	PREMICRON G/W 2/0 (3) 8X90 2XDR18 CV MLP
M0027775	PREMICRON G/W 2/0 (3) 8X75 2XHR17 CV MLP
M0027776	PREMICRON G/W 2/0 8X75 2XHR17 CV PF6 MLP
M0027791	PREMICRON G/W 2/0(3) 8X75 2XHRC26 CV MLP
M0027792	PREMICRN G/W 2/0 8X75 2XHRC26 CV PF6 MLP
M0027801	PREMICRON G/W 2/0(3)8X75 2XHRC26CVPF3MLP
M0027802	PREMICRON G/W 2/0 8X75 2XHR17 CV PF3 MLP
M0027811	PREMICRN G/W 2/0 8X75 2XHRC17 CV PS6 MLP
M0027813	PREMICRN G/W 2/0 4X75 2XHRC17 CV PS6 MLP
M0027819	PREMICRN G/W 2/0 8X75 2XHRC17 CV PS3 MLP
M0027824	PREMICRN G/W 2/0 4X75 2XHRC17 CV PS3 MLP
M0027836	PREMICRN G/W 2/0 4X75 2XHRC22 CV PS6 MLP
M0027841	PREMICRN G/W 2/0 8X75 2XHRC22 CV PS6 MLP
M0027843	PREMICRN G/W 2/0 4X75 2XHRC26 CV PS6 MLP
M0027845	PREMICRN G/W 2/0 8X75 2XHRC26 CV PS6 MLP
M0027858	PREMICRON GW0(3,5)4X75CM2XHRC26BCVPF6MLP
M0027859	PREMICRON G/W0(3,5)4X75CM2XHRC26CVPF6MLP
M0027863	PREMICRN G/W 2/0 8X90 2XHRC17 CV PF6 MLP
M0027867	PREMICRON G/W2/0(3)8X75CM2XHR22 CVPS6MLP
M0027905	PREMICRON G/W 2/0 4X75 2XHR20 CV PF6 MLP
M0027915	PREMICRON G/W 2/0(3)8X75CM HRC26 CV MLP
M0027921	PREMICRON G/W 2/0(3) 4X75 2XHRC17 CV MLP
M0027925	PREMICRON G/W0 (3,5)4X75CM 2XHRC26 CVMLP
M0027926	PREMICRON G/W 0(3,5)4X75CM 2XHRC26BCVMLP
M0027929	PREMICRON G/W 2/0(3)8X90CM 2XHRC17CV MLP
M0027935	PREMICRN G/W 2/0 8X90 2XDRC18 CV PF6 MLP

M0027941	PREMICRON G/W 2/0(3) 8X90 2XDRC18 CV MLP
M0027945	PREMICRON G/W2/0(3)8X75CM 2XHR20CVPS3MLP
M0027949	PREMICRN G/W 2/0 4X75 2XHRC22 CV PF6 MLP
M0027953	PREMICRON G/W2/0(3)8X75CM 2XHR20CVPF3MLP
M0027961	PREMICRON G/W3/0(2)4X90CM2XHRC17MPS3MLPP
M0027969	PREMICRON G/W3/0(2)4X90CM2XHRC17MPS6MLP
M0027971	PREMICRON GW0(3,5)4X90CM2XHRC17BCVPF6MLP
M0027972	PREMICRON G/W 0(3,5)4X90CM 2XHRC17BCVMLP
M0027975	PREMICRON G/W0(3,5)4X90CM2XHRC17PF6CVMLP
M0027976	PREMICRON G/W 0(3,5)4X90CM 2XHRC17 CVMLP
M0027979	PREMICRON G/W 2/0(3)8X75CM 2XFRC22(M)MLP
X0026006	PREMICRON GREEN 2/0 (3) 75CM HR17 CVDDP
X0026008	PREMICRON GREEN 6/0(0,7)45CM HR10 CVDDP
X0026016	PREMICRON GREEN 2/0 (3) 75CM HR26 CVDDP
X0026036	PREMICRON GREEN 2/0 (3) 75CM HR37 (M)DDP
X0026121	PREMICRON GREEN 5/0 (1) 45CM DR12 CVDDP
X0026155	PREMICRON GREEN 5/0 (1)45CM 2XVLM8(M)DDP
X0026156	PREMICRON GREEN 4/0(1,5)75CM HRC22CV DDP
X0026173	PREMICRON GREEN 4/0 (1,5)45CM DS12(M)DDP
X0026190	PREMICRON GREEN 2/0 (3) 45CM DS24 (M)DDP
X0026199	PREMICRON GREEN 4/0 (1,5)45CM HS21(M)DDP
X0026208	PREMICRON GREEN 5(7)50CM DS76(M)PEEL CE
X0026223	PREMICRON GRN 5(7)75CM 2XDS76(M)PEEL CE
X0026278	PREMICRON GREEN 2 (5) 75CM DS60(M) DDP
X0026331	PREMICRON GREEN 1 (4) 75CM DS39(M) DDP
X0026335	PREMICRON GREEN 3/0 (2) 75CM HRT26(M)DDP
X0026336	PREMICRON GREEN 2/0 (3) 75CM HRT26(M)DDP
X0026337	PREMICRON GREEN 0 (3,5) 75CM HRT26(M)DDP
X0026352	PREMICRON GREEN 6/0(0,7)75CM DRC12CVDDP
X0026394	PREMICRON GREEN 2 (5) 75CM HRT37 (M) DDP
X0026404	PREMICRON GREEN 1 (4) 75CM HS23SS(M) DDP
X0026420	PREMICRON GREEN 3/0 (2) 75CM HS26(M) DDP
X0026472	PREMICRON GREEN 1 (4) 75CM HS37(M) DDP
X0026514	PREMICRON GREEN 4/0 (1,5)75CM DS24(M)DDP
X0026531	PREMICRON GREEN 2/0 (3) 75CM DS30(M) DDP
X0026542	PREMICRON GREEN 0 (3,5) 75CM HS30(M) DDP
X0026543	PREMICRON GREEN 1 (4) 75CM HS30(M) DDP
X0026854	PREMICRON GRN 4/0(1,5)75CM 2XDR12 CVDDP
X0026903	PREMICRON GREEN 5/0(1)90CM 2XHR17 CVDDP
X0027115	PREMICRON WHITE 5/0 (1)45CM 2XVLM8 (MDDP
X0102025	PREMICRON WHITE 2 (5) 150CM DDP
X0120006	PREMICRON GREEN 4/0 (1,5) 10X45CM DDP
X0120020	PREMICRON GREEN 3/0 (2) 10X75CM DDP
X0120021	PREMICRON GREEN 2/0 (3) 10X75CM DDP
X0120027	PREMICRON GREEN 2/0 (3) 150CM DDP
X0120063	PREMICRON GREEN 1 (4) 250CM ARO DDP
X0120080	PREMICRON GREEN 3/0 (2) 2X75CM DDP
X0120081	PREMICRON GREEN 2/0 (3) 2X75CM DDP
X0120082	PREMICRON GREEN 0 (3,5) 2X75CM DDP
X0120083	PREMICRON GREEN 1 (4) 2X75CM DDP
X0120084	PREMICRON GREEN 2 (5) 2X75CM DDP
X0120085	PREMICRON GREEN 3 (6) 2X75CM DDP
X0120090	PREMICRON GREEN 5/0 (1) 3X45CM DDP
X0120091	PREMICRON GREEN 4/0 (1,5) 3X45CM DDP
X0120102	PREMICRON GREEN 3/0 (2) 5X45CM DDP
B0026043	PREMICRON GREEN 2 (5) 2X75CM HRT48 (M)
B0026188	PREMICRON GREEN 5 (7) 75CM 2XDS80MPEELCE
B0026207	PREMICRON GREEN 1 (4) 90CM HRC30SS(M)DDP
B0026219	PREMICRON GREEN 2(5) 4X75CM HRC37S(M)DDP
B0026221	PREMICRON GREEN 2(5) 4X75CM HRC40S(M)DDP
B0026247	PREMICRON GREEN 2 (5) 4X75CM HRC43(M)DDP
B0026258	PREMICRON GREEN 2 (5) 75CM HRC48 (M) DDP
B0026808	PREMICRON GREEN 1(4)75 2XHR76 (M)PEEL CE
B0027014	PREMICRON G/W 2/0(3)2X752XHRC26 CVPF6DDP
B0027027	PREMICRON G/W 2/0(3)2X752XHRC17 CVPF6DDP
C0026097	PREMICRON GREEN 2 (5) 75CM HRT48 (M)
C0026182	PREMICRON GREEN 2 (5) 75CM HR48 (M)
C0026206	PREMICRON GREEN 0(3,5)75CM HRC30SS(M)DDP
C0026216	PREMICRON GREEN 2 (5) 75CM HRC37S (M)DDP
C0026218	PREMICRON GREEN 2/0(3)90CM2XHRC37S(M)DDP
C0026224	PREMICRON GREEN 0(3,5) 75CM HRC40S(M)DDP
C0026229	PREMICRON GREEN 2 (5) 75CM HRC40SS(M)DDP
C0026243	PREMICRON GREEN 4/0 (1,5)75CM DS12(M)DDP
C0026256	PREMICRON GREEN 2 (5) 75CM HRC43SS(M)DDP
C0026264	PREMICRON GREEN 5/0 (1) 75CM DRC12CV DDP
C0026273	PREMICRON GREEN 5/0 (1) 75CM DRC18CV DDP
C0026274	PREMICRON GREEN 3/0 (2) 75CM DRC18CV DDP
C0026276	PREMICRON GREEN 2/0(3)90CM 2XDRC18CV DDP

C0026278	PREMICRON GREEN 2 (5) 75CM DS60(M)
C0026311	PREMICRON GREEN 4/0(1.5)90CM 2XHRC26 (M)
C0026355	PREMICRON GREEN 2 (5) 75CM DS48(M) DDP
C0026364	PREMICRON GREEN 4/0 (1.5) 75CM DRC18 (M)
C0026376	PREMICRON GREEN 2 (5) 75CM HS48(M) DDP
C0026400	PREMICRON GREEN 2/0 (3) 75CM DRC26CV DDP
C0026410	PREMICRON GRN 2/0(3)90CM2XDRC26 CVPF6DDP
C0026424	PREMICRON GREEN 2/0 (3) 75CM HRC17CV DDP
C0026456	PREMICRON GREEN 5/0(1)90CM 2XHRC17CV DDP
C0026459	PREMICRON GREEN 0 (3,5) 75CM HRC26CV DDP
C0026474	PREMICRON GREEN 4/0 (1.5)75CM 2XDRC18(M)
C0026475	PREMICRON GREEN 1 (4) 90CM HRC26 CV DDP
C0026478	PREMICRON GREEN 0(3,5) 75CM HRC37S(M)DDP
C0026492	PREMICRON GREEN 2/0 (3) 45CM DS30 (M)DDP
C0026497	PREMICRON GREEN 1 (4) 75CM HRC37S (M)DDP
C0026509	PREMICRON GREEN 2 (5) 75CM GS51(M)
C0026720	PREMICRON GREEN 2 (5) 75CM 2XHR48 (M)
C0026759	PREMICRON GREEN 2 (5) 90CM 2XHR37S (M)
C0026808	PREMICRON GREEN 1(4) 75CM 2XHR76 (M)PEEL
C0026998	PREMICRON GRN 4/0(1,5)45CM 2XGRC24(M)DDP
C0027021	PREMICRON WHITE 3/0 (2) 75CM DRC18CV DDP
C0027022	PREMICRON WHITE 4/0(1,5)75CM DRC18CV DDP
C0027023	PREMICRON WHITE 3/0 (2) 75CM DRC26CV DDP
C0027038	PREMICRON WHITE 3/0(2)90CM 2XHRC17CV DDP
C0027041	PREMICRON WHITE 2/0(3)75CM 2XHRC26CV DDP
C0027323	PREMICRON WHT 2/0(3)75CM2XHRC17 CV PS6
C0102035	PREMICRON WHITE 3 (6) 150CM
C0120009	PREMICRON GREEN 0 (3,5) 10X45CM
C0120010	PREMICRON GREEN 1 (4) 10X45CM
C0120031	PREMICRON GREEN 3 (6) 150CM
C0120032	PREMICRON GREEN 4 (6) 150CM
C0120085	PREMICRON GREEN 3 (6) 2X75CM DDP
C0120122	PREMICRON GREEN 1 (4) 5X75CM
C2026312	PREMICRON GRN 3/0(2)90CM 2XHRC26 CV
C2026456	PREMICRON GREEN 5/0(1)90CM2XHRC17CV .RCP
C2026474	PREMICRON GRN 4/0(1,5)75CM 2XDRC18 CV
C2026494	PREMICRON GRN 0(3,5)752XHRC17CV PF6 .RCP
C2027038	PREMICRON WHITE 3/0(2)90CM2XHRC17CV .RCP
C2027039	PREMICRON WHT 2/0(3)752XHRC17CV PF6 .RCP
C2027041	PREMICRON WHITE 2/0(3)75CM2XHRC26CV .RCP
C2027316	PREMICRON WHITE2/0(3)90CM2XHRC17CVPF3RCP
C2027336	PREMICRON WHITE2/0(3)90CM2XHRC26CVPF6RCP
G0026166	PREMICRON GREEN 5(7)4X75CM HRC48 (M)
G0026222	PREMICRON GREEN 5(7) 4X75CM HRC40S(M)DDP
G0026259	PREMICRON GREEN 5 (7) 75CM HRC48 (M) DDP
G0026262	PREMICRON GREEN 5 (7) 75CM HRC55 (M) DDP
M0027026	PREMICRON G/W 2/0(3)8X902XDRC26 CVPF6MLP
M0027029	PREMICRON G/W 2/0(3)4X902XHRC17 CVPF6MLP
M0027030	PREMICRON G/W 2/0(3)8X902XHRC17CV PF3MLP
M0027032	PREMICRON G/W 4/0 8X75 2XHRC17 CV PF3MLP
M0027033	PREMICRON G/W 4/0 8X75 2XHRC17 CV PF6MLP
M0027037	PREMICRON G/W 4/0 8X75 2XHRC22 CV PF6MLP
M0027125	PREMICRON G/W 2/0(3)8X75CM 2XHR22MPF6MLP
M0027127	PREMICRON G/W 2/0(3)8X75CM 2XHR22(M) MLP
M0027128	PREMICRON G/W 2/0(3)8X75CM2XHRC17MPF6MLP
M0027372	PREMICRON G/W 2/0(3)4X75 2XHRC26CVPF3MLP
M0027616	PREMICRON G/W 2/0(3)4X90CM 2XHRC22(M)MLP
M0027647	PREMICRON G/W 2/0 4X90 2XDRC18 CV PF6MLP
M0027910	PREMICRON G/W 2/0(3)8X90CM 2XHRC26(M)MLP
X0120095	PREMICRON GREEN 1 (4) 3X45CM