

RT330 Optiflow Tubing Kit

RT330 CIRCUIT, CHAMBER AND PRESSURE MANIFOLD

PERFORMANCE SPECS

Max Flow Rate (L/Min)	0.5-25
Length Of Tubing	1.75m/5.6ft
Ambient Range °C	18-26 °C
Connection	Unique Swivel connection to Optiflow Cannula (OPT31X)
Compatible Humidifier	MR850
Humidifier Mode	Invasive
Duration Of Use	7 days
Carton Size	Box of 10
Recommended Gas Source	Air/Oxygen
Pressure Manifold	40cm H ₂ O pressure pop off

COMPONENTS & COMPOSITION

Predominant materials	High Density Polyethylene; Polypropylene; PTFE
Materials Not Present	Contains no natural rubber Latex, PVC or Phthalates (DEHP, DBP BBP)
Manufacturing Mode	Noninvasive device; Produced in CWE
Compatible With	OPT31X Range
Destruction	Incineration, or according to hospital protocol

REGULATORY

Classification	AU-IIa, EU-IIa, Canada-II, USA-II
Country Of Origin	New Zealand
Intended Use	Single patient use; Max 7 days
EU Notified Body	TÜV SÜD Product Service GmbH CE0123

