

DOW CONFIDENTIAL

Do not reproduce or distribute without written permission from The Dow Chemical Company.

Dow Microbial Control
Dow Europe GmbH
Bachtobelstrasse 3
CH – 8810 Horgen
Switzerland



Microbial Control

Microbiology Laboratory Report

Bacterial and Fungal In-can Challenge test according to ISO 11930:2012.

EVLY PHARMA KOZMETİK SANAYİ VE TİCARET LİMİTED ŞİRKETİ
ŞERİFALİ MAH. SÖYLEŞİ SK. -
SİTESİ - BLOK NO: 41A
ÜMRANIYE / İSTANBUL
Turkey

Contact Person: Mr. Senol Kocabas

Report No.: **BIO 21-PL0883**

Date:
07 January 2021

Laboratory Microbiologists:
Independent Laboratory

Dow TS&D Contact: Paul
Wood

Summary / Recommendations.

The sample both preserved was all free from microbial contamination as received.

The sample labelled

THE PUREST SOLUTIONS EXFOLIATING FACIAL PEELING - preserved by
manufacturer

All achieved a pass in the ISO 11930-2012 test procedure.

Did not exhibit sufficient protection to meet the pass criteria of the ISO 11930-2012 test.

DOW CONFIDENTIAL - Do not share without permission

Do not reproduce or distribute without written permission from The Dow Chemical Company.

1. Introduction

The aim of the conducted test was to evaluate antimicrobial protection of a skin cosmetic product produced by EVLY PHARMA KOZMETİK from Turkey. The test was performed according to ISO 11930:2012.

2. Materials tested

Sample of skin product produced by EVLY PHARMA KOZMETİK, Istanbul Turkey were received for testing. Sample was unpacked in the laboratory and internal laboratory number was assigned to sample.

Sample identification:

THE PUREST SOLUTIONS EXFOLIATING FACIAL PEELING - Preserved by manufacturer

Sterility Test

To assess the initial contamination a Z-streak was applied with sterile cotton-tipped applicator on Tryptic Soy Agar (TSA) petri dishes to check for bacterial contamination (incubation 30°C for upto 7 days) and on Malt Extract Agar (MEA + Gentamycin) to check for fungal contamination (incubation 28°C for upto 7 days). Gentamycin is an antibiotic to suppress bacterial growth.

Table 1

B = bacteria, F = filamentous fungi (mould), Y = yeast

| pH | Sample Description | B | F | Y |
|------|---|---|---|---|
| 5.34 | THE PUREST SOLUTIONS EXFOLIATING FACIAL PEELING | 0 | 0 | 0 |

Rating Legend:

| Plating Results | Score |
|----------------------|-------|
| No detectable growth | 0 |
| 1–10 colonies | 1 |
| 11–100 colonies | 2 |
| 101–1000 colonies | 3 |
| >1000 colonies | 4 |

Test methodology

The test was run separately for each micro-organism. For each product 20 ml of each formulation was dispensed into five sterile containers, one per each strain of micro-organism. In accordance with ISO standard 11930:2012 requirements the following strains of test organisms were used:

- *Pseudomonas aeruginosa* (ATCC 9027)
- *Staphylococcus aureus* (ATCC 6538)
- *Escherichia coli* (ATCC 8739)
- *Candida albicans* (ATCC 10231)
- *Aspergillus brasiliensis* (ATCC 16404)

The initial concentration of inoculum per micro-organism was as follows:

- *Pseudomonas aeruginosa* – 2.45×10^5 CFU/ml
- *Staphylococcus aureus* – 1.30×10^5 CFU/ml
- *Escherichia coli* – 1.75×10^5 CFU/ml
- *Candida albicans* – 6.05×10^4 CFU/ml
- *Aspergillus brasiliensis* – 6.50×10^4 CFU/ml

At T0, T7, T14 and T28 1 g of each tested product was taken for assessment of the survival of test microorganisms. 1 g of each sample was mixed with 9 ml of Dey-Eagle Neutralizing broth to neutralize the activity of the biocides. Neutralized samples were then serially diluted (10^0 to 10^{-4}) and plated on TSA for bacteria, PDA for *A. brasiliensis*, and SDA for *C. albicans*. The plates were then incubated at $(32.50 \pm 2.50)^\circ\text{C}$ for 48h to 72h for the bacteria and *C. albicans* and at $(22.50 \pm 2.50)^\circ\text{C}$ for 3 to 5 days for *A. brasiliensis*.

3. Enumeration of the microorganism and interpretation of the results

Plates with the growth of 30 to 300 colonies were taken for enumeration of bacteria and *C. albicans* and 15 to 150 colonies for *A. brasiliensis*. Then the numbers of surviving microorganisms in test samples were calculated and demonstrated as CFU/ml.

For each sample the log of reduction of total viable count was calculated. Based on the log reduction the preservation efficacy of each tested sample was evaluated according to the following criteria:

Log reduction values

| Microorganisms | Bacteria | | | <i>C. albicans</i> | | | <i>A. brasiliensis</i> | |
|----------------|------------------|------------------------|------------------------|--------------------|-----------------------|-----------------------|------------------------|-----------------------|
| Sampling time | T7 | T14 | T28 | T7 | T14 | T28 | T14 | T28 |
| Criteria A | ≥ 3 | ≥ 3 and NI* | ≥ 3 and NI | ≥ 1 | ≥ 1 and NI | ≥ 1 and NI | ≥ 0 | ≥ 1 and NI |
| Criteria B | Not performed | ≥ 3 | ≥ 3 and NI* | Not performed | ≥ 1 | ≥ 1 and NI | ≥ 0 | ≥ 0 and NI |

*No increase in the count from the previous contact time

The ISO 11930:2012 standard describes two criteria of product protection:

- Criterion A, whereby the formulation is protected against microbial proliferation that may present a potential risk for the user and no additional factors are considered;
- Criterion B, whereby the level of protection is acceptable if the risk analysis demonstrates the existence of control factors not related to the formulation indicating that the microbiological risk is tolerable for the cosmetic product.

Test results presented below include the log reduction values for each sample and microorganism and also the pass/fail data for both criteria A and B.

| Bacterial strains | Initial TVC (CFU/ml) | Log reduction values | | | Criteria A result | Criteria B result | Final result |
|---------------------------------|----------------------|----------------------|-----|-----|-------------------|-------------------|--------------|
| | | T7 | T14 | T28 | | | |
| <i>Escherichia coli</i> | 2.45×10^5 | >5 | >5 | >5 | Passed | Passed | Passed |
| <i>Staphylococcus aureus</i> | 1.30×10^5 | >5 | >5 | >5 | Passed | Passed | |
| <i>Pseudomonas aeruginosa</i> | 1.75×10^5 | >5 | >5 | >5 | Passed | Passed | |
| <i>Candida albicans</i> | 6.50×10^4 | >5 | >5 | >5 | Passed | Passed | |
| <i>Aspergillus brasiliensis</i> | 6.50×10^4 | >5 | >5 | >5 | Passed | Passed | |

NOTICE: No freedom from any patent owned by Seller or others is to be inferred. Because use conditions and applicable laws may differ from one location to another and may change with time, Customer is responsible for determining whether products and the information in this document are appropriate for Customer's use and ensuring that Customer's workplace and disposal practices are in compliance with applicable laws and other governmental enactments. Seller assumes no obligation or liability for the information in this document. NO WARRANTIES ARE GIVEN; ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED.