


phenox

pRESET Thrombectomy Device

PRES-5-40

5 mm x 40 mm

ID 0.021" 2 mm



A	B	C	D	E	F
Shaft Diameter (mm)	Shaft Length (mm)	Insertion Wire Length (mm)	Marker Distance (mm)	Working Length (mm)	Min. Vessel Diameter (mm)
5	40	2000	1000	40	2

REF PRES-5-40 LOT 481888 2026-04-27 1 x pRESET STANDARD

pRESET STANDARD PRES-5-40 LOT 481888

(01)04280123760888 1726042710481888

phenox STANDARD phenox GmbH www.phenox.net
 LOT PRES-5-40 481888 2026-04-27
 (01)04280123760888 1726042710481888
 phenox STANDARD phenox GmbH www.phenox.net
 LOT PRES-5-40 481888 2026-04-27
 (01)04280123760888 1726042710481888



IMPORTADO Y DISTRIBUIDO POR:
 BIOSOL COLOMBIA SAS
 CR. 80 # 136-54
 BOGOTÁ DC
 TEL. 38011113
 REGISTRO SANITARIO
 (DVMA 01504-012281)




pRESET LITE Thrombectomy Device

PRE-LT-4-20

4 mm x 20 mm ID 0.0165" / 0.017" 1.5 mm

LITE



A	B	C	D	E	F
Shaft Diameter (mm)	Shaft Length (mm)	Insertion Wire Length (mm)	Marker Distance (mm)	Working Length (mm)	Min. Vessel Diameter (mm)
4	20	2000	1000	40	1.5

REF PRE-LT-4-20 LOT DV4153 2023-12-31 1 x pRESET LITE

pRESET LITE PRE-LT-4-20 LOT DV4153

(01)04280123760888 1726042710481888

phenox LITE phenox GmbH www.phenox.net

LOT PRE-LT-4-20 481888 2023-12-31

(01)04280123760888 1726042710481888

phenox LITE phenox GmbH www.phenox.net

LOT PRE-LT-4-20 481888 2023-12-31

(01)04280123760888 1726042710481888

pRESET LITE Thrombectomy Device PRE-LT-4-20 4 mm x 20 mm LOT DV4153 2023-12-31

SAMPLE Not for human use!