


Bonnets

Bonnets · Disposable




REF see table >>

The bonnets are the most important connection to the baby.

Only with these bonnets you can realize a perfect fixation of the prongs and masks. There are different sizes available. The correct size is measured by our measuring tape (REF 2150). Material: 95% Polyamid, 5% Elastan

technical data & facts

Size	XXS	XS	S	M	L	XL	XXL	XXXL
REF	1213-10	1214-10	1215-10	1216-10	1217-10	1218-10	1219-10	1220-10
[cm] Head circumference	17 - 19	19 - 21	21 - 23	23 - 25,5	25,5 - 28	28 - 30	30 - 33	33 - 36
Packaging unit	10	10	10	10	10	10	10	10



REF see table >>

Fixation for disposable bonnets: Medifoam or Minifoam

20 pieces of Medifoam respectively Minifoam are included in REF1000 (Medijet) respectively REF4000 (Miniflow).

Optionally you can order both of them separately:

Medifoam or Minifoam

Fixes Medijet or Miniflow on the bonnet

Medifoam	REF 1030	Packaging unit	10
Minifoam	REF 4030	Packaging unit	10

Bonnets · Reusable



REF see table >>

technical data & facts

Size	XS	S	M	L	XL	XXL	XXXL	XXXXL
REF	1214	1215	1216	1217	1218	1219	1220	1221
[cm] Head circumference	19 - 21	21 - 23	23 - 25,5	25,5 - 28	28 - 30	30 - 33	33 - 36	36 - 40
Packaging unit	1	1	1	1	1	1	1	1

Fixation pillow

Fixes the generator more stable on the cap

REF 2018	Packaging unit	10
----------	----------------	----

Measuring tape



REF 2150

The measuring tape is the best way to choose the correct size.

Packaging unit	20
----------------	----